

LZI Methadone Metabolite (EDDP) Enzyme Immunoassay IVD For In Vitro Diagnostic Use Only



REF 0190 (100/37.5 mL R₁/R₂ Kit)
0191 (1000/375 mL R₁/R₂ Kit)



Lin-Zhi International, Inc.

Intended Use

The Lin-Zhi International, Inc. (LZI) Methadone Metabolite (EDDP) Enzyme Immunoassay is intended for the qualitative and semi-quantitative determination of methadone metabolite in human urine at a cutoff value of 300 ng/mL. The assay is designed for prescription use with a number of automated clinical chemistry analyzers.

The assay provides only a preliminary analytical result. A more specific alternative analytical chemistry method must be used in order to obtain a confirmed analytical result. Gas or Liquid Chromatography/Mass Spectrometry (GC/MS or LC/MS) are the preferred confirmatory methods (1, 2). Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary test result is positive.

Summary and Explanation of Test

Methadone is a synthetic diphenylheptanonylamine opioid that has similar analgesic activity and potency as morphine when administered parenterally. However, unlike morphine, it reliably retains its effectiveness when given orally, and tolerance and physical dependency develop slowly (3, 4). Although methadone is prescribed to relieve chronic pain, its primary, medical application, however, is the detoxification and/or maintenance treatment of narcotic or heroin addiction (3-6). The abuse potential of methadone is comparable to that of morphine due to its similar pharmacological activity (3, 5, 7).

Methadone is readily absorbed from the gastrointestinal tract when ingested, and metabolized extensively in the liver. Initial N-demethylation results in normethadone metabolite, which rapidly undergoes cyclization followed by dehydration to form 2-ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine, commonly known as EDDP. Further N-demethylation yields a secondary metabolite, 2-ethyl-5-methyl-3, 3-diphenyl-1-pyrroline (EMDP). The metabolites are secreted in urine or bile along with unchanged drug (9). Current methadone specific immunoassays can only detect the parent drug and are subject to false positives from adulterated samples for drug of abuse testing or false negatives from urine samples from persons that are fast metabolizers of methadone. As a result, confirmation of the presence of EDDP by GC/MS is often required. EDDP (methadone metabolite) immunoassay detection will ease compliance testing and rule out the possibility of urine adulteration during unsupervised collections.

Assay Principle

The LZI Methadone Metabolite assay is a homogeneous enzyme immunoassay ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent (9). Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity. In the absence of drug in the sample, methadone metabolite-labeled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when drug is present in the sample, antibody binds to the free drug; the unbound methadone metabolite-labeled G6PDH then exhibits its maximal enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at a 340 nm primary wavelength.

Reagents Provided

Antibody/Substrate Reagent (R₁): Contains mouse monoclonal anti-methadone metabolite antibody, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide (0.09 %) as a preservative.


Enzyme-drug Conjugate Reagent (R₂): Contains glucose-6-phosphate dehydrogenase (G6PDH) labeled with methadone metabolite in buffer with sodium azide (0.09 %) as a preservative.

Calibrators and Controls*

*Calibrators and controls are sold separately and contain negative human urine with sodium azide as a preservative.

METHADONE METABOLITE (EDDP) Calibrators		REF
Negative Calibrator		0001
Low Calibrator: Contains 150 ng/mL methadone metabolite		0192
Cutoff Calibrator: Contains 300 ng/mL methadone metabolite		0193
Intermediate Calibrator: Contains 600 ng/mL methadone metabolite		0194
High Calibrator: Contains 1000 ng/mL methadone metabolite		0195
METHADONE METABOLITE (EDDP) Controls		REF
Level 1 Control: Contains 225 ng/mL methadone metabolite		0197
Level 2 Control: Contains 375 ng/mL methadone metabolite		0198

Precautions and Warning

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagent contains sodium azide as a preservative, which may form explosive compounds in metal drain lines. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up. See National Institute for Occupational Safety and Health Bulletin: Explosive Azide Hazards (10).
- Do not use the reagents beyond their expiration dates.
-  For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

Reagent Preparation and Storage

The reagents are ready-to-use. No reagent preparation is required. All assay components should be refrigerated at 2-8°C when not in use.

Specimen Collection and Handling

Urine samples may be collected in plastic or glass containers. Some plastics may absorb drugs. Use of plastics such as polyethylene is recommended (11). Use fresh urine specimens for the test. If the sample cannot be analyzed immediately, it may be refrigerated at 2-8°C for up to seven days. For longer storage, keep sample frozen at -20°C and then thaw before use. Studies have shown EDDP analytes in urine are stable at -20°C up to 6 months (12). Samples should be at room temperature (18-25°C) for testing. Samples with high turbidity should be centrifuged before analysis. Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and forward both samples to the laboratory for testing. Handle all urine specimens as if they are potentially infectious.

Instrument

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting samples, mixing reagents, measuring enzyme rates at a 340 nm primary wavelength and timing the reaction accurately can be used to perform this homogeneous immunoassay. Performance characteristics presented in this package insert have been validated on the Hitachi 717. If other instruments are used, performance will need to be validated by the laboratory (13, 14).

Assay Procedure

Analyzers with the specifications indicated above are suitable for performing this homogeneous enzyme immunoassay. Refer to the specific parameters used for each analyzer before performing the assay. Typical assay parameters used for the Hitachi 717 analyzer include a 20 µL sample, 200 µL of antibody reagent (R₁), and 75 µL of enzyme conjugate reagent (R₂) at 37°C incubation temperature, 30-35 reading frames, and a 340 nm primary wavelength. For qualitative analysis, use the 300 ng/mL as the cutoff calibrator. For semi-quantitative analysis, use all five calibrators. Recalibration should be performed after reagent bottle change or if there is a change in calibrators or reagent lot. Two levels of controls are also available for monitoring the cutoff level: use the 225 ng/mL and 375 ng/mL controls for the 300 ng/mL cutoff.

Calibration and Quality Control

Good laboratory practices recommend the use of at least two levels of control specimens (one positive and one negative control near the cutoff) to ensure proper assay performance. Controls should be run with each new calibration and after specific maintenance or troubleshooting procedures as detailed in the instrument system manual. Each laboratory should establish its own control frequency. If any trends or sudden change in control value are observed, review all operating parameters, or contact LZI technical support for further assistance. Laboratories should comply with all federal, state, and local laws, as well as all guidelines and regulations.

Results

Note: A preliminary positive test result does not necessarily mean a person took illegal drugs and a negative test result does not necessarily mean a person did not take illegal drugs. There are a number of factors that influence the reliability of drug tests.

Qualitative: The cutoff calibrator, which contains 300 ng/mL of methadone metabolite, is used as a reference for distinguishing preliminary positive from negative samples. A sample with a change in absorbance (ΔA/min) equal to or greater than that obtained with the cutoff calibrator is considered positive. A sample with a change in absorbance (ΔA/min) lower than that obtained with the cutoff calibrator is considered negative.

Semi-Quantitative: The semi-quantitative mode is for purposes of (1) enabling laboratories to determine an appropriate dilution of the specimen for verification by a confirmatory method such as GC/MS, LC/MS or (2) permitting laboratories to establish quality control procedures. When an approximation of concentration is required, a calibration curve can be established with five calibrators. The concentration of methadone metabolite in the sample may then be estimated from the calibration curve.

Limitations

1. A preliminary positive result from the assay indicates only the presence of methadone metabolite. The test is not intended for quantifying this single analyte in samples.
2. A preliminary positive result does not necessarily indicate drug abuse.
3. A negative result does not necessarily mean a person did not take illegal drugs.
4. Care should be taken when reporting results as numerous factors (e.g., fluid intake, endogenous or exogenous interferents) may influence the urine test result.
5. Preliminary positive results should be confirmed by other affirmative, analytical chemistry methods (e.g., chromatography), preferably GC/MS or LC/MS.
6. The test is designed for use with human urine only.
7. The test is not for therapeutic drug monitoring.

Typical Performance Characteristics

The results shown below were performed with a single Hitachi 717 automated clinical chemistry analyzer.

Precision:

Qualitative analysis: The three calibrators and two levels of controls were evaluated. Typical results (mA/min) are as follows:

Concentration	Within Run (N=21)			Run-to-Run (N=12)		
	Mean	SD	% CV	Mean	SD	% CV
0 ng/mL	297.9	2.2	0.7 %	296.7	2.2	0.7 %
225 ng/mL	354.0	2.4	0.7 %	357.3	2.9	0.8 %
300 ng/mL	374.9	2.7	0.7 %	377.2	2.8	0.5 %
375 ng/mL	390.9	3.0	0.8 %	390.6	2.7	0.6 %
1000 ng/mL	449.4	3.4	0.8 %	452.5	4.2	0.9 %

Semi-quantitative analysis: The concentrations of the cutoff level and the two levels of controls were determined with reference curves from five calibrators. Typical results (ng/mL) are as follows:

Concentration	Within Run (N=21)			Run-to-Run* (N=12)		
	Mean	SD	% CV	Mean	SD	% CV
225 ng/mL	227.9	8.8	3.9 %	227.2	8.2	3.6 %
300 ng/mL	304.9	13.7	4.5 %	298.4	10.8	3.6 %
375 ng/mL	382.8	13.0	3.4 %	371.7	11.8	3.2 %

*Run-to-run data taken in 3 weeks

Sensitivity: Sensitivity, defined as the lowest concentration that can be differentiated from negative urine with 95 % confidence, was tested to be 15 ng/mL.

Accuracy: One hundred and thirty nine (139) clinical urine specimens were tested with the LZI EDDP EIA; 48 samples were found positive, and 91 samples were negative. All positive samples were confirmed with GC/MS.

Cutoff Value (300 ng/mL)	GC/MS	LZI EDDP EIA	% Agreement with Predicate
# Positive Samples	48	48	100 %
# Negative Samples	91	91	100 %
Total # of Samples	139	139	N/A

In addition to the above study, 20 diluted clinical samples with methadone metabolite GC/MS concentration ranging from 125 ng/mL to 1360 ng/mL were evaluated with the current EIA. Of the 13 samples with methadone metabolite GC/MS values greater than the cutoff (ranging from 328 to 1360 ng/mL), 12 were found positive. One discrepant sample was found at the cut-off borderline. The seven samples with methadone metabolite GC/MS values below the cutoff (125 ng/mL to 269 ng/mL) had negative EIA results.

Analytical Recovery: Analytical recovery was evaluated by spiking known concentrations of EDDP to negative urine samples.

In qualitative analysis, the assay correctly identified spiked samples containing more than 300 ng/mL of methadone metabolite (n=25, spiked levels equal to or higher than Level 2 Control, 375 ng/mL) as positive, and those containing less than 300 ng/mL of methadone metabolite (n=25, spiked levels equal to or less than the Level 1 Control, 225 ng/mL) as negative.

For semi-quantitative analysis, the average recovery for samples at each concentration, ranging from 30 ng/mL to 900 ng/mL (five samples at each concentration) of methadone metabolite is summarized in the following table:

Expected Value (ng/mL)	Observed Value (ng/mL)	% Recovery	SD	% CV
30	32.6	108.8 %	3.0	9.2 %
60	61.3	102.1 %	5.7	9.3 %
120	124.4	103.7 %	6.7	5.4 %
180	171.6	95.3 %	15.5	9.1 %
225	213.8	95.0 %	16.1	7.5 %
375	395.2	105.4 %	12.7	3.2 %
400	476.6	95.3 %	6.2	1.3 %
600	568.8	94.8 %	15.1	3.7 %
750	699.0	93.2 %	38.8	5.5 %
900	863.2	95.9 %	35.0	4.3 %

Specificity: Various potentially interfering substances were tested for cross-reactivity with the assay. Test compounds were spiked into the drug-free urine calibrator matrix to various concentrations and evaluated against the cutoff calibrator.

The table below lists the concentration of each test compound that gave a response approximately equivalent to that of the cutoff calibrator (as positive) or the maximal concentration of the compound tested that gave a response below the response of the cutoff calibrator (as negative).

Structurally Related Methadone Metabolite Compounds:

Compound	Target [] (µg/mL)	% Cross-Reactivity
EDDP	0.3	Positive
EMDP	550	Positive
(-) α-Methadol	60	Positive
Methadone	200	Positive
LAAM	100	Negative
nor-LAAM	10	Negative

Structurally Unrelated Pharmacological Compounds:

Compound	Target [] (µg/mL)	% Cross-Reactivity
Acetaminophen	1000	Negative
Acetylsalicylic Acid	1000	Negative
Amitriptyline	50	Negative
Amobarbital	1000	Negative
Amphetamine	1000	Negative
Benzoyllecgonine	1000	Negative
Bupropion	1000	Negative
Caffeine	1000	Negative
Chlorpheniramine	15	Negative
Chlorpromazine	20	Negative
Cocaine	1000	Negative
Codeine	1000	Negative
Dextromethorphan	50	Negative
Ecgonine	1000	Negative
Ephedrine	1000	Negative
Imipramine	20	Negative
Lidocaine	1000	Negative
Meperidine	50	Negative
Methamphetamine	1000	Negative
Methaqualone	1000	Negative
Morphine	1000	Negative
Nortriptyline	50	Negative
Oxazepam	1000	Negative
Phencyclidine	1.5	Negative
Ranitidine	1000	Negative
Secobarbital	1000	Negative
Valproic Acid	1000	Negative

It is possible that other substances and/or factors not listed above may interfere with the test and cause false positive results.

Bibliography:

1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph **73**, (1986).
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, **53**(69): 11970 (1988).
3. Goodman, L.S., and Gilman, A., The Pharmacological Basis of Therapeutics, 4th edition, 380, The MacMillan Co., (1970).
4. Roper, J.D., Garside, D., and Goldberger, B.A., "Opiates", in Contemporary Practice in Clinical Toxicology, 2nd edition, Leslie M. Shaw, editor-in-chief, AACCC, (2000).
5. Lacy, C., Armstrong, L.L., Lippy, R.J., and Lance, L.L., Drug Information Handbook. Hudson, OH: Lexi-Comp, (1993).

Bibliography, continued:

6. Katzung, B.G., Basic and Clinical Pharmacology, 6th edition, Appleton & Lange, Norwalk, CT. (1995).
7. Physician's Desk Reference, 54th edition, 2711-2713, Medical Economics Company, Montvale, NJ. (2000).
8. Baselt, R.C., and Carvey, R.H., Disposition of Toxic Drugs and Chemicals in Man, 4th edition, Chemical Toxicology Institute, Foster City, CA. (1995).
9. Rubenstein, K.E., Schneider, R.S., and Ullman, E.F., Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun.*, **47**: 846 (1972).
10. Sodium Azide. National Institute for Occupational Safety (NIOSH). Pocket Guide to Chemical Hazards. Third Printing, September 2007. Available online at: <https://www.cdc.gov/niosh/npg/default.html>
11. Yahya, A.M., McElnay, J.C., and D'Arcy, P.F., Drug absorption to glass and plastics, *Drug Metabol Drug Interact.*, **6**(1):1-45 (1988).
12. Gonzales, E., et al., Stability of pain-related medications, metabolites, and illicit substances in urine, *Clinica Chimica Acta*, **416**:80-85 (2013).
13. Nichols, J. Instrument Validation: The Road to Success. CLN's Lab 2004: From Basic to Advanced Series. 14-16 (2004).
14. CDRH Guidance for Industry and FDA Staff: Replacement Reagent and Instrument Family Policy (2003).

Additions, deletions, or changes are indicated by a change bar in the margin.
For technical assistance please call: (408) 970-8811

Manufacturer:

Lin-Zhi International, Inc.

2945 Oakmead Village Court

Santa Clara, CA 95051

USA

Tel: (408) 970-8811

Fax: (408) 970-9030

www.lin-zhi.com

EC REP

Authorized European

Rep. within the EU:

CEpartner4U

Esdoornlaan 13

3951 DB Maarn

The Netherlands

www.cepartner4u.eu



© June 2019 Rev. 10

Printed in USA