



Fentanyl (qualitative), AU5800

System Reagent: C68808

The information provided in this application sheet is intended as a supplement to the package insert.

Refer to the package insert for information on intended use, reagent storage, and additional performance data.

Parameters		Specific Test Parameters				
General	LIH	ISE	HbA1c		Calculated Test	Range
Test Name:		FEN	<	>	Type: Urine	Operation: Yes
Sample Volume	15	μL	Dilution	0	μL	OD Limit
Pre-Dilution Rate	1		Diluent Bottle	#		Min.OD: -2.0000 Max.OD: 3.0000
Rgt. Volume	R1(R1-1)	120	μL	Dilution	0	μL
	R1-2		μL	Dilution		μL
	R2(R2-1)	45	μL	Dilution	10	μL
Common Rgt. Type	None		Name	None		Dynamic Range Low: -999999.9 High: 999999.9
Wavelength	Pri	340	nm	Sec.	410	nm
Method	Fixed		Correlation Factor A	1.0		B: 0
Reaction Slope	+		Factor for Maker A	1		B: 0
Measuring Point 1 1 st	14		Last	19		Onboard Stability Period: Day Hour
Measuring Point 2 1 st			Last			LIH Influence Check: Lipemia Icterus Hemolysis
Linearity Limit						
Lag Time Check						

Parameters		Specific Test Parameters				
General	LIH	ISE	HbA1c		Calculated Test	Range
Test Name:		FEN	<	>	Type: Urine	
Value/Flag:	Flag					
Specific Ranges:	From		Level To		Low	High
	Sex	Year	Month	Year	Month	Low High
<input type="checkbox"/> 1.	#	#	#	#	#	# #
<input type="checkbox"/> 2.	#	#	#	#	#	# #
<input type="checkbox"/> 3.	#	#	#	#	#	# #
<input type="checkbox"/> 4.	#	#	#	#	#	# #
<input type="checkbox"/> 5.	#	#	#	#	#	# #
<input type="checkbox"/> 6.	#	#	#	#	#	# #
7.	Standard demographics					# #
8.	Not within expected values					# #
Panic Value	Low	#	High	#	Unit	Decimal Places: 1

Parameters		Calibration Parameters				
Calibrators	Calibration Specific					
General	ISE					
Test Name:		FEN	<	>	Type: Urine	Cuvette:
		<input type="checkbox"/> Use Serum Cal.				
Calibration Type:	AB		Formula:	Y=AX+B		Counts: 2
<Calibrator Parameters>		Range				Slope Check: None
Calibrator	OD	Conc	Low	High	Allowance Range Check	
Point 1:	#	100.0*	-999999	9999999	<input type="checkbox"/> Reagent Blank	
Point 2:					<input type="checkbox"/> Calibration	
Point 3:					Advanced Calibration Operation: Yes	
Point 4:					Interval (RB/ACAL): Lot/Lot	
Point 5:					<input type="checkbox"/> Lot Calibration	
Point 6:					Stability	
Point 7:					Reagent Blank	14 Day 0 Hour
Point 8:					Calibration	14 Day 0 Hour
Point 9:					MB Type Factor: 1-Point Calibration Point	
Point 10:					<input type="checkbox"/> with Conc-0	

User Defined

* The cutoff is normalized to 100. Positive samples are ≥ 100 and are flagged with a (P). LZI Norfentanyl Calibrator Ref No.: C68815.

Lin-Zhi International, Inc.

Santa Clara, CA 95051

Tel: (408) 970-8811

Fax: (408) 970-9030

www.lin-zhi.com



Fentanyl (wash setting), AU5800

System Reagent: C68808

The information provided in this application sheet is intended as a supplement to the package insert.

Refer to the package insert for information on intended use, reagent storage, and additional performance data.

Miscellaneous						
Checked Tests	Contamination Parameters			Data Check Parameters		
Contamination Prevention	Carry-Over Prevention (Type Changes)			Carry-Over Prevention (Test)		
Test Name	Pre-Dispense Wash Count			Post-Dispense Wash Count		
	Detergent-1	Detergent-2	Water	Detergent-1	Detergent-2	Water
1. Fentanyl	0	0	0	0	0	3
2. X	x	x	x	x	x	x
3. X	x	x	x	x	x	x
4. X	x	x	x	x	x	x
5. X	x	x	x	x	x	x
6. X	x	x	x	x	x	x
7. X	x	x	x	x	x	x
8. X	x	x	x	x	x	x
9. X	x	x	x	x	x	x
10. X	x	x	x	x	x	x
11. X	x	x	x	x	x	x
12. X	x	x	x	x	x	x
13. X	x	x	x	x	x	x
14. X	x	x	x	x	x	x
15. X	x	x	x	x	x	x
16. X	x	x	x	x	x	x
17. X	x	x	x	x	x	x
18. X	x	x	x	x	x	x
19. X	x	x	x	x	x	x
20. X	x	x	x	x	x	x