



MONOCLONAL BECKMAN BUPRENORPHINE ENZYME IMMUNOASSAY - QUALITATIVE

Catalog # A53684

Application	AU400	AU400*	AU480
	AU600	AU640	AU680

Please note that these are preliminary parameters only and validation studies have not been completed by the supplier of this assay. Analytical laboratories are solely responsible for the validation of these parameters for their diagnostic and screening applications. These parameters are for informational use only.

Specific Test Parameters									
General									
LIH	ISE	Range							
Test Name:	BUP	<	>	Type:	URINE	Operation:	Yes		
Sample:	Volume:	12.5	Dilution:		μL	Pre-Dilution Rate:	1		
Reagents:	R ₁ Volume:	100	Dilution:		μL	Min OD:		Max OD:	
	R ₂ Volume:	25	Dilution:	20	μL	L:	-2.00	H:	2.50
Wavelength:	Pri:	340	Sec:	410		Reagent OD Limit:			
Method:	FIXED 1					First L:	-2.00	First H:	2.50
Reaction Slope:	+					Last L:	-2.00	Last H:	2.50
Measuring Point 1:	First:	13	Last:	16		Dynamic Range:			
Measuring Point 2:	First:		Last:			L:	-9999999.9	H:	9999999.9
Linearity:						Correlation Factor:			
No Lag Time:						A:	1.0000	B:	§
						On-board stability period:	#		

Specific Test Parameters									
General									
LIH	ISE	Range							
Test Name:	BUP	<	>	Type:	URINE				
Value/Flag:	FLAG		Level L:	-99999.9		Level H:	¥		
Normal Ranges:									
	Sex	Year	Month	Year	Month	L	H		
<input type="checkbox"/>	1.	#	#	#	#	#	#		
<input type="checkbox"/>	2.	#	#	#	#	#	#		
<input type="checkbox"/>	3.	#	#	#	#	#	#		
<input type="checkbox"/>	4.	#	#	#	#	#	#		
<input type="checkbox"/>	5.	#	#	#	#	#	#		
<input type="checkbox"/>	6.	#	#	#	#	#	#		
<input type="checkbox"/>	7.	None Selected					-9999999.9	9999999.9	
<input type="checkbox"/>	8.	Out of Range					#	#	
Panic Value	L	H				Unit:	Decimal Places: 1		
	#	#							

Calibration Specific									
General									
ISE									
Test Name:	BUP	<	>	Type:	URINE				
Calibration Type:	AB	Formula:	Y = AX + B		Counts:	2	Process:	OD	
	Cal. No.	OD	CONC		Factor/OD-L	Factor/OD-H			
Point 1:	*		100**		-9999999	9999999			
Point 2:									
Point 3:									
Point 4:									
Point 5:									
Point 6:									
Point 7:									
1-Point Cal. Point:				<input type="checkbox"/> With CONC-0	Advanced Calibration:		No		
MB Type Factor:	N/A				Calibration Stability Period:	#			

* Cutoff Calibrator
 ** % of Cutoff Calibrator
 # User Defined
 ¥ Specific Test Parameters, Level H: VALUE, Option 1 enter 99999
 FLAG, Option 2 enter 99.9
 FLAG, Option 3 enter -0.1
 § Correlation Factor B for results:
 Option 1, enter 0
 Option 2, enter 0
 Option 3, enter -100



MONOCLONAL BECKMAN BUPRENORPHINE ENZYME IMMUNOASSAY - SEMI-QUANTITATIVE

Catalog # A53684

Application	AU400	AU400*	AU480
	AU600	AU640	AU680

Please note that these are preliminary parameters only and validation studies have not been completed by the supplier of this assay. Analytical laboratories are solely responsible for the validation of these parameters for their diagnostic and screening applications. These parameters are for informational use only.

Specific Test Parameters									
General									
LIH	ISE	Range							
Test Name:	BUP	<	>	Type:	URINE	Operation:	Yes		
Sample:	Volume	12.5	Dilution		μL	Pre-Dilution Rate:	1		
Reagents:	R ₁ Volume	100	Dilution		μL	Min OD		Max OD	
	R ₂ Volume	25	Dilution	20	μL	L	-2.00	H	2.50
Wavelength:	Pri	340	Sec	410		Reagent OD Limit:			
Method:	FIXED 1					First L	-2.00	First H	2.50
Reaction Slope:	+					Last L	-2.00	Last H	2.50
Measuring Point 1:	First	13	Last	16		Dynamic Range:			
Measuring Point 2:	First		Last			L	#	H	#
Linearity:						Correlation Factor:			
No Lag Time:						A	1.000	B	0.000
						On-board stability period:	#		

Specific Test Parameters									
General									
LIH	ISE	Range							
Test Name:	BUP	<	>	Type:	URINE				
Value/Flag:	VALUE		Level L:	-99999.9		Level H:	99999.9		
Normal Ranges:									
<input type="checkbox"/>	Sex	Year	Month	Year	Month	L	H		
<input type="checkbox"/>	1.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	2.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	3.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	4.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	5.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	6.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	7.	None Selected					-99999.9	99999.9	
<input type="checkbox"/>	8.	Out of Range					#	#	
Panic Value	L	H				Unit:	Decimal Places: 1		

Calibration Specific									
General									
LIH	ISE	Range							
Test Name:	BUP	<	>	Type:	URINE				
Calibration Type:	6AB	Formula:	Polygonal		Counts:	2	Process:	CONC	
	Cal. No.	OD	CONC		Factor/OD-L	Factor/OD-H			
Point 1:	#		0		-2.000	2.500			
Point 2:	#		5		-2.000	2.500			
Point 3:	#		10		-2.000	2.500			
Point 4:	#		20		-2.000	2.500			
Point 5:	#		40		-2.000	2.500			
Point 6:	#		100		-2.000	2.500			
Point 7:									
1-Point Cal. Point:				<input type="checkbox"/> With CONC-0	Advanced Calibration:	No			
MB Type Factor:	N/A				Calibration Stability Period:	#			

User Defined

¥ Specific Test Parameters, Level H:
 VALUE, Option 1 enter 99999
 FLAG, Option 2 enter 99.9
 FLAG, Option 3 enter -0.1

§ Correlation Factor B for results:
 Option 1, enter 0
 Option 2, enter 0
 Option 3, enter -100