

LZI Fentanyl II Enzyme Immunoassay

IVD For In Vitro Diagnostic Use Only



REF 0570 (100/37.5 mL R₁/R₂ Kit)
0571 (1000/375 mL R₁/R₂ Kit)



Lin-Zhi International, Inc.

Intended Use

The LZI Fentanyl II Enzyme Immunoassay is intended for the qualitative determination of norfentanyl in human urine at a cutoff value of 5 ng/mL when calibrated against norfentanyl. The assay is designed for prescription use with a number of automated clinical chemistry analyzers.

The assay provides only a preliminary analytical result. A more specific alternative chemical confirmatory method (e.g., gas or liquid chromatography and mass spectrometry) must be used to obtain a confirmed analytical result (1, 2). Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary test result is positive.

Summary and Explanation of Test

Fentanyl is an important opioid analgesic used widely in surgical operations and is a controlled substance (3). Fentanyl is most commonly encountered in the form of patches applied to the skin, as “lollipops” which can be dissolved in the mouth through the mucous membrane, or can be administered intravenously. It is 50-100 times stronger than morphine (4, 5) and cases of fentanyl abuse via intravenous injection, inhalation, oral, or nasal applications have been previously reported (6). Fentanyl is used in the treatment of acute and chronic pain, usually in patients who no longer respond to high doses of less potent opioids such as morphine or oxycodone. Due to its potency and wide availability as a prescribed drug, fentanyl has been abused and misused by health professionals, pain management patients, and recreational abusers (7).

Due to its short elimination half-life and approximately 90 % metabolism, fentanyl is difficult to detect in urine (8). Fentanyl undergoes extensive hepatic biotransformation to metabolites coming from hydrolysis, N-dealkylation, or hydroxylation reactions (9). In an intravenous dose of fentanyl, up to 85 % is excreted in urine over a three- to four- day period with 0.4-6 % eliminated as unchanged fentanyl and 26-55 % eliminated as the norfentanyl metabolite (10).

Fentanyl analogs also have high potency analgesic activities. Numerous reports have been published with modified fentanyl-related compounds abused as designer drugs (11-13).

Other recently available fentanyl analogs associated with abuse and severe intoxication include butyryl fentanyl and 4-fluorobutyryl fentanyl (14-18).

Assay Principle

The LZI Fentanyl II Enzyme Immunoassay is a homogeneous enzyme immunoassay ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent (19). The drug-labeled G6PDH conjugate is traceable to a commercially available fentanyl standard and referred to as fentanyl-labeled G6PDH conjugate. Enzyme activity decreases upon binding to the antibody, and the fentanyl concentration in the sample is measured in terms of enzyme activity. In the absence of norfentanyl in the sample, fentanyl-labeled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when norfentanyl is present in the sample, antibody would bind to free norfentanyl; the unbound fentanyl-labeled G6PDH then exhibits its maximal enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm.

Reagents Provided

Antibody/Substrate Reagent (R₁): Contains a mouse monoclonal anti-fentanyl antibody, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide (0.09 %) as a preservative.

Enzyme-drug Conjugate Reagent (R₂): Contains fentanyl-labeled glucose-6-phosphate dehydrogenase (G6PDH) in buffer with sodium azide (0.09 %) as a preservative.

Calibrators and Controls*

*Calibrators and Controls are sold separately and contain negative human urine with sodium azide as a preservative.

| NORFENTANYL Calibrators | REF |
|---|------|
| Cutoff Calibrator: Contains 5 ng/mL norfentanyl | 0313 |

| NORFENTANYL Controls | REF |
|--|------|
| Level 1 Control: Contains 3.75 ng/mL norfentanyl | 0317 |
| Level 2 Control: Contains 6.25 ng/mL norfentanyl | 0318 |

Precautions and Warning

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagent contains sodium azide as a preservative, which may form explosive compounds in metal drain lines. When disposing such reagents or wastes, always flush with a large volume of water to prevent azide build-up. See National Institute for Occupational Safety and Health Bulletin: Explosive Azide Hazards (20).
- Do not use the reagents beyond their expiration dates.
- For USA: Federal law restricts this device to sale by or on the order of a physician.

Reagent Preparation and Storage

The reagents are ready to use. No reagent preparation is required. All assay components should be refrigerated at 2-8°C when not in use.

Specimen Collection and Handling

Use fresh urine specimens for the test. If the sample cannot be analyzed immediately, it may be refrigerated at 2-8°C for up to four weeks (21) or at room temperature for up to four weeks (21, 22). For longer storage, keep sample frozen at -20°C and then thaw before use. Studies have shown norfentanyl samples in urine are stable at -20°C for up to six months (23). Sample stability claims were established by experimental data by the manufacturer or based on reference literature and only for the temperatures/time frames as stated in the method sheet. It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory.

Samples should be equilibrated to room temperature (18-25°C) for testing. Samples with high turbidity should be centrifuged before analysis.

Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and both samples should be forwarded to a laboratory for testing.

Handle all urine specimens as if they are potentially infectious.

Instrument

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting sample, mixing reagents, measuring enzyme rates at 340 nm and timing the reaction accurately can be used to perform this homogeneous immunoassay.

Performance characteristics presented in this package insert have been validated on the Beckman Coulter AU480 automated clinical analyzer.

Assay Procedure

Typical assay parameters used for the Beckman Coulter AU480 analyzer include a 20 µL sample, 120 µL of antibody reagent (R₁), 45 µL of enzyme conjugate reagent (R₂), 12-16 reading frame, FIXED method, and 340 nm primary wavelength.

For qualitative analysis, use the 5 ng/mL as the cutoff calibrator.

Recalibration should be performed after reagent bottle change or a change in calibrators or reagent lot. Two levels of controls are also available for monitoring the cutoff level: 3.75 and 6.25 ng/mL.

Calibration and Quality Control

Good laboratory practices recommend the use of at least two levels of control specimens (one positive and one negative control near the cutoff) to ensure proper assay performance. Controls should be run with each new calibration and after specific maintenance or troubleshooting procedures as detailed in the instrument system manual. Each laboratory should establish its own control frequency. If any trends or sudden change in control value are observed, review all operating parameters, or contact LZI technical support for further assistance. Laboratories should comply with all federal, state, and local laws, as well as all guidelines and regulations.

Results

Note: A preliminary positive test result does not necessarily mean a person took a specific drug and a negative test result does not necessarily mean a person did not take a specific drug. There are a number of factors that influence the reliability of drug tests.

Qualitative: The cutoff calibrator, which contains 5 ng/mL of norfentanyl, is used as a reference for distinguishing positive from negative samples. A sample with a change in absorbance (Δ MAU) equal to or greater than that obtained with the cutoff calibrator is considered positive. A sample with a change in absorbance (Δ MAU) lower than that obtained with the cutoff calibrator is considered negative.

Limitations

1. Boric Acid at 1% w/v may cause false negative results.
2. Dextromethorphan may cause false positive results at concentrations greater than 40,000 ng/mL.
3. A preliminary positive result from this assay indicates only the presence of norfentanyl. The test is not intended for quantifying this single analyte in samples.
4. A preliminary positive result does not necessarily indicate drug abuse.
5. A negative result does not necessarily mean a person did not take illegal drugs.
6. Care should be taken when reporting results, as numerous factors (e.g., fluid intake, endogenous or exogenous interferents) may influence the urine test result.
7. Preliminary positive results must be confirmed by other affirmative, analytical methods (e.g., chromatography), preferably GC/MS or LC/MS.
8. The test is designed for use with human urine only.
9. This test should not be used for therapeutic drug monitoring.

Typical Performance Characteristics

The results shown below were performed with a single Beckman Coulter AU480 automated chemistry analyzer.

Precision:

Qualitative analysis: The following concentrations were evaluated. Typical qualitative results (measured by Δ OD, mAU) are as follows:

| 5 ng/mL Cutoff | | Within Run (N = 22) | | Run-to-Run (N = 88) | |
|----------------|-------------|---------------------|------------------|---------------------|-------------------|
| Concentration | % of Cutoff | # Samples | EIA Result | # Samples | EIA Result |
| 0 ng/mL | 0 % | 22 | 22 Neg | 88 | 88 Neg |
| 1.25 ng/mL | 25 % | 22 | 22 Neg | 88 | 88 Neg |
| 2.5 ng/mL | 50 % | 22 | 22 Neg | 88 | 88 Neg |
| 3.75 ng/mL | 75 % | 22 | 22 Neg | 88 | 88 Neg |
| 5 ng/mL | 100 % | 22 | 13 Neg/ 9 Pos | 88 | 60 Neg/ 28 Pos |
| 6.25 ng/mL | 125 % | 22 | 22 Pos | 88 | 88 Pos |
| 7.5 ng/mL | 150 % | 22 | 22 Pos | 88 | 88 Pos |
| 8.75 ng/mL | 175 % | 22 | 22 Pos | 88 | 88 Pos |
| 10 ng/mL | 200 % | 22 | 22 Pos | 88 | 88 Pos |

Accuracy: One hundred (100) unaltered clinical urine specimens were tested with the LZI Fentanyl II Enzyme Immunoassay and confirmed with LC/MS. Specimens with a norfentanyl concentration greater than or equal to 5 ng/mL by LC/MS are defined as positive, and specimens with a norfentanyl concentration below 5 ng/mL by LC/MS are defined as negative in the table below. Near cutoff samples are defined as \pm 50 % of the cutoff value. The correlation results are summarized as follows:

Qualitative Accuracy Study:

| NFEN Results 5 ng/mL Cutoff | Neg | < 50 % below the cutoff | Near Cutoff Neg | Near Cutoff Pos | > 50 % above the cutoff |
|--------------------------------|-----|----------------------------|--------------------|--------------------|----------------------------|
| Positive | 0 | 1* | 8* | 10 | 40 |
| Negative | 20 | 19 | 2 | 0 | 0 |

The following table summarizes the results for the qualitative discordant samples:

| Sample # | Norfentanyl LC/MS (ng/mL) | LC/MS Pos/Neg Result | AU480 EIA Qualitative Result (mAU) | AU480 EIA Qualitative Cutoff Rate (mAU) | LZI FEN II EIA Pos/Neg Result |
|----------|---------------------------------|----------------------------|---|--|-------------------------------------|
| 37* | 1.5 | - | 85.9 | 83.0 | + |
| 41* | 2.7 | - | 111.3 | 83.0 | + |
| 43* | 3.0 | - | 207.9 | 83.0 | + |
| 44* | 3.0 | - | 107.7 | 83.0 | + |
| 45* | 3.3 | - | 124.7 | 83.0 | + |
| 46* | 3.5 | - | 169.6 | 83.0 | + |
| 47* | 3.8 | - | 204.6 | 83.0 | + |
| 48* | 3.9 | - | 113.6 | 83.0 | + |
| 49* | 4.2 | - | 263.1 | 83.0 | + |

* Discrepant below the cutoff concentration (0 ng/mL – 4.9 ng/mL)

These samples contained levels of fentanyl that contributed to the false positive result.

Specificity: Various potentially interfering substances were tested for cross-reactivity with the assay. Test compounds were spiked into a drug free–urine pool to various concentrations and evaluated against the cutoff calibrator.

The following table lists the concentration of each test compound that gave a response approximately equivalent to that of the cutoff calibrator (as positive) or the maximal concentration of the compound tested that gave a response below the response of the cutoff calibrator (as negative).

Compounds tested at high concentration (100,000 ng/mL) with results below the cutoff value were listed as Not Detected (ND).

Fentanyl and Metabolites:

| Compound | Concentration Tested (ng/mL) | % Cross- Reactivity | Result |
|-------------|------------------------------------|------------------------|----------|
| Fentanyl | 3.8 | 131.58 % | Positive |
| Norfentanyl | 5 | 100.00 % | Positive |

Structurally Related Compounds:

| Compound | Concentration Tested (ng/mL) | % Cross-Reactivity |
|--------------------------------------|---------------------------------|--------------------|
| 4-Fluoro-Isobutyryl Fentanyl | 20.0 | 25.00 % |
| 9-Hydroxy Risperidone | 100,000 | ND |
| Acetyl Fentanyl | 7.0 | 71.43 % |
| Acetyl Norfentanyl | 100.0 | 5.00 % |
| Acryl Fentanyl | 4.0 | 125.00 % |
| Alfentanil | 100,000 | ND |
| Butyryl Fentanyl | 6.0 | 83.33 % |
| Butyryl Norfentanyl | 40.0 | 12.50 % |
| Carfentanil Oxalate | 100,000 | ND |
| Cis- d,l 3-Methyl Fentanyl | 8.0 | 62.50 % |
| Cyclopropyl Fentanyl | 3.2 | 156.25 % |
| Cyclopropyl Norfentanyl | 25.0 | 20.00 % |
| Despropionyl Fentanyl (4-ANPP) | 100,000 | ND |
| Furanyl Fentanyl | 5.5 | 90.91 % |
| Furanyl Norfentanyl | 180.0 | 2.78 % |
| (±) β -Hydroxy ThioFentanyl | 5.0 | 100.00 % |
| Isobutyryl Fentanyl | 15.0 | 33.33 % |
| Isobutyryl Norfentanyl | 500.0 | 1.00 % |
| Labetalol Hydrochloride | 100,000 | ND |
| Methoxyacetyl Fentanyl | 3.5 | 142.86 % |
| MT-45 | 100,000 | ND |
| N-benzyl Furanyl Norfentanyl | 11.0 | 45.45 % |
| N-benzyl Para-fluoro Norfentanyl | 4.0 | 125.00 % |
| Norcarfentanil Oxalate | 100,000 | ND |
| Ocfentanil | 3.8 | 131.58 % |
| Para-fluoro Butyryl Fentanyl (p-FBF) | 4.5 | 111.11 % |
| Para-fluoro Fentanyl | 3.2 | 156.25 % |
| Remifentanil | 100,000 | ND |
| Risperidone | 100,000 | ND |
| Sufentanil | 100,000 | ND |
| Thienyl Fentanyl | 4.0 | 125.00 % |
| Thiofentanyl | 3.2 | 156.25 % |
| Trans- d,l 3-Methyl Fentanyl | 6.0 | 83.33 % |
| Trazodone | 100,000 | ND |
| U-47700 | 100,000 | ND |
| Valeryl Fentanyl | 70.0 | 7.14 % |
| ω -1-Hydroxy Fentanyl | 300.0 | 1.67 % |

Structurally Unrelated Compounds:

| Compound | Spiked [] (ng/mL) | Spiked Norfentanyl Concentration | | |
|----------------------|-----------------------|----------------------------------|-----------------------|-----------------------|
| | | 0 ng/mL | 3.75 ng/mL Control | 6.25 ng/mL Control |
| Acetaminophen | 100,000 | ND | Neg | Pos |
| 6-Acetylmorphine | 100,000 | ND | Neg | Pos |
| Acetylsalicylic Acid | 100,000 | ND | Neg | Pos |
| Amitriptyline | 100,000 | ND | Neg | Pos |
| Amlodipine Besylate | 100,000 | ND | Neg | Pos |
| Amoxicillin | 100,000 | ND | Neg | Pos |
| d-Amphetamine | 100,000 | ND | Neg | Pos |
| Atorvastatin | 100,000 | ND | Neg | Pos |
| Benzoyllecgonine | 100,000 | ND | Neg | Pos |
| Buprenorphine | 100,000 | ND | Neg | Pos |
| Bupropion | 100,000 | ND | Neg | Pos |
| Caffeine | 100,000 | ND | Neg | Pos |
| Carbamazepine | 100,000 | ND | Neg | Pos |
| Cetirizine | 100,000 | ND | Neg | Pos |
| Chlorpheniramine | 100,000 | ND | Neg | Pos |
| Chlorpromazine | 100,000 | ND | Neg | Pos |
| Clomipramine | 100,000 | ND | Neg | Pos |

Structurally Unrelated Compounds, continued:

| Compound | Spiked [] (ng/mL) | Spiked Norfentanyl Concentration | | |
|---|--------------------|----------------------------------|----------------------------|----------------------------|
| | | 0 ng/mL (ng/mL) | 3.75 ng/mL Control (ng/mL) | 6.25 ng/mL Control (ng/mL) |
| Codeine | 100,000 | ND | Neg | Pos |
| Desipramine | 100,000 | ND | Neg | Pos |
| Dextromethorphan | 40,000 | 0.01 % | Pos | Pos |
| Diphenhydramine | 100,000 | ND | Neg | Pos |
| Duloxetine | 100,000 | ND | Pos | Pos |
| Fluoxetine | 100,000 | ND | Neg | Pos |
| Fluphenazine | 100,000 | ND | Neg | Pos |
| Gabapentin | 100,000 | ND | Neg | Pos |
| Hydrocodone | 100,000 | ND | Neg | Pos |
| Hydromorphone | 100,000 | ND | Neg | Pos |
| Ibuprofen | 100,000 | ND | Neg | Pos |
| Imipramine | 100,000 | ND | Neg | Pos |
| Lisinopril | 100,000 | ND | Neg | Pos |
| Losartan | 100,000 | ND | Neg | Pos |
| Loratadine | 100,000 | ND | Neg | Pos |
| MDA (3,4-methylene-dioxyamphetamine) | 100,000 | ND | Neg | Pos |
| MDEA | 100,000 | ND | Neg | Pos |
| MDMA (3,4-methylene-dioxyamphetamine) | 100,000 | ND | Neg | Pos |
| Meperidine | 100,000 | ND | Neg | Pos |
| Metformin | 100,000 | ND | Neg | Pos |
| Metoprolol | 100,000 | ND | Neg | Pos |
| Methadone | 100,000 | ND | Neg | Pos |
| <i>d</i> -Methamphetamine | 100,000 | ND | Neg | Pos |
| Morphine | 100,000 | ND | Neg | Pos |
| Nalmefene | 100,000 | ND | Neg | Pos |
| Nicotine | 100,000 | ND | Neg | Pos |
| Nortriptyline | 100,000 | ND | Neg | Pos |
| Omeprazole | 100,000 | ND | Neg | Pos |
| Oxazepam | 100,000 | ND | Neg | Pos |
| Oxycodone | 100,000 | ND | Neg | Pos |
| Oxymorphone | 100,000 | ND | Neg | Pos |
| Phenobarbital | 100,000 | ND | Neg | Pos |
| (1 <i>S</i> ,2 <i>S</i>)-(+)-Pseudoephedrine | 100,000 | ND | Neg | Pos |
| Quetiapine | 100,000 | ND | Neg | Pos |
| Ranitidine | 100,000 | ND | Neg | Pos |
| Salbutamol (Albuterol) | 100,000 | ND | Neg | Pos |
| Sertraline | 100,000 | ND | Neg | Pos |
| THC-COOH (11-Nor-Delta-9-THC-9-carboxylic acid) | 100,000 | ND | Neg | Pos |
| <i>l</i> -Thyroxine | 100,000 | ND | Neg | Pos |
| Tramadol | 100,000 | ND | Neg | Pos |
| Zolpidem | 100,000 | ND | Neg | Pos |
| Phencyclidine | 100,000 | ND | Neg | Pos |

It is possible that other substances and/or factors not listed above may interfere with the test and cause false positive results.

Endogenous and Preservatives Compound Interference Study:

Various potentially interfering endogenous and preservative substances were tested for interference with the assay. Test compounds were split into three portions each and either left un-spiked or spiked to a norfentanyl concentration of either 3.75 or 6.25 ng/mL (the negative and positive control concentrations, respectively). These samples were then evaluated in qualitative mode.

Boric Acid (1 % w/v) were found to cause interference with the assay. Results are summarized in the following table:

Endogenous and Preservatives Compound Interference Study

| Endogenous or Preservative Substance | Spiked [] (mg/dL) | Spiked Norfentanyl Concentration | | |
|--------------------------------------|--------------------|----------------------------------|--------------------|--------------------|
| | | 0 ng/mL | 3.75 ng/mL Control | 6.25 ng/mL Control |
| Acetone | 1000 | Neg | Neg | Pos |
| Ascorbic Acid | 500 | Neg | Neg | Pos |
| Bilirubin | 2 | Neg | Neg | Pos |
| Biotin | 0.5 | Neg | Neg | Pos |
| Boric Acid | 1000 | Neg | Neg | Neg |
| Calcium Chloride (CaCl2) | 300 | Neg | Neg | Pos |
| Citric Acid (pH 3) | 200 | Neg | Neg | Pos |
| Creatinine | 500 | Neg | Neg | Pos |
| Ethanol | 1000 | Neg | Neg | Pos |
| Galactose | 10 | Neg | Neg | Pos |
| γ-Globulin | 500 | Neg | Neg | Pos |
| Glucose | 3000 | Neg | Neg | Pos |
| Hemoglobin | 300 | Neg | Neg | Pos |
| HSA | 500 | Neg | Neg | Pos |
| Human Urine (pooled) | N/A | Neg | Neg | Pos |
| β-hydroxybutyric Acid | 100 | Neg | Neg | Pos |
| Oxalic Acid | 100 | Neg | Neg | Pos |
| Potassium Chloride | 1000 | Neg | Neg | Pos |

Endogenous and Preservatives Compound Interference Study, continued:

| Endogenous or Preservative Substance | Spiked [] (mg/dL) | Spiked Norfentanyl Concentration | | |
|--------------------------------------|--------------------|----------------------------------|--------------------|--------------------|
| | | 0 ng/mL | 3.75 ng/mL Control | 6.25 ng/mL Control |
| Riboflavin | 7.5 | Neg | Neg | Pos |
| Sodium Azide | 1000 | Neg | Neg | Pos |
| Sodium Chloride | 1000 | Neg | Neg | Pos |
| Sodium Fluoride | 1000 | Neg | Neg | Pos |
| Sodium Phosphate | 300 | Neg | Neg | Pos |
| Urea | 6000 | Neg | Neg | Pos |
| Uric Acid | 10 | Neg | Neg | Pos |
| LZI Urine-Based Calibrator Buffer | N/A | Neg | Neg | Pos |

The following endogenous compounds which showed interference at ±25 % of cutoff concentrations were then spiked into negative urine and at ±50 % of cutoff concentrations (2.5 ng/mL and 7.5 ng/mL) for the assay. Interference was still observed with Boric Acid at 1 % w/v. Results are summarized in the following table:

| Endogenous Substance | Spiked [] (mg/dL) | Spiked Norfentanyl Concentration | | |
|----------------------|--------------------|----------------------------------|---------------------------|---------------------------|
| | | 0 ng/mL (ng/mL) | 2.5 ng/mL Control (ng/mL) | 7.5 ng/mL Control (ng/mL) |
| Boric Acid | 1000 | Neg | Neg | Neg |

pH Interference Study: pH 3 to pH 11 was tested for interference with the assay. Each pH level was split into three portions each and either left un-spiked or spiked to a norfentanyl concentration of either 3.75 or 6.25 ng/mL (the negative and positive control concentrations, respectively). These samples were then evaluated in qualitative mode. No pH interference was observed.

| pH | Spiked Norfentanyl Concentration | | |
|-------|----------------------------------|--------------------|--------------------|
| | 0 ng/mL | 3.75 ng/mL Control | 6.25 ng/mL Control |
| pH 3 | Neg | Neg | Pos |
| pH 4 | Neg | Neg | Pos |
| pH 5 | Neg | Neg | Pos |
| pH 6 | Neg | Neg | Pos |
| pH 7 | Neg | Neg | Pos |
| pH 8 | Neg | Neg | Pos |
| pH 9 | Neg | Neg | Pos |
| pH 10 | Neg | Neg | Pos |
| pH 11 | Neg | Neg | Pos |
















Specific Gravity: Samples ranging in specific gravity from 1.000 to 1.027 were split into three portions each and either left un-spiked or spiked to a norfentanyl concentration of either 3.75 or 6.25 ng/mL (the negative and positive control concentrations, respectively). These samples were then evaluated in qualitative mode. No interference was observed.

| Specific Gravity | Spiked Norfentanyl Concentration | | |
|------------------|----------------------------------|--------------------|--------------------|
| | 0 ng/mL | 3.75 ng/mL Control | 6.25 ng/mL Control |
| 1.000 | Neg | Neg | Pos |
| 1.003 | Neg | Neg | Pos |
| 1.005 | Neg | Neg | Pos |
| 1.008 | Neg | Neg | Pos |
| 1.010 | Neg | Neg | Pos |
| 1.012 | Neg | Neg | Pos |
| 1.015 | Neg | Neg | Pos |
| 1.018 | Neg | Neg | Pos |
| 1.020 | Neg | Neg | Pos |
| 1.022 | Neg | Neg | Pos |
| 1.025 | Neg | Neg | Pos |
| 1.027 | Neg | Neg | Pos |

Open-Vial Reagent and Calibrator/Control Stability: Real-time data for open-vial reagent and calibrator/control stability studies at Cold Temperature (2-8°C) have been carried out up to Day 377. Results from open-vial studies indicate that degradation is minimal up to Day 377, and, based on the real-time data, suggests an open-vial stability of up to 12 months. Open-vial reagents and calibrators/controls should be stored at 2-8°C for maximum shelf life.

Closed-Vial Calibrator/Control Stability: Real-time data for closed-vial calibrator/control stability studies at Cold Temperature (2-8°C) have been carried out up to Day 377. Results from closed-vial studies indicate that degradation is minimal at Cold Temperature (2-8°C) up to Day 377 in comparison to Day 1. Closed-vial calibrators/controls should be stored at 2-8°C for maximum shelf life.

Symbols Used

| | | | |
|---|--|---|---|
|  | Authorized Representative |  | Manufacturer |
|  | Biological Risks |  | R ₁ , Antibody/ Substrate Reagent |
|  | CE Mark |  | R ₂ , Enzyme- Drug Conjugate Reagent |
|  | Consult Instructions for Use |  | Reference Number |
|  | Contents |  | Safety Data Sheet |
|  | Global Trade Item Number |  | Temperature Limits |
|  | <i>In Vitro</i> Diagnostic medical device | T.K. | Test Kit Number |
|  | Lot Number |  | Use-by Date |

Bibliography:

- Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
- Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, **23**(82):7920-7970 (2017).
- Ikeda R., Pelton C. Diversion programs for impaired physicians. *West J Med*, **152**: 617-621 (1990).
- Poklis A., Fentanyl: a review for clinical and analytical toxicologists. *Clinical Toxicology*, **33**, 439-447 (1995).
- Jaffe J.H., Martin W.R., Opioid analgesics and antagonists. In: Pharmacological Basis of Therapeutics, 7th Ed. Gilman AG, Goodman LS, Rall TW, Murad F, eds. New York: MacMillan, 517 (1985).
- Ellenhorn M.D., Schonwald S., Ordog G., Wasserberger J, Fentanyl, in: Ellenhorn's Medical Toxicology, Williams & Wilkins, Baltimore, Philadelphia, London, Paris, Bangkok, Buenos Aires, Hong Kong, Munich, Sydney, Tokyo, Wroclaw, 416-420 (1997).
- Atluri S., Sudarshan G., Manchikanti L., Assessment of the trends in medical use and misuse of opioid analgesics from 2004 to 2011. *Pain Physician*, **17**, E119-E128 (2014).
- McClain DA, Hug CC Jr., Intravenous fentanyl kinetics. *Clin Pharmacol Ther* **28**:106-114 (1980).
- Hug C.C., Murphy M.R., Anesthesiology, **55** (1981) 369. T. Tateishi, A.J.J. Wood, F.P. Guengerich and M. Wood, *Biochem. Pharmacol.*, **50** (1995) 1921.
- Baselt R.C., Disposition of Toxic Drugs and Chemicals in Man, 10th edition. Biomedical Publications: Seal Beach, California; 846-849 (2014).
- Baum R.M., New variety of street drugs poses growing problem. *Chem. Eng News* Sep **9**:7-16 (1985).
- Suzuki S., Inoue, T., New abused drug "China White". *Forensic Toxicol. News* **5**:5-12 (1987).
- Heagy J.A., Identification of 4-propionoxy-4-phenyl-N-methylpiperidine. *Microgram* **15**: 181-186 (1982).
- Poklis J., Poklis A., Wolf C., Hathaway C., Arbefeville E., Chrostowski L., Devers K., Hair L., Mainland M., Merves M., and Pearson J., *Journal of Analytical Toxicology*, **40**:703-708 (2016).
- Backberg M., Beck O., Jonsson K.H., Helander A., Opioid intoxications involving butyryl fentanyl, 4-fluorobutyrfentanyl, and fentanyl from the Swedish STRIDA project. *Clinical toxicology*, **53**:609-617 (2015).
- EMCDDA. European Drug Report 2014: Trends and developments 2014. Available at: http://www.emcdda.europa.eu/attachements.cfm/att_228272_EN_TDAT14001ENN.pdf.
- EMCDDA. (2015) New psychoactive substances in Europe. An update from the EU Early Warning System (March 2015). Available at: http://www.emcdda.europa.eu/attachements.cfm/att_235958_EN_TD0415135ENN.pdf.
- Helander A., Backberg M., Hulten P., Al-Saffer Y., Beck O., Detection of new psychoactive substance use among emergency room patients: results from Swedish STRIDA project. *Forensic Science International*, **243**:23-29(2014).
- Rubenstein K.E., Schneider R.S., Ullman E.F., Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, **47**: 846 (1972).

Bibliography, continued:

- Sodium Azide. National Institute for Occupational Safety (NIOSH) Pocket Guide to Chemical Hazards Third Printing, September 2007. Available online at: <https://www.cdc.gov/niosh/npg/default.html>
- Hammargren, W.R. and G.L. Henderson. Analyzing normetabolites of the Fentanyl by Gas Chromatography/Electron Capture Detection. *J. of Anal. Tox.* **12**:183-191 (1988).
- Thevis, M., Geyer, H., Bahr, D., and W. Schänzer. Identification of fentanyl, alfentanil, sufentanil, remifentanyl and their major metabolites in human urine by liquid chromatography/tandem mass spectrometry for doping control purposes. *Eur. J. Mass Spectrom.* **11**:419-427 (2005).
- Gonzales, E., Ng, G., Pesce, A., West, C., West, R., Mikel, C., Latyshev, S., and P. Almazan. Stability of pain-related medications, metabolites, and illicit substances in urine. *Clinica Chimica Acta* **416**:80-85 (2013).

|Additions, deletions, or changes are indicated by a change bar in the margin.

For technical assistance please call: (408) 970-8811

 **Manufacturer:**
Lin-Zhi International, Inc.
2945 Oakmead Village Court
Santa Clara, CA 95051 USA
Tel: (408) 970-8811
Fax: (408) 970-9030
www.lin-zhi.com

© December 2020 Rev. 2