



Reagent Name Hydrocodone 300 C/O REF C68823 DxC 500 AU Urine Settings

Calibrator Name LZI Hydrocodone 300 Qualitative Calibrator Ref C68830

The information provided in this application sheet is intended as a supplement to the package insert.

Refer to the package insert for information on intended use, reagent storage, and additional performance data.

TEST CONFIGURATION & CHEMISTRY DETAILS

Assay Name	Test	Rev	Discipline	Chemistry
Test ID	HYD300		Calculated Result	<input type="checkbox"/>
LIS Code	HYD300		Result Type	Qualitative
UNITS AND RANGE SETTINGS				
Use Settings from	None	Units	None	Decimal Places
Test Kind	General	Revision	01	<input checked="" type="checkbox"/> Multi Reagent Switch
Reagent Name	HYD	Reagent ID	249	<input type="checkbox"/> FSE Test
Region	<input checked="" type="checkbox"/> US	<input checked="" type="checkbox"/> OUS	<input checked="" type="checkbox"/> AP	<input type="checkbox"/> JP <input checked="" type="checkbox"/> EU <input type="checkbox"/> Other
ABB Name	HYDIG	Parameter Long Name	Hydrocodone 300 C/O C68823 HYDIG Urine	

GENERAL PARAMETERS

SAMPLE VOLUME		Sample Volume	12	µL	Dilution	0	µL	REACTION OD LIMIT					
REAGENT VOLUME		Pre-Dilution Rate	1		REACTION BLANK OD LIMIT		Min. OD	-2.0000	Max OD	3.0000			
R1 (R1-1)	120	µL	Dilution	0	µL	First: Low	-2.0000	High	3.0000				
R1-2		µL	Dilution		µL	Last: Low	-2.0000	High	3.0000				
R2 (R2-1)	45	µL	Dilution	10	µL	ANALYTICAL MEASURING RANGE		Low	-999999.9	High	999999.9		
WAVELENGTH		Primary	340	nm	Secondary	410	nm	MANUFACTURER FACTOR		A	1	B	0
METHOD		FIXED		REAGENT ONBOARD STABILITY		8		Days	0		Hours		
REACTION SLOPE		+		LIH INFLUENCE CHECK		<input type="checkbox"/> Perform LIH Check		Lipemia		+			
MEASURING POINT		Point 1: First	14	Last	18	Icterus		+		+			
Linearity Limit				Last		Hemolysis		+		+			
Lag Time Check		<input type="checkbox"/> Perform Lag Time Check											

CALIBRATION PARAMETERS

Base Unit	Decimal Place	Unit1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4
None	1	None	0	None	0	None	0	None	0
CALIBRATOR SPECIFIC					CALIBRATION OD AND CONCENTRATION PARAMETERS				
Calibration Type		AB	Counts	2	<input type="checkbox"/> Use highest calibrator for Upper AMR				
Formula		Y=AX+B	MB Factor		Calibrator Name	Conc.	Factor Range Low	Factor Range High	
Calibrator Name		Add #	Positive Cutoff	100.0	#	100.0*	-9999999	9999999	
<input type="checkbox"/> SLOPE CHECK		Slope Check	Number of Levels	5	Point 1				
STABILITY AND INTERVAL		Reagent Blank Stability	8	Days	0	Hours	Point 2		
Calibration Stability		8	Days	0	Hours	Point 3			
Interval		Lot	Interval	Lot	Point 4				
<input type="checkbox"/> OD DELTA CHECK		<input type="checkbox"/> Reagent Blank	0.0000		Point 5				
		<input type="checkbox"/> Calibration	0.0000		Point 6				
					Point 7				

User Defined

PROZONE CHECK PARAMETERS



Reagent Name **Hydrocodone 300 C/O REF C68823 Dx C 500 AU Urine Settings**

Calibrator Name LZI Hydrocodone 300 Qualitative Calibrator Ref C68830

The information provided in this application sheet is intended as a supplement to the package insert.

Refer to the package insert for information on intended use, reagent storage, and additional performance data.

<input type="checkbox"/> Logic Check 1		<input type="checkbox"/> Logic Check 2		<input type="checkbox"/> Logic Check 3	
Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values
Point 1	<input type="text" value="0"/>	Point 1	<input type="text" value="0"/>	Point 1	<input type="text" value="0"/>
Point 2	<input type="text" value="0"/>	Interval	<input type="text" value="1"/>	Interval	<input type="text" value="1"/>
Point 3	<input type="text" value="0"/>	Value 1	<input type="text" value="0"/>	Value 1	<input type="text" value="0"/>
		Value 2	<input type="text" value="0"/>	Value 2	<input type="text" value="0"/>
		Value 3	<input type="text" value="0"/>		
Limit Points		Limit Points		Limit Points	
Limit 1	<input type="text" value="0"/>	Limit 1	<input type="text" value="0"/>	Limit 1	<input type="text" value="0"/>
Limit 2	<input type="text" value="27"/>	Limit 2	<input type="text" value="27"/>	Limit 2	<input type="text" value="27"/>
Check Pattern					
Pattern	<input type="text" value="Pattern 1"/>				



Reagent Name Hydrocodone 300 S/Q REF C68823 DxC 500 AU Urine Settings

Calibrator Name LZI Hydrocodone 300 Semi-Quantitative Calibrator

Ref C68807, C68831

The information provided in this application sheet is intended as a supplement to the package insert.

Refer to the package insert for information on intended use, reagent storage, and additional performance data.

TEST CONFIGURATION & CHEMISTRY DETAILS

Assay Name	Test	Rev	Discipline	Chemistry
Test ID	HYD300-		Calculated Result	<input type="checkbox"/>
LIS Code	HYD300-		Result Type	Semi-Quantitative
UNITS AND RANGE SETTINGS				
Use Settings from	None	Units	ng/mL	Decimal Places
Test Kind	General	Revision	01	<input checked="" type="checkbox"/> Multi Reagent Switch
Reagent Name	HYD	Reagent ID	249	<input type="checkbox"/> FSE Test
Region	<input checked="" type="checkbox"/> US	<input checked="" type="checkbox"/> OUS	<input checked="" type="checkbox"/> AP	<input type="checkbox"/> JP
	<input checked="" type="checkbox"/> EU	<input type="checkbox"/> Other		
ABB Name	HYD2G	Parameter Long Name	Hydrocodone 300 S/Q C68823 HYD2G Urine	

GENERAL PARAMETERS

SAMPLE VOLUME	Sample Volume	12	µL	Dilution	0	µL	REACTION OD LIMIT	Min. OD	-2.0000	Max OD	3.0000
	Pre-Dilution Rate	1					REACTION BLANK OD LIMIT	First: Low	-2.0000	High	3.0000
REAGENT VOLUME	R1 (R1-1)	120	µL	Dilution	0	µL		Last: Low	-2.0000	High	3.0000
	R1-2		µL	Dilution		µL	ANALYTICAL MEASURING RANGE	Low	150.0	High	800.0
	R2 (R2-1)	45	µL	Dilution	10	µL	MANUFACTURER FACTOR	A	1	B	0
WAVELENGTH	Primary	340	nm	Secondary	410	nm	REAAGENT ONBOARD STABILITY		8	Days	0
							LIH INFLUENCE CHECK	<input type="checkbox"/> Perform LIH Check			
METHOD		FIXED					Lipemia	+			
REACTION SLOPE		+					Icterus	+			
MEASURING POINT	Point 1: First	14		Last	18		Hemolysis	+			
	Point 2: First			Last							
Linearity Limit			%								
Lag Time Check	<input type="checkbox"/> Perform Lag Time Check										

CALIBRATION PARAMETERS

Base Unit	Decimal Place	Unit1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4	
ng/mL	1	None	0	None	0	None	0	None	0	
CALIBRATOR SPECIFIC					CALIBRATION OD AND CONCENTRATION PARAMETERS					
Calibration Type	5AB	Counts	2	<input type="checkbox"/> Use highest calibrator for Upper AMR						
Formula	Polygonal	MB Factor		Point 1	Calibrator Name	Conc.	Factor Range Low	Factor Range High		
Calibrator Name	Add #	Positive Cutoff	300.0	Point 2	#	0.0	-2.0000	3.0000		
<input checked="" type="checkbox"/> SLOPE CHECK	Slope Check	Number of Levels	5	Point 3	#	150.0	-2.0000	3.0000		
STABILITY AND INTERVAL	Reagent Blank Stability	8	Days	0	Hours	Point 4	#	300.0	-2.0000	3.0000
	Calibration Stability	8	Days	0	Hours	Point 5	#	500.0	-2.0000	3.0000
	Interval	Lot		Point 6						
	Interval	Lot		Point 7						
	OD DELTA CHECK									
	<input type="checkbox"/> Reagent Blank	0.0000								
	<input type="checkbox"/> Calibration	0.0000								

User Defined

PROZONE CHECK PARAMETERS



Reagent Name Hydrocodone 300 S/Q REF C68823 DxC 500 AU Urine Settings

Calibrator Name LZI Hydrocodone 300 Semi-Quantitative Calibrator Ref C68807, C68831

The information provided in this application sheet is intended as a supplement to the package insert.

Refer to the package insert for information on intended use, reagent storage, and additional performance data.

<input type="checkbox"/> Logic Check 1		<input type="checkbox"/> Logic Check 2		<input type="checkbox"/> Logic Check 3	
Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values
Point 1	<input type="text" value="0"/>	Point 1	<input type="text" value="0"/>	Point 1	<input type="text" value="0"/>
Point 2	<input type="text" value="0"/>	Interval	<input type="text" value="1"/>	Interval	<input type="text" value="1"/>
Point 3	<input type="text" value="0"/>	Value 1	<input type="text" value="0"/>	Value 1	<input type="text" value="0"/>
		Value 2	<input type="text" value="0"/>	Value 2	<input type="text" value="0"/>
		Value 3	<input type="text" value="0"/>		
Limit Points		Limit Points		Limit Points	
Limit 1	<input type="text" value="0"/>	Limit 1	<input type="text" value="0"/>	Limit 1	<input type="text" value="0"/>
Limit 2	<input type="text" value="27"/>	Limit 2	<input type="text" value="27"/>	Limit 2	<input type="text" value="27"/>
Check Pattern					
Pattern	<input type="text" value="Pattern 1"/>				