



URINE LZI Fentanyl III (qualitative), AU480

System Reagent: C82916 or D29739

The information provided in this application sheet is intended as a supplement to the package insert.
 Refer to the package insert for information on intended use, reagent storage, and additional performance data.

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	FEN3#	<	>	Type:	Urine	Operation:	Yes			
Sample Volume	15	μL	Dilution	0	μL	OD Limit				
Pre-Dilution Rate	1					Min OD	-2.0000	Max	2.5000	
Reagents Volume: R1 (R1-1)	65	μL	Dilution	0	μL	Reagent OD Limit:				
						First Low	-2.0000	High	2.5000	
						Last Low	-2.0000	High	2.5000	
R2 (R2-1)	25	μL	Dilution	0	μL	Dynamic Range Low	-999999.9	High	999999.9	
Wavelength: Pri.	340	▼	Sec	410	▼	Correlation Factor A	1.0	B	0	
Method:	FIXED1	▼				Factor for Maker A	1	B	0	
Reaction Slope:	+	▼				On-board stability	90	Day	0	Hour
Measuring Point 1:	First	14	Last	18		LIH Influence Check	No			
Measuring Point 2:	First		Last			Lipemia				
Linearity:		%				Icterus				
No Lag Time:						Hemolysis				

Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:	FEN3#	<	>	Type:	Urine				
Value/Flag:	Flag#	▼	Level L:	99.9#	Level H:	100.0#			
Specific Ranges:	From	To							
	Sex	Year	Month	Year	Month	Low	High	Panic Value	
								Low	High
<input type="checkbox"/>	1.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	2.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	3.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	4.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	5.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	6.	#	#	#	#	#	#	#	#
	7.	No demographics			#	#	#	#	
	8.	Not within expected values			#	#	#	#	
Unit					Decimal Places	1			

Calibration Specific												
General		ISE										
Test Name:	FEN3#	<	>	Type:	Urine					<input type="checkbox"/> Use Serum Cal.		
Calibration Type:	AB	▼	Formula:	Y = AX + B	▼	Counts:	2	Slope Check	None	▼		
<Calibrator Parameters>												
	Calibrator	OD	Conc	Factor Range	Low	High	Allowable Range Check					
Point 1:	FEN3 Cutoff Cal#		100.0*	#	#		<input type="checkbox"/> Reagent Blank					
Point 2:							<input type="checkbox"/> Calibration					
Point 3:							Advanced Calibration					
Point 4:							Operation	No	▼			
Point 5:							Interval (RB/ACAL)					
Point 6:												
Point 7:												
Point 8:												
Point 9:												
Point 10:												
<Point Cal. For Master Curve>												
	Calibrator	OD	Conc	No. of Correction Points		Use Master Curve		<input type="checkbox"/> Lot Calibration				
						OD Range	Stability					
						Low	High	Reagent Blanks	14	Day	0	Hour
Point 1								Calibration	14	Day	0	Hour
Point 2												
MB Type Factor			1-Point Calibration Point				<input type="checkbox"/> With CONC-O					

User Defined
 * The cutoff is normalized to 100. Positive samples are ≥ 100. LZI Fentanyl Qualitative Calibrator Ref No.: D26784.



URINE LZI Fentanyl III (qualitative), (wash setting), AU480

System Reagent: C82916 or D29739

The information provided in this application sheet is intended as a supplement to the package insert.

Refer to the package insert for information on intended use, reagent storage, and additional performance data.

Miscellaneous						
Checked Tests		Contamination Parameters			Data Check Parameters	
Contamination Prevention		Carry-Over Prevention (Type Changes)			Carry-Over Prevention (Test)	
Test Name	Pre-Dispense Wash Count			Post-Dispense Wash Count		
	Detergent-1	Detergent-2	Water	Detergent-1	Detergent-2	Water
1. Fentanyl	0	0	0	0	0	4
2. X	x	x	x	x	x	x
3. X	x	x	x	x	x	x
4. X	x	x	x	x	x	x
5. X	x	x	x	x	x	x
6. X	x	x	x	x	x	x
7. X	x	x	x	x	x	x
8. X	x	x	x	x	x	x
9. X	x	x	x	x	x	x
10. X	x	x	x	x	x	x
11. X	x	x	x	x	x	x
12. X	x	x	x	x	x	x
13. X	x	x	x	x	x	x
14. X	x	x	x	x	x	x
15. X	x	x	x	x	x	x
16. X	x	x	x	x	x	x
17. X	x	x	x	x	x	x
18. X	x	x	x	x	x	x
19. X	x	x	x	x	x	x
20. X	x	x	x	x	x	x