



PRODUCT CATALOG: Drugs of Abuse Urinalysis by Catalog Number

<u>REF #</u>	<u>Product Name</u>	<u>REF #</u>	<u>Product Name</u>
0001	Universal Negative Calibrator, 5 mL	0100	^{510(k)} _{Cleared} Amphetamines 500 EIA Kit, small test kit
0010	^{510(k)} _{Cleared} Phencyclidine EIA Kit, small test kit	0101	Amphetamines 500 EIA Kit, large test kit
0011	Phencyclidine EIA Kit, large test kit	0102	Amphetamines 500 Low Calibrator (MAMP 250 ng/mL)
0012	Phencyclidine Low Calibrator (12.5 ng/mL)	0103	Amphetamines 500 Cutoff Calibrator (MAMP 500 ng/mL)
0013	Phencyclidine Cutoff Calibrator (25 ng/mL)	0104	Amphetamines 500 Intermediate Calibrator (MAMP 1000 ng/mL)
0014	Phencyclidine Intermediate Calibrator (50 ng/mL)	0105	Amphetamines 500 High Calibrator (MAMP 2000 ng/mL)
0015	Phencyclidine High Calibrator (100 ng/mL)	0107	Amphetamines 500 Level 1 Control (MAMP 375 ng/mL)
0017	Phencyclidine Level 1 Control (18 ng/mL)	0108	Amphetamines 500 Level 2 Control (MAMP 625 ng/mL)
0018	Phencyclidine Level 2 Control (32 ng/mL)	0110	^{510(k)} _{Cleared} Methadone EIA Kit, small test kit
0020b	^{510(k)} _{Cleared} Opiate EIA Kit, small test kit	0111	Methadone EIA Kit, large test kit
0021b	Opiate EIA Kit, large test kit	0112	Methadone Low Calibrator (150 ng/mL)
0022	Opiate Low Calibrator (Morphine 150 ng/mL)	0113	Methadone Cutoff Calibrator (300 ng/mL)
0023	Opiate Cutoff Calibrator (Morphine 300 ng/mL)	0114	Methadone Intermediate Calibrator (600 ng/mL)
0024	Opiate Intermediate Calibrator (Morphine 600 ng/mL)	0115	Methadone High Calibrator (1000 ng/mL)
0025	Opiate High Calibrator (Morphine 1000 ng/mL)	0117	Methadone Level 1 Control (225 ng/mL)
0027	Opiate Level 1 Control (Morphine 225 ng/mL)	0118	Methadone Level 2 Control (375 ng/mL)
0028	Opiate Level 2 Control (Morphine 375 ng/mL)	0120	^{510(k)} _{Cleared} Propoxyphene EIA Kit, small test kit
0030	^{510(k)} _{Cleared} Cocaine Metabolite EIA Kit, small test kit	0121	Propoxyphene EIA Kit, large test kit
0031	Cocaine Metabolite EIA Kit, large test kit	0122	Propoxyphene Low Calibrator (150 ng/mL)
0032	Cocaine Metabolite Low Calibrator (BE 150 ng/mL)	0123	Propoxyphene Cutoff Calibrator (300 ng/mL)
0033	Cocaine Metabolite Cutoff Calibrator (BE 300 ng/mL)	0124	Propoxyphene Intermediate Calibrator (600 ng/mL)
0034	Cocaine Metabolite Intermediate Calibrator (BE 1000 ng/mL)	0125	Propoxyphene High Calibrator (1000 ng/mL)
0035	Cocaine Metabolite High Calibrator (BE 3000 ng/mL)	0127	Propoxyphene Level 1 Control (225 ng/mL)
0037	Cocaine Metabolite Level 1 Control (BE 225 ng/mL)	0128	Propoxyphene Level 2 Control (375 ng/mL)
0038	Cocaine Metabolite Level 2 Control (BE 375 ng/mL)	0130	^{510(k)} _{Cleared} Benzodiazepine EIA Kit, small test kit
0040	^{510(k)} _{Cleared} Amphetamines EIA Kit, small test kit	0131	Benzodiazepine EIA Kit, large test kit
0041	Amphetamines EIA Kit, large test kit	BZO 200	
0042	Amphetamines Low Calibrator (MAMP 500 ng/mL)	0132	Benzodiazepine Calibrator #2/Control I (Oxazepam 100 ng/mL)
0043	Amphetamines Cutoff Calibrator (MAMP 1000 ng/mL)	0133	Benzodiazepine Calibrator #3/Cutoff A (Oxazepam 200 ng/mL)
0044	Amphetamines Intermediate Calibrator (MAMP 1500 ng/mL)	0134	Benzodiazepine Calibrator #4/Cutoff III (Oxazepam 300 ng/mL)
0045	Amphetamines High Calibrator (MAMP 2000 ng/mL)	0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)
0047	Amphetamines Level 1 Control (MAMP 750 ng/mL)	BZO 300	
0048	Amphetamines Level 2 Control (MAMP 1250 ng/mL)	0132	Benzodiazepine Calibrator #2 (Oxazepam 100 ng/mL)
0070c	^{510(k)} _{Cleared} Cannabinoids (cTHC) EIA Kit, small test kit	0133	Benzodiazepine Calibrator #3/Control II (Oxazepam 200 ng/mL)
0071c	Cannabinoids (cTHC) EIA Kit, large test kit	0134	Benzodiazepine Calibrator #4/Cutoff B (Oxazepam 300 ng/mL)
0002c	Cannabinoids (cTHC) Negative Calibrator (0 ng/mL)	0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)
	cTHC 25	0136	Benzodiazepine Control IV (Oxazepam 400 ng/mL)
0072c	Cannabinoids (cTHC) Low Calibrator (12.5 ng/mL)	0140	^{510(k)} _{Cleared} Barbiturate EIA Kit, small test kit
0073c	Cannabinoids (cTHC) Cutoff Calibrator (25 ng/mL)	0141	Barbiturate EIA Kit, large test kit
0007c	Cannabinoids (cTHC) Intermediate Calibrator (37.5 ng/mL)	BARB 200	
0075c	Cannabinoids (cTHC) High Calibrator (50 ng/mL)	0142	Barbiturate Calibrator #2/Control I (Secobarbital 100 ng/mL)
0006c	Cannabinoids (cTHC) Level 1 Control (18.75 ng/mL)	0143	Barbiturate Calibrator #3/Cutoff A (Secobarbital 200 ng/mL)
0074c	Cannabinoids (cTHC) Level 2 Control (31.25 ng/mL)	0144	Barbiturate Calibrator #4/Control III (Secobarbital 300 ng/mL)
	cTHC 50	0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)
0073c	Cannabinoids (cTHC) Low Calibrator (25 ng/mL)	BARB 300	
0075c	Cannabinoids (cTHC) Cutoff Calibrator (50 ng/mL)	0142	Barbiturate Calibrator #2 (Secobarbital 100 ng/mL)
0076c	Cannabinoids (cTHC) Intermediate Calibrator (75 ng/mL)	0143	Barbiturate Calibrator #3/Control II (Secobarbital 200 ng/mL)
0077c	Cannabinoids (cTHC) High Calibrator (100 ng/mL)	0144	Barbiturate Calibrator #4/Cutoff B (Secobarbital 300 ng/mL)
0007c	Cannabinoids (cTHC) Level 1 Control (37.5 ng/mL)	0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)
0008c	Cannabinoids (cTHC) Level 2 Control (62.5 ng/mL)	0146	Barbiturate Control IV (Secobarbital 400 ng/mL)
	cTHC 100	0160	^{510(k)} _{Cleared} Ecstasy (MDMA) EIA Kit, small test kit
0075c	Cannabinoids (cTHC) Low Calibrator (50 ng/mL)	0161	Ecstasy (MDMA) EIA Kit, large test kit
0077c	Cannabinoids (cTHC) Cutoff Calibrator (100 ng/mL)	0162	Ecstasy Low Calibrator (100 ng/mL)
0078c	Cannabinoids (cTHC) Intermediate Calibrator (150 ng/mL)	0163	Ecstasy Cutoff Calibrator (500 ng/mL)
0079c	Cannabinoids (cTHC) High Calibrator (200 ng/mL)	0164	Ecstasy Intermediate Calibrator (750 ng/mL)
0076c	Cannabinoids (cTHC) Level 1 Control (75 ng/mL)	0165	Ecstasy High Calibrator (1000 ng/mL)
0009c	Cannabinoids (cTHC) Level 2 Control (125 ng/mL)	0167	Ecstasy Level 1 Control (375 ng/mL)
		0168	Ecstasy Level 2 Control (625 ng/mL)

^{510(k)}_{Cleared} = FDA 510(k) Cleared Assay, Calibrators, & Controls.
[†] Forensic Use Only
[‡] Employment and Insurance Testing Use Only
[§] Product is made to order. Please allow 10 to 12 weeks of production time

* EIA = Enzyme Immunoassay
 * All single analyte calibrators and controls are 5 mL per bottle.
 * A small test kit is 100/37.5 mL R₁/R₂
 * A large test kit is 1000/375 mL R₁/R₂



PRODUCT CATALOG: Drugs of Abuse Urinalysis by Catalog Number

<u>REF #</u>	<u>Product Name</u>	<u>REF #</u>	<u>Product Name</u>
0190b	^{510(k) cleared} Methadone Metabolite (EDDP) EIA Kit, small test kit		BUP 10
0191b	Methadone Metabolite (EDDP) EIA Kit, large test kit	0273	Norbuprenorphine Low Calibrator (5 ng/mL)
	EDDPb 100	0274	Norbuprenorphine Level 1 Control (7 ng/mL)
0202b	Methadone Metabolite 100 Low Calibrator (50 ng/mL)	0275	Norbuprenorphine Cutoff Calibrator (10 ng/mL)
0203b	Methadone Metabolite 100 Cutoff Calibrator (100 ng/mL)	0276	Norbuprenorphine Level 2 Control (13 ng/mL)
0204b	Methadone Metabolite 100 Intermediate Calibrator (250 ng/mL)	0277	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)
0205b	Methadone Metabolite 100 High Calibrator (500 ng/mL)	0278	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)
0207b	Methadone Metabolite 100 Level 1 Control (75 ng/mL)	0279	Norbuprenorphine High Calibrator (75 ng/mL)
0208b	Methadone Metabolite 100 Level 2 Control (125 ng/mL)		
	EDDPb 300	0290	^{510(k) cleared} 6-Acetylmorphine EIA Kit, small test kit
0192b	Methadone Metabolite 300 Low Calibrator (150 ng/mL)	0291	6-Acetylmorphine EIA Kit, large test kit
0193b	Methadone Metabolite 300 Cutoff Calibrator (300 ng/mL)	0299	6-Acetylmorphine Negative Calibrator (0 ng/mL)
0194b	Methadone Metabolite 300 Intermediate Calibrator (600 ng/mL)	0292	6-Acetylmorphine Low Calibrator (5 ng/mL)
0195b	Methadone Metabolite 300 High Calibrator (1000 ng/mL)	0293	6-Acetylmorphine Cutoff Calibrator (10 ng/mL)
0197b	Methadone Metabolite 300 Level 1 Control (225 ng/mL)	0294	6-Acetylmorphine Intermediate Calibrator (20 ng/mL)
0198b	Methadone Metabolite 300 Level 2 Control (375 ng/mL)	0295	6-Acetylmorphine High Calibrator (40 ng/mL)
		0297	6-Acetylmorphine Level 1 Control (7.5 ng/mL)
0220	^{510(k) cleared} Ethyl Alcohol Enzymatic Assay Kit, small test kit	0298	6-Acetylmorphine Level 2 Control (12.5 ng/mL)
0221	Ethyl Alcohol Enzymatic Assay Kit, large test kit		
0223	Ethyl Alcohol Calibrator (100 mg/dL)	0310	^{510(k) cleared} Fentanyl EIA Kit (Q), small test kit
0224	Ethyl Alcohol Level 1 Control (50 mg/dL)	0311	Fentanyl EIA Kit (Q), large test kit
0225	Ethyl Alcohol Level 2 Control (300 mg/dL)	0313	Norfentanyl Cutoff Calibrator (5 ng/mL)
		0317	Norfentanyl Level 1 Control (3.75 ng/mL)
0230	Cotinine EIA Kit, small test kit [‡]	0318	Norfentanyl Level 2 Control (6.25 ng/mL)
0231	Cotinine EIA Kit, large test kit [‡]		
0232	Cotinine Low Calibrator/Level 1 Control (250 ng/mL) [‡]	0320b	Ethyl Glucuronide EIA Kit, small test kit [†]
0233	Cotinine Cutoff Calibrator (500 ng/mL) [‡]	0321b	Ethyl Glucuronide EIA Kit, large test kit [†]
0234	Cotinine Intermediate Calibrator (1000 ng/mL) [‡]	0322b	Ethyl Glucuronide Low Calibrator (250 ng/mL) [†]
0235	Cotinine High Calibrator (2000 ng/mL) [‡]	0323b	Ethyl Glucuronide Cutoff Calibrator (1000 ng/mL) [†]
0238	Cotinine Level 2 Control (750 ng/mL) [‡]	0324b	Ethyl Glucuronide Intermediate Calibrator #1 (2000 ng/mL) [†]
		0325b	Ethyl Glucuronide Intermediate Calibrator #2 (3000 ng/mL) [†]
0240b	^{510(k) cleared} Oxycodone EIA Kit, small test kit	0326b	Ethyl Glucuronide High Calibrator (4000 ng/mL) [†]
0241b	Oxycodone EIA Kit, large test kit	0327b	Ethyl Glucuronide Level 1 Control (500 ng/mL) [†]
	OXYb 100	0328b	Ethyl Glucuronide Level 2 Control (1500 ng/mL) [†]
0252b	Oxycodone Low Calibrator (50 ng/mL)		
0242b	Oxycodone Level 1 Control (75 ng/mL)	0330	^{510(k) cleared} Opiate 2000 EIA Kit, small test kit
0243b	Oxycodone Cutoff Calibrator #1 (100 ng/mL)	0331	Opiate 2000 EIA Kit, large test kit
0244b	Oxycodone Level 2 Control (125 ng/mL)	0332	Opiate 2000 Low Calibrator (Morphine 1000 ng/mL)
0246b	Oxycodone Cutoff Calibrator #2 (300 ng/mL)	0333	Opiate 2000 Cutoff Calibrator (Morphine 2000 ng/mL)
0248b	Oxycodone Intermediate Calibrator (500 ng/mL)	0334	Opiate 2000 Intermediate Calibrator (Morphine 4000 ng/mL)
0249b	Oxycodone High Calibrator (800 ng/mL)	0335	Opiate 2000 High Calibrator (Morphine 6000 ng/mL)
	OXYb 300	0337	Opiate 2000 Level 1 Control (Morphine 1500 ng/mL)
0252b	Oxycodone Low Calibrator (50 ng/mL)	0338	Opiate 2000 Level 2 Control (Morphine 2500 ng/mL)
0243b	Oxycodone Cutoff Calibrator #1 (100 ng/mL)		
0245b	Oxycodone Level 1 Control (225 ng/mL)	0340	^{510(k) cleared} Cocaine Metabolite 150 EIA Kit, small test kit
0246b	Oxycodone Cutoff Calibrator #2 (300 ng/mL)	0341	Cocaine Metabolite 150 EIA Kit, large test kit
0247b	Oxycodone Level 2 Control (375 ng/mL)	0342	Cocaine Metabolite 150 Low Calibrator (BE 75 ng/mL)
0248b	Oxycodone Intermediate Calibrator (500 ng/mL)	0343	Cocaine Metabolite 150 Cutoff Calibrator (BE 150 ng/mL)
0249b	Oxycodone High Calibrator (800 ng/mL)	0344	Cocaine Metabolite 150 Intermediate Calibrator (BE 300 ng/mL)
		0345	Cocaine Metabolite 150 High Calibrator (BE 1000 ng/mL)
0270	^{510(k) cleared} Buprenorphine EIA Kit, small test kit	0347	Cocaine Metabolite 150 Level 1 Control (BE 112.5 ng/mL)
0271	Buprenorphine EIA Kit, large test kit	0348	Cocaine Metabolite 150 Level 2 Control (BE 187.5 ng/mL)
	BUP 5		
0272	Norbuprenorphine Level 1 Control (3 ng/mL)	0350	^{510(k) cleared} Methamphetamine EIA Kit, small test kit
0273	Norbuprenorphine Cutoff Calibrator #1 (5 ng/mL)	0351	Methamphetamine EIA Kit, large test kit
0274	Norbuprenorphine Level 2 Control (7 ng/mL)	0352	Methamphetamine Low Calibrator (250 ng/mL)
0275	Norbuprenorphine Cutoff Calibrator #2 (10 ng/mL)	0353	Methamphetamine Cutoff Calibrator (500 ng/mL)
0277	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)	0354	Methamphetamine Intermediate Calibrator (1000 ng/mL)
0278	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)	0355	Methamphetamine High Calibrator (2000 ng/mL)
0279	Norbuprenorphine High Calibrator (75 ng/mL)	0357	Methamphetamine Level 1 Control (375 ng/mL)
		0358	Methamphetamine Level 2 Control (625 ng/mL)

^{510(k) cleared} = FDA 510(k) Cleared Assay, Calibrators, & Controls.

[†] Forensic Use Only

[‡] Employment and Insurance Testing Use Only

[§] Product is made to order. Please allow 10 to 12 weeks of production time

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂

* A large test kit is 1000/375 mL R₁/R₂



PRODUCT CATALOG: Drugs of Abuse Urinalysis by Catalog Number

<u>REF #</u>	<u>Product Name</u>
0370	^{510(k) cleared} Carisoprodol Metabolite (SOMA) EIA Kit, small test kit
0371	Carisoprodol Metabolite (SOMA) EIA Kit, large test kit
0372	Meprobamate Low Calibrator (50 ng/mL)
0373	Meprobamate Cutoff Calibrator (100 ng/mL)
0374	Meprobamate Intermediate Calibrator (200 ng/mL)
0375	Meprobamate High Calibrator (400 ng/mL)
0377	Meprobamate Level 1 Control (75 ng/mL)
0378	Meprobamate Level 2 Control (125 ng/mL)
0380	^{510(k) cleared} Hydrocodone 100 EIA Kit, small test kit
0381	Hydrocodone 100 EIA Kit, large test kit
0382	Hydrocodone 100 Low Calibrator (50 ng/mL)
0383	Hydrocodone 100 Cutoff Calibrator (100 ng/mL)
0384	Hydrocodone 100 Intermediate Calibrator (150 ng/mL)
0385	Hydrocodone 100 High Calibrator (300 ng/mL)
0387	Hydrocodone 100 Level 1 Control (75 ng/mL)
0388	Hydrocodone 100 Level 2 Control (125 ng/mL)
0390	^{510(k) cleared} Hydrocodone 300 EIA Kit, small test kit
0391	Hydrocodone 300 EIA Kit, large test kit
0392	Hydrocodone 300 Low Calibrator (150 ng/mL)
0393	Hydrocodone 300 Cutoff Calibrator (300 ng/mL)
0394	Hydrocodone 300 Intermediate Calibrator (500 ng/mL)
0395	Hydrocodone 300 High Calibrator (800 ng/mL)
0397	Hydrocodone 300 Level 1 Control (225 ng/mL)
0398	Hydrocodone 300 Level 2 Control (375 ng/mL)
0400	SPICE I (JWH-018) EIA Kit, small test kit [†]
0401	SPICE I (JWH-018) EIA Kit, large test kit [†]
0002c	THC Negative Calibrator (0 ng/mL)
0402	SPICE I (JWH-018) Low Calibrator (10 ng/mL) [†]
0403	SPICE I (JWH-018) Cutoff Calibrator (20 ng/mL) [†]
0404	SPICE I (JWH-018) Intermediate Calibrator (35 ng/mL) [†]
0405	SPICE I (JWH-018) High Calibrator (50 ng/mL) [†]
0407	SPICE I (JWH-018) Level 1 Control (15 ng/mL) [†]
0408	SPICE I (JWH-018) Level 2 Control (25 ng/mL) [†]
0420	SPICE II (UR-144/XLR-11) EIA Kit, small test kit [†]
0421	SPICE II (UR-144/XLR-11) EIA Kit, large test kit [†]
0002c	THC Negative Calibrator (0 ng/mL)
0422	SPICE II (UR-144/XLR-11) Low Calibrator (10 ng/mL) [†]
0423	SPICE II (UR-144/XLR-11) Cutoff Calibrator (20 ng/mL) [†]
0424	SPICE II (UR-144/XLR-11) Intermediate Cal (35 ng/mL) [†]
0425	SPICE II (UR-144/XLR-11) High Calibrator (50 ng/mL) [†]
0427	SPICE II (UR-144/XLR-11) Level 1 Control (15 ng/mL) [†]
0428	SPICE II (UR-144/XLR-11) Level 2 Control (25 ng/mL) [†]
0450	Ethyl Glucuronide III EIA Kit, small test kit [†]
0451	Ethyl Glucuronide III EIA Kit, large test kit [†]
	EtG 200 for EtG III EIA Kit
0452	Ethyl Glucuronide 200 Low Calibrator (100 ng/mL) [†]
0453	Ethyl Glucuronide 200 Cutoff Calibrator (200 ng/mL) [†]
0454	Ethyl Glucuronide 200 Intermediate Calibrator (500 ng/mL) [†]
0455	Ethyl Glucuronide 200 High Calibrator (1000 ng/mL) [†]
0457	Ethyl Glucuronide 200 Level 1 Control (150 ng/mL) [†]
0458	Ethyl Glucuronide 200 Level 2 Control (250 ng/mL) [†]
	EtG 500 for EtG III EIA Kit
0462	Ethyl Glucuronide 500 Low Calibrator (200 ng/mL) [†]
0463	Ethyl Glucuronide 500 Cutoff Calibrator (500 ng/mL) [†]
0464	Ethyl Glucuronide 500 Intermediate Calibrator (1000 ng/mL) [†]
0465	Ethyl Glucuronide 500 High Calibrator (2000 ng/mL) [†]
0467	Ethyl Glucuronide 500 Level 1 Control (375 ng/mL) [†]
0468	Ethyl Glucuronide 500 Level 2 Control (625 ng/mL) [†]

<u>REF #</u>	<u>Product Name</u>
0520	^{510(k) cleared} Cotinine II EIA Kit, small test kit
0521	Cotinine II EIA Kit, large test kit
0522	Cotinine II Low Calibrator (100 ng/mL)
0523	Cotinine II Cutoff Calibrator (200 ng/mL)
0524	Cotinine II Intermediate Calibrator (400 ng/mL)
0525	Cotinine II High Calibrator (1000 ng/mL)
0527	Cotinine II Level 1 Control (150 ng/mL)
0528	Cotinine II Level 2 Control (250 ng/mL)
0550	Fentanyl EIA Kit (SQ), small test kit [†]
0551	Fentanyl EIA Kit (SQ), large test kit [†]
0552	Norfentanyl Low Calibrator (2.5 ng/mL) [†]
0553	Norfentanyl Cutoff Calibrator (5 ng/mL) [†]
0554	Norfentanyl Intermediate Calibrator (10 ng/mL) [†]
0555	Norfentanyl High Calibrator (20 ng/mL) [†]
0557	Norfentanyl Level 1 Control (3.75 ng/mL) [†]
0558	Norfentanyl Level 2 Control (6.25 ng/mL) [†]
0570	^{510(k) cleared} Fentanyl II EIA Kit (Q), small test kit
0571	Fentanyl II EIA Kit (Q), large test kit
0313	Norfentanyl Cutoff Calibrator (5 ng/mL)
0317	Norfentanyl Level 1 Control (3.75 ng/mL)
0318	Norfentanyl Level 2 Control (6.25 ng/mL)
0580	Fentanyl II EIA Kit (SQ), small test kit [†]
0581	Fentanyl II EIA Kit (SQ), large test kit [†]
0552	Norfentanyl Low Calibrator (2.5 ng/mL) [†]
0553	Norfentanyl Cutoff Calibrator (5 ng/mL) [†]
0554	Norfentanyl Intermediate Calibrator (10 ng/mL) [†]
0555	Norfentanyl High Calibrator (20 ng/mL) [†]
0557	Norfentanyl Level 1 Control (3.75 ng/mL) [†]
0558	Norfentanyl Level 2 Control (6.25 ng/mL) [†]
0610	^{510(k) cleared} Oxycodone III EIA Kit, small test kit
0611	Oxycodone III EIA Kit, large test kit
	OXY 100 for OXY III EIA Kit
0612	Oxycodone 100 Low Calibrator (50 ng/mL)
0613	Oxycodone 100 Cutoff Calibrator (100 ng/mL)
0614	Oxycodone 100 Intermediate Calibrator (150 ng/mL)
0615	Oxycodone 100 High Calibrator (300 ng/mL)
0617	Oxycodone 100 Level 1 Control (75 ng/mL)
0618	Oxycodone 100 Level 2 Control (125 ng/mL)
	OXY 300 for OXY III EIA Kit
0622	Oxycodone 300 Low Calibrator (150 ng/mL)
0623	Oxycodone 300 Cutoff Calibrator (300 ng/mL)
0624	Oxycodone 300 Intermediate Calibrator (500 ng/mL)
0625	Oxycodone 300 High Calibrator (800 ng/mL)
0627	Oxycodone 300 Level 1 Control (225 ng/mL)
0628	Oxycodone 300 Level 2 Control (375 ng/mL)
0630	^{510(k) cleared} Methadone II EIA Kit, small test kit
0631	Methadone II EIA Kit, large test kit
0112	Methadone Low Calibrator (150 ng/mL)
0113	Methadone Cutoff Calibrator (300 ng/mL)
0114	Methadone Intermediate Calibrator (600 ng/mL)
0115	Methadone High Calibrator (1000 ng/mL)
0117	Methadone Level 1 Control (225 ng/mL)
0118	Methadone Level 2 Control (375 ng/mL)
0852	^{510(k) cleared} Multi-Analyte Set C Low Calibrator, 15 mL [§]
0853	Multi-Analyte Set C Cutoff Calibrator, 15 mL [§]
0854	Multi-Analyte Set C Intermediate Calibrator, 15 mL [§]
0855	Multi-Analyte Set C High Calibrator, 15 mL [§]
0857	Multi-Analyte Set C Level 1 Control, 15 mL [§]
0858	Multi-Analyte Set C Level 2 Control, 15 mL [§]

Tests COC 150, MAMP, AMP 500, MTD, OPI2K, BZO 200, and BARB 200

<sup>510(k)
cleared</sup> = FDA 510(k) Cleared Assay, Calibrators, & Controls.

[†] Forensic Use Only

[‡] Employment and Insurance Testing Use Only

[§] Product is made to order. Please allow 10 to 12 weeks of production time

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂

* A large test kit is 1000/375 mL R₁/R₂



PRODUCT CATALOG: Drugs of Abuse Urinalysis by Catalog Number

<u>REF #</u>	<u>Product Name</u>
0871	Multi-Analyte Set D Negative Calibrator, 15 mL [§]
0872	Multi-Analyte Set D Low Calibrator, 15 mL [§]
0873	Multi-Analyte Set D Cutoff Calibrator, 15 mL [§]
0874	Multi-Analyte Set D Intermediate Calibrator, 15 mL [§]
0875	Multi-Analyte Set D High Calibrator, 15 mL [§]
0877	Multi-Analyte Set D Level 1 Control, 15 mL [§]
0878	Multi-Analyte Set D Level 2 Control, 15 mL [§]

Tests AMP 1000, BZO 200, COC 300, OPI 300, and OXY 100

0891	Multi-Analyte Set E Negative Calibrator, 15 mL [§]
0892	Multi-Analyte Set E Low Calibrator, 15 mL [§]
0893	Multi-Analyte Set E Cutoff Calibrator, 15 mL [§]
0894	Multi-Analyte Set E Intermediate Calibrator, 15 mL [§]
0895	Multi-Analyte Set E High Calibrator, 15 mL [§]
0897	Multi-Analyte Set E Level 1 Control, 15 mL [§]
0898	Multi-Analyte Set E Level 2 Control, 15 mL [§]

Tests BARB 200, BUP 10, MTD 300, and PCP 25

[§] = FDA 510(k) Cleared Assay, Calibrators, & Controls.

[†] Forensic Use Only

[‡] Employment and Insurance Testing Use Only

[§] Product is made to order. Please allow 10 to 12 weeks of production time

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂

* A large test kit is 1000/375 mL R₁/R₂