



PRODUCT CATALOG: Drugs of Abuse Urinalysis – EU ONLY by Alphabetical Order

<u>REF #</u>	<u>Product Name</u>	<u>REF #</u>	<u>Product Name</u>
0001	Universal Negative Calibrator, 5 mL	0270	^{510(k)} _{Cleared} Buprenorphine EIA Kit, small test kit
0290	^{510(k)} _{Cleared} 6-Acetylmorphine EIA Kit, small test kit	0271	Buprenorphine EIA Kit, large test kit
0291	6-Acetylmorphine EIA Kit, large test kit		BUP 5
0299	6-Acetylmorphine Negative Calibrator (0 ng/mL)	0272	Norbuprenorphine Level 1 Control (3 ng/mL)
0292	6-Acetylmorphine Low Calibrator (5 ng/mL)	0273	Norbuprenorphine Cutoff Calibrator #1 (5 ng/mL)
0293	6-Acetylmorphine Cutoff Calibrator (10 ng/mL)	0274	Norbuprenorphine Level 2 Control (7 ng/mL)
0294	6-Acetylmorphine Intermediate Calibrator (20 ng/mL)	0275	Norbuprenorphine Cutoff Calibrator #2 (10 ng/mL)
0295	6-Acetylmorphine High Calibrator (40 ng/mL)	0277	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)
0297	6-Acetylmorphine Level 1 Control (7.5 ng/mL)	0278	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)
0298	6-Acetylmorphine Level 2 Control (12.5 ng/mL)	0279	Norbuprenorphine High Calibrator (75 ng/mL)
			BUP 10
0100	^{510(k)} _{Cleared} Amphetamines 500 EIA Kit, small test kit	0273	Norbuprenorphine Low Calibrator (5 ng/mL)
0101	Amphetamines 500 EIA Kit, large test kit	0274	Norbuprenorphine Level 1 Control (7 ng/mL)
0102	Amphetamines 500 Low Calibrator (MAMP 250 ng/mL)	0275	Norbuprenorphine Cutoff Calibrator (10 ng/mL)
0103	Amphetamines 500 Cutoff Calibrator (MAMP 500 ng/mL)	0276	Norbuprenorphine Level 2 Control (13 ng/mL)
0104	Amphetamines 500 Intermediate Calibrator (MAMP 1000 ng/mL)	0277	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)
0105	Amphetamines 500 High Calibrator (MAMP 2000 ng/mL)	0278	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)
0107	Amphetamines 500 Level 1 Control (MAMP 375 ng/mL)	0279	Norbuprenorphine High Calibrator (75 ng/mL)
0108	Amphetamines 500 Level 2 Control (MAMP 625 ng/mL)		
0040	^{510(k)} _{Cleared} Amphetamines EIA Kit, small test kit	0070c	^{510(k)} _{Cleared} Cannabinoids (cTHC) EIA Kit, small test kit
0041	Amphetamines EIA Kit, large test kit	0071c	Cannabinoids (cTHC) EIA Kit, large test kit
0042	Amphetamines Low Calibrator (MAMP 500 ng/mL)	0002c	Cannabinoids (cTHC) Negative Calibrator (0 ng/mL)
0043	Amphetamines Cutoff Calibrator (MAMP 1000 ng/mL)		cTHC 25
0044	Amphetamines Intermediate Calibrator (MAMP 1500 ng/mL)	0072c	Cannabinoids (cTHC) Low Calibrator (12.5 ng/mL)
0045	Amphetamines High Calibrator (MAMP 2000 ng/mL)	0073c	Cannabinoids (cTHC) Cutoff Calibrator (25 ng/mL)
0047	Amphetamines Level 1 Control (MAMP 750 ng/mL)	0007c	Cannabinoids (cTHC) Intermediate Calibrator (37.5 ng/mL)
0048	Amphetamines Level 2 Control (MAMP 1250 ng/mL)	0075c	Cannabinoids (cTHC) High Calibrator (50 ng/mL)
		0006c	Cannabinoids (cTHC) Level 1 Control (18.75 ng/mL)
0140	^{510(k)} _{Cleared} Barbiturate EIA Kit, small test kit	0074c	Cannabinoids (cTHC) Level 2 Control (31.25 ng/mL)
0141	Barbiturate EIA Kit, large test kit		cTHC 50
	BARB 200	0073c	Cannabinoids (cTHC) Low Calibrator (25 ng/mL)
0142	Barbiturate Calibrator #2/Control I (Secobarbital 100 ng/mL)	0075c	Cannabinoids (cTHC) Cutoff Calibrator (50 ng/mL)
0143	Barbiturate Calibrator #3/Cutoff A (Secobarbital 200 ng/mL)	0076c	Cannabinoids (cTHC) Intermediate Calibrator (75 ng/mL)
0144	Barbiturate Calibrator #4/Control III (Secobarbital 300 ng/mL)	0077c	Cannabinoids (cTHC) High Calibrator (100 ng/mL)
0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)	0007c	Cannabinoids (cTHC) Level 1 Control (37.5 ng/mL)
	BARB 300	0008c	Cannabinoids (cTHC) Level 2 Control (62.5 ng/mL)
0142	Barbiturate Calibrator #2 (Secobarbital 100 ng/mL)		cTHC 100
0143	Barbiturate Calibrator #3/Control II (Secobarbital 200 ng/mL)	0075c	Cannabinoids (cTHC) Low Calibrator (50 ng/mL)
0144	Barbiturate Calibrator #4/Cutoff B (Secobarbital 300 ng/mL)	0077c	Cannabinoids (cTHC) Cutoff Calibrator (100 ng/mL)
0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)	0078c	Cannabinoids (cTHC) Intermediate Calibrator (150 ng/mL)
0146	Barbiturate Control IV (Secobarbital 400 ng/mL)	0079c	Cannabinoids (cTHC) High Calibrator (200 ng/mL)
		0076c	Cannabinoids (cTHC) Level 1 Control (75 ng/mL)
0130	^{510(k)} _{Cleared} Benzodiazepine EIA Kit, small test kit	0009c	Cannabinoids (cTHC) Level 2 Control (125 ng/mL)
0131	Benzodiazepine EIA Kit, large test kit		
	BZO 200	0370	^{510(k)} _{Cleared} Carisoprodol Metabolite (SOMA) EIA Kit, small test kit
0132	Benzodiazepine Calibrator #2/Control I (Oxazepam 100 ng/mL)	0371	Carisoprodol Metabolite (SOMA) EIA Kit, large test kit
0133	Benzodiazepine Calibrator #3/Cutoff A (Oxazepam 200 ng/mL)	0372	Meprobamate Low Calibrator (50 ng/mL)
0134	Benzodiazepine Calibrator #4/Control III (Oxazepam 300 ng/mL)	0373	Meprobamate Cutoff Calibrator (100 ng/mL)
0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)	0374	Meprobamate Intermediate Calibrator (200 ng/mL)
	BZO 300	0375	Meprobamate High Calibrator (400 ng/mL)
0132	Benzodiazepine Calibrator #2 (Oxazepam 100 ng/mL)	0377	Meprobamate Level 1 Control (75 ng/mL)
0133	Benzodiazepine Calibrator #3/Control II (Oxazepam 200 ng/mL)	0378	Meprobamate Level 2 Control (125 ng/mL)
0134	Benzodiazepine Calibrator #4/Cutoff B (Oxazepam 300 ng/mL)		
0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)	0340	^{510(k)} _{Cleared} Cocaine Metabolite 150 EIA Kit, small test kit
0136	Benzodiazepine Control IV (Oxazepam 400 ng/mL)	0341	Cocaine Metabolite 150 EIA Kit, large test kit
		0342	Cocaine Metabolite 150 Low Calibrator (BE 75 ng/mL)
		0343	Cocaine Metabolite 150 Cutoff Calibrator (BE 150 ng/mL)
		0344	Cocaine Metabolite 150 Intermediate Calibrator (BE 300 ng/mL)
		0345	Cocaine Metabolite 150 High Calibrator (BE 1000 ng/mL)
		0347	Cocaine Metabolite 150 Level 1 Control (BE 112.5 ng/mL)
		0348	Cocaine Metabolite 150 Level 2 Control (BE 187.5 ng/mL)

^{510(k)}_{Cleared} = U.S. FDA 510(k) cleared Assay, Calibrators, & Controls.
 § Product is made to order. Please allow 10 to 12 weeks of production time

* EIA = Enzyme Immunoassay
 * All single analyte calibrators and controls are 5 mL per bottle.
 * A small test kit is 100/37.5 mL R₁/R₂
 * A large test kit is 1000/375 mL R₁/R₂

**PRODUCT CATALOG: Drugs of Abuse Urinalysis – EU ONLY by Alphabetical Order**

REF #	Product Name	REF #	Product Name
0030	^{510(k)} _{Cleared} Cocaine Metabolite EIA Kit, small test kit	0560	Fentanyl EIA Kit, small test kit
0031	Cocaine Metabolite EIA Kit, large test kit	0561	Fentanyl EIA Kit, large test kit
0032	Cocaine Metabolite Low Calibrator (BE 150 ng/mL)	0562	Norfentanyl Low Calibrator (2.5 ng/mL)
0033	Cocaine Metabolite Cutoff Calibrator (BE 300 ng/mL)	0563	Norfentanyl Cutoff Calibrator (5 ng/mL)
0034	Cocaine Metabolite Intermediate Cal (BE 1000 ng/mL)	0564	Norfentanyl Intermediate Calibrator (10 ng/mL)
0035	Cocaine Metabolite High Calibrator (BE 3000 ng/mL)	0565	Norfentanyl High Calibrator (20 ng/mL)
0037	Cocaine Metabolite Level 1 Control (BE 225 ng/mL)	0567	Norfentanyl Level 1 Control (3.75 ng/mL)
0038	Cocaine Metabolite Level 2 Control (BE 375 ng/mL)	0568	Norfentanyl Level 2 Control (6.25 ng/mL)
0480	Cotinine EIA Kit, small test kit	0590	Fentanyl II EIA Kit, small test kit
0481	Cotinine EIA Kit, large test kit	0591	Fentanyl II EIA Kit, large test kit
0482	Cotinine Low Calibrator/Level 1 Control (250 ng/mL)	0562	Norfentanyl Low Calibrator (2.5 ng/mL)
0483	Cotinine Cutoff Calibrator (500 ng/mL)	0563	Norfentanyl Cutoff Calibrator (5 ng/mL)
0484	Cotinine Intermediate Calibrator (1000 ng/mL)	0564	Norfentanyl Intermediate Calibrator (10 ng/mL)
0485	Cotinine High Calibrator (2000 ng/mL)	0565	Norfentanyl High Calibrator (20 ng/mL)
0488	Cotinine Level 2 Control (750 ng/mL)	0567	Norfentanyl Level 1 Control (3.75 ng/mL)
0520	^{510(k)} _{Cleared} Cotinine II EIA Kit, small test kit	0568	Norfentanyl Level 2 Control (6.25 ng/mL)
0521	Cotinine II EIA Kit, large test kit	0380	^{510(k)} _{Cleared} Hydrocodone 100 EIA Kit, small test kit
0522	Cotinine II Low Calibrator (100 ng/mL)	0381	Hydrocodone 100 EIA Kit, large test kit
0523	Cotinine II Cutoff Calibrator (200 ng/mL)	0382	Hydrocodone 100 Low Calibrator (50 ng/mL)
0524	Cotinine II Intermediate Calibrator (400 ng/mL)	0383	Hydrocodone 100 Cutoff Calibrator (100 ng/mL)
0525	Cotinine II High Calibrator (1000 ng/mL)	0384	Hydrocodone 100 Intermediate Calibrator (150 ng/mL)
0527	Cotinine II Level 1 Control (150 ng/mL)	0385	Hydrocodone 100 High Calibrator (300 ng/mL)
0528	Cotinine II Level 2 Control (250 ng/mL)	0387	Hydrocodone 100 Level 1 Control (75 ng/mL)
0160	^{510(k)} _{Cleared} Ecstasy (MDMA) EIA Kit, small test kit	0388	Hydrocodone 100 Level 2 Control (125 ng/mL)
0161	Ecstasy (MDMA) EIA Kit, large test kit	0390	^{510(k)} _{Cleared} Hydrocodone 300 EIA Kit, small test kit
0162	Ecstasy Low Calibrator (100 ng/mL)	0391	Hydrocodone 300 EIA Kit, large test kit
0163	Ecstasy Cutoff Calibrator (500 ng/mL)	0392	Hydrocodone 300 Low Calibrator (150 ng/mL)
0164	Ecstasy Intermediate Calibrator (750 ng/mL)	0393	Hydrocodone 300 Cutoff Calibrator (300 ng/mL)
0165	Ecstasy High Calibrator (1000 ng/mL)	0394	Hydrocodone 300 Intermediate Calibrator (500 ng/mL)
0167	Ecstasy Level 1 Control (375 ng/mL)	0395	Hydrocodone 300 High Calibrator (800 ng/mL)
0168	Ecstasy Level 2 Control (625 ng/mL)	0397	Hydrocodone 300 Level 1 Control (225 ng/mL)
0220	^{510(k)} _{Cleared} Ethyl Alcohol Enzymatic Assay Kit, small test kit	0398	Hydrocodone 300 Level 2 Control (375 ng/mL)
0221	Ethyl Alcohol Enzymatic Assay Kit, large test kit	0110	^{510(k)} _{Cleared} Methadone EIA Kit, small test kit
0223	Ethyl Alcohol Calibrator (100 mg/dL)	0111	Methadone EIA Kit, large test kit
0224	Ethyl Alcohol Level 1 Control (50 mg/dL)	0112	Methadone Low Calibrator (150 ng/mL)
0225	Ethyl Alcohol Level 2 Control (300 mg/dL)	0113	Methadone Cutoff Calibrator (300 ng/mL)
0490b	Ethyl Glucuronide EIA Kit, small test kit	0114	Methadone Intermediate Calibrator (600 ng/mL)
0491b	Ethyl Glucuronide EIA Kit, large test kit	0115	Methadone High Calibrator (1000 ng/mL)
0492b	Ethyl Glucuronide Low Calibrator (250 ng/mL)	0117	Methadone Level 1 Control (225 ng/mL)
0493b	Ethyl Glucuronide Cutoff Calibrator (1000 ng/mL)	0118	Methadone Level 2 Control (375 ng/mL)
0494b	Ethyl Glucuronide Intermediate Calibrator #1 (2000 ng/mL)	0630	^{510(k)} _{Cleared} Methadone II EIA Kit, small test kit
0495b	Ethyl Glucuronide Intermediate Calibrator #2 (3000 ng/mL)	0631	Methadone II EIA Kit, large test kit
0496b	Ethyl Glucuronide High Calibrator (4000 ng/mL)	0112	Methadone Low Calibrator (150 ng/mL)
0497b	Ethyl Glucuronide Level 1 Control (500 ng/mL)	0113	Methadone Cutoff Calibrator (300 ng/mL)
0498b	Ethyl Glucuronide Level 2 Control (1500 ng/mL)	0114	Methadone Intermediate Calibrator (600 ng/mL)
0530	Ethyl Glucuronide III EIA Kit, small test kit	0115	Methadone High Calibrator (1000 ng/mL)
0531	Ethyl Glucuronide III EIA Kit, large test kit	0117	Methadone Level 1 Control (225 ng/mL)
	EtG 200 for EtG III EIA Kit	0118	Methadone Level 2 Control (375 ng/mL)
0532	Ethyl Glucuronide 200 Low Calibrator (100 ng/mL)	0190b	^{510(k)} _{Cleared} Methadone Metabolite (EDDP) EIA Kit, small test kit
0533	Ethyl Glucuronide 200 Cutoff Calibrator (200 ng/mL)	0191b	Methadone Metabolite (EDDP) EIA Kit, large test kit
0534	Ethyl Glucuronide 200 Intermediate Calibrator (500 ng/mL)		EDDPb 100
0535	Ethyl Glucuronide 200 High Calibrator (1000 ng/mL)	0202b	Methadone Metabolite 100 Low Calibrator (50 ng/mL)
0537	Ethyl Glucuronide 200 Level 1 Control (150 ng/mL)	0203b	Methadone Metabolite 100 Cutoff Calibrator (100 ng/mL)
0538	Ethyl Glucuronide 200 Level 2 Control (250 ng/mL)	0204b	Methadone Metabolite 100 Intermediate Calibrator (250 ng/mL)
	EtG 500 for EtG III EIA Kit	0205b	Methadone Metabolite 100 High Calibrator (500 ng/mL)
0542	Ethyl Glucuronide 500 Low Calibrator (200 ng/mL)	0207b	Methadone Metabolite 100 Level 1 Control (75 ng/mL)
0543	Ethyl Glucuronide 500 Cutoff Calibrator (500 ng/mL)	0208b	Methadone Metabolite 100 Level 2 Control (125 ng/mL)
0544	Ethyl Glucuronide 500 Intermediate Calibrator (1000 ng/mL)		
0545	Ethyl Glucuronide 500 High Calibrator (2000 ng/mL)		
0547	Ethyl Glucuronide 500 Level 1 Control (375 ng/mL)		
0548	Ethyl Glucuronide 500 Level 2 Control (625 ng/mL)		

^{510(k)}_{Cleared} = U.S. FDA 510(k) cleared Assay, Calibrators, & Controls.
 § Product is made to order. Please allow 10 to 12 weeks of production time

* EIA = Enzyme Immunoassay
 * All single analyte calibrators and controls are 5 mL per bottle.
 * A small test kit is 100/37.5 mL R₁/R₂
 * A large test kit is 1000/375 mL R₁/R₂



PRODUCT CATALOG: Drugs of Abuse Urinalysis – EU ONLY by Alphabetical Order

<u>REF #</u>	<u>Product Name</u>
EDDPb 300	
0192b	Methadone Metabolite 300 Low Calibrator (150 ng/mL)
0193b	Methadone Metabolite 300 Cutoff Calibrator (300 ng/mL)
0194b	Methadone Metabolite 300 Intermediate Calibrator (600 ng/mL)
0195b	Methadone Metabolite 300 High Calibrator (1000 ng/mL)
0197b	Methadone Metabolite 300 Level 1 Control (225 ng/mL)
0198b	Methadone Metabolite 300 Level 2 Control (375 ng/mL)
0350	^{510(k)} _{Cleared} Methamphetamine EIA Kit, small test kit
0351	Methamphetamine EIA Kit, large test kit
0352	Methamphetamine Low Calibrator (250 ng/mL)
0353	Methamphetamine Cutoff Calibrator (500 ng/mL)
0354	Methamphetamine Intermediate Calibrator (1000 ng/mL)
0355	Methamphetamine High Calibrator (2000 ng/mL)
0357	Methamphetamine Level 1 Control (375 ng/mL)
0358	Methamphetamine Level 2 Control (625 ng/mL)
0852	^{510(k)} _{Cleared} Multi-Analyte Set C Low Calibrator, 15 mL [§]
0853	Multi-Analyte Set C Cutoff Calibrator, 15 mL [§]
0854	Multi-Analyte Set C Intermediate Calibrator, 15 mL [§]
0855	Multi-Analyte Set C High Calibrator, 15 mL [§]
0857	Multi-Analyte Set C Level 1 Control, 15 mL [§]
0858	Multi-Analyte Set C Level 2 Control, 15 mL [§]
<i>Tests COC 150, MAMP, AMP 500, MTD, OPI2K, BZO 200, and BARB 200</i>	
0871	Multi-Analyte Set D Negative Calibrator, 15 mL [§]
0872	Multi-Analyte Set D Low Calibrator, 15 mL [§]
0873	Multi-Analyte Set D Cutoff Calibrator, 15 mL [§]
0874	Multi-Analyte Set D Intermediate Calibrator, 15 mL [§]
0875	Multi-Analyte Set D High Calibrator, 15 mL [§]
0877	Multi-Analyte Set D Level 1 Control, 15 mL [§]
0878	Multi-Analyte Set D Level 2 Control, 15 mL [§]
<i>Tests AMP 1000, BZO 200, COC 300, OPI 300, and OXY 100</i>	
0891	Multi-Analyte Set E Negative Calibrator, 15 mL [§]
0892	Multi-Analyte Set E Low Calibrator, 15 mL [§]
0893	Multi-Analyte Set E Cutoff Calibrator, 15 mL [§]
0894	Multi-Analyte Set E Intermediate Calibrator, 15 mL [§]
0895	Multi-Analyte Set E High Calibrator, 15 mL [§]
0897	Multi-Analyte Set E Level 1 Control, 15 mL [§]
0898	Multi-Analyte Set E Level 2 Control, 15 mL [§]
<i>Tests BARB 200, BUP 10, MTD 300, and PCP 25</i>	
0020b	^{510(k)} _{Cleared} Opiate EIA Kit, small test kit
0021b	Opiate EIA Kit, large test kit
0022	Opiate Low Calibrator (Morphine 150 ng/mL)
0023	Opiate Cutoff Calibrator (Morphine 300 ng/mL)
0024	Opiate Intermediate Calibrator (Morphine 600 ng/mL)
0025	Opiate High Calibrator (Morphine 1000 ng/mL)
0027	Opiate Level 1 Control (Morphine 225 ng/mL)
0028	Opiate Level 2 Control (Morphine 375 ng/mL)
0330	^{510(k)} _{Cleared} Opiate 2000 EIA Kit, small test kit
0331	Opiate 2000 EIA Kit, large test kit
0332	Opiate 2000 Low Calibrator (Morphine 1000 ng/mL)
0333	Opiate 2000 Cutoff Calibrator (Morphine 2000 ng/mL)
0334	Opiate 2000 Intermediate Calibrator (Morphine 4000 ng/mL)
0335	Opiate 2000 High Calibrator (Morphine 6000 ng/mL)
0337	Opiate 2000 Level 1 Control (Morphine 1500 ng/mL)
0338	Opiate 2000 Level 2 Control (Morphine 2500 ng/mL)

<u>REF #</u>	<u>Product Name</u>
0240b	^{510(k)} _{Cleared} Oxycodone EIA Kit, small test kit
0241b	Oxycodone EIA Kit, large test kit
OXY 100 "b"	
0252b	Oxycodone Low Calibrator (50 ng/mL)
0242b	Oxycodone Level 1 Control (75 ng/mL)
0243b	Oxycodone Cutoff Calibrator #1 (100 ng/mL)
0244b	Oxycodone Level 2 Control (125 ng/mL)
0246b	Oxycodone Cutoff Calibrator #2 (300 ng/mL)
0248b	Oxycodone Intermediate Calibrator (500 ng/mL)
0249b	Oxycodone High Calibrator (800 ng/mL)
OXY 300 "b"	
0252b	Oxycodone Low Calibrator (50 ng/mL)
0243b	Oxycodone Cutoff Calibrator #1 (100 ng/mL)
0245b	Oxycodone Level 1 Control (225 ng/mL)
0246b	Oxycodone Cutoff Calibrator #2 (300 ng/mL)
0247b	Oxycodone Level 2 Control (375 ng/mL)
0248b	Oxycodone Intermediate Calibrator (500 ng/mL)
0249b	Oxycodone High Calibrator (800 ng/mL)
0610	^{510(k)} _{Cleared} Oxycodone III EIA Kit, small test kit
0611	Oxycodone III EIA Kit, large test kit
OXY 100 for OXY III EIA Kit	
0612	Oxycodone 100 Low Calibrator (50 ng/mL)
0613	Oxycodone 100 Cutoff Calibrator (100 ng/mL)
0614	Oxycodone 100 Intermediate Calibrator (150 ng/mL)
0615	Oxycodone 100 High Calibrator (300 ng/mL)
0617	Oxycodone 100 Level 1 Control (75 ng/mL)
0618	Oxycodone 100 Level 2 Control (125 ng/mL)
OXY 300 for OXY III EIA Kit	
0622	Oxycodone 300 Low Calibrator (150 ng/mL)
0623	Oxycodone 300 Cutoff Calibrator (300 ng/mL)
0624	Oxycodone 300 Intermediate Calibrator (500 ng/mL)
0625	Oxycodone 300 High Calibrator (800 ng/mL)
0627	Oxycodone 300 Level 1 Control (225 ng/mL)
0628	Oxycodone 300 Level 2 Control (375 ng/mL)
0010	^{510(k)} _{Cleared} Phencyclidine EIA Kit, small test kit
0011	Phencyclidine EIA Kit, large test kit
0012	Phencyclidine Low Calibrator (12.5 ng/mL)
0013	Phencyclidine Cutoff Calibrator (25 ng/mL)
0014	Phencyclidine Intermediate Calibrator (50 ng/mL)
0015	Phencyclidine High Calibrator (100 ng/mL)
0017	Phencyclidine Level 1 Control (18 ng/mL)
0018	Phencyclidine Level 2 Control (32 ng/mL)
0120	^{510(k)} _{Cleared} Propoxyphene EIA Kit, small test kit
0121	Propoxyphene EIA Kit, large test kit
0122	Propoxyphene Low Calibrator (150 ng/mL)
0123	Propoxyphene Cutoff Calibrator (300 ng/mL)
0124	Propoxyphene Intermediate Calibrator (600 ng/mL)
0125	Propoxyphene High Calibrator (1000 ng/mL)
0127	Propoxyphene Level 1 Control (225 ng/mL)
0128	Propoxyphene Level 2 Control (375 ng/mL)
0500	SPICE I (JWH-018) EIA Kit, small test kit
0501	SPICE I (JWH-018) EIA Kit, large test kit
0002c	THC Negative Calibrator (0 ng/mL)
0502	SPICE I (JWH-018) Low Calibrator (10 ng/mL)
0503	SPICE I (JWH-018) Cutoff Calibrator (20 ng/mL)
0504	SPICE I (JWH-018) Intermediate Calibrator (35 ng/mL)
0505	SPICE I (JWH-018) High Calibrator (50 ng/mL)
0507	SPICE I (JWH-018) Level 1 Control (15 ng/mL)
0508	SPICE I (JWH-018) Level 2 Control (25 ng/mL)

^{510(k)}_{Cleared} = U.S. FDA 510(k) cleared Assay, Calibrators, & Controls.
[§] Product is made to order. Please allow 10 to 12 weeks of production time

* EIA = Enzyme Immunoassay
 * All single analyte calibrators and controls are 5 mL per bottle.
 * A small test kit is 100/37.5 mL R₁/R₂
 * A large test kit is 1000/375 mL R₁/R₂



PRODUCT CATALOG: Drugs of Abuse Urinalysis – EU ONLY by Alphabetical Order

REF #	Product Name
0510	SPICE II (UR-144/XLR-11) EIA Kit, small test kit
0511	SPICE II (UR-144/XLR-11) EIA Kit, large test kit
0002c	THC Negative Calibrator (0 ng/mL)
0512	SPICE II (UR-144/XLR-11) Low Calibrator (10 ng/mL)
0513	SPICE II (UR-144/XLR-11) Cutoff Calibrator (20 ng/mL)
0514	SPICE II (UR-144/XLR-11) Intermediate Cal (35 ng/mL)
0515	SPICE II (UR-144/XLR-11) High Calibrator (50 ng/mL)
0517	SPICE II (UR-144/XLR-11) Level 1 Control (15 ng/mL)
0518	SPICE II (UR-144/XLR-11) Level 2 Control (25 ng/mL)
0410	^{510(k)} Tramadol EIA Kit, small test kit
0411	Tramadol EIA Kit, large test kit
0412	Tramadol Low Calibrator (50 ng/mL)
0413	Tramadol Cutoff Calibrator (100 ng/mL)
0414	Tramadol Intermediate Calibrator (225 ng/mL)
0415	Tramadol High Calibrator (400 ng/mL)
0417	Tramadol Level 1 Control (75 ng/mL)
0418	Tramadol Level 2 Control (125 ng/mL)

^{510(k)} = U.S. FDA 510(k) cleared Assay, Calibrators, & Controls.
§ Product is made to order. Please allow 10 to 12 weeks of production time

* EIA = Enzyme Immunoassay
* All single analyte calibrators and controls are 5 mL per bottle.
* A small test kit is 100/37.5 mL R₁/R₂
* A large test kit is 1000/375 mL R₁/R₂