

## Multi-Analyte Calibrators/Controls Summary

## Multi-Analyte Calibrators and Controls Set C

<i>Works with 7 EIAs</i>	<b>#0001</b> Universal Neg Calibrator	<b>#0852</b> Low Calibrator	<b>#0857</b> Level I Control	<b>#0853</b> Cutoff Calibrator	<b>#0858</b> Level II Control	<b>#0854</b> Intermediate Cal	<b>#0855</b> High Calibrator
<b>Benzoyllecgonine</b> (Cocaine 150 EIA)	0 ng/mL	75 ng/mL	112.5 ng/mL	150 ng/mL	187.5 ng/mL	300 ng/mL	1000 ng/mL
<b>d -Methamphetamine</b> (Amphetamines 500 EIA & Methamphetamine EIA)	0 ng/mL	250 ng/mL	375 ng/mL	500 ng/mL	625 ng/mL	750 ng/mL	1000 ng/mL
<b>Methadone</b> (Methadone EIA)	0 ng/mL	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	600 ng/mL	1000 ng/mL
<b>Morphine</b> (Opiate 2000 EIA)	0 ng/mL	1000 ng/mL	1500 ng/mL	2000 ng/mL	2500 ng/mL	4000 ng/mL	6000 ng/mL
<b>Oxazepam</b> (Benzodiazepine 200 EIA)	0 ng/mL	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL
<b>Secobarbital</b> (Barbiturate 200 EIA)	0 ng/mL	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL

## Multi-Analyte Calibrators and Controls Set D

<i>Works with 6 EIAs</i>	<b>#0871</b> Neg Calibrator	<b>#0872</b> Low Calibrator	<b>#0877</b> Level I Control	<b>#0873</b> Cutoff Calibrator	<b>#0878</b> Level II Control	<b>#0874</b> Intermediate Cal	<b>#0875</b> High Calibrator
<b>Benzoyllecgonine</b> (Cocaine 150 EIA)	0 ng/mL	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	1000 ng/mL	3000 ng/mL
<b>d -Methamphetamine</b> (Amphetamines 1000 EIA)	0 ng/mL	500 ng/mL	750 ng/mL	1000 ng/mL	1250 ng/mL	1500 ng/mL	2000 ng/mL
<b>Morphine</b> (Opiate 300 EIA)	0 ng/mL	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	600 ng/mL	1000 ng/mL
<b>Oxazepam</b> (Benzodiazepine 200 EIA)	0 ng/mL	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL
<b>Oxycodone*</b> (Oxycodone 100 EIA)	0 ng/mL	50 ng/mL	75 ng/mL	100 ng/mL	125 ng/mL	300 ng/mL	500 ng/mL

\* Semi-quantitative calibration of the Oxycodone EIA 100 ng/mL cutoff requires additional purchase of 800 ng/mL high calibrator (Ref#0249b)

## Multi-Analyte Calibrators and Controls Set E

<i>Works with 4 EIAs</i>	<b>#0891</b> Neg Calibrator	<b>#0892</b> Low Calibrator	<b>#0897</b> Level I Control	<b>#0893</b> Cutoff Calibrator	<b>#0898</b> Level II Control	<b>#0894</b> Intermediate Cal	<b>#0895</b> High Calibrator
<b>Secobarbital</b> (Barbiturate 200 EIA)	0 ng/mL	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL
<b>Norbuprenorphine**</b> (Buprenorphine 10 EIA)	0 ng/mL	5 ng/mL	7 ng/mL	10 ng/mL	13 ng/mL	20 ng/mL	40 ng/mL
<b>Methadone</b> (Methadone 300 EIA)	0 ng/mL	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	600 ng/mL	1000 ng/mL
<b>Phencyclidine</b> (Phencyclidine 25 EIA)	0 ng/mL	12.5 ng/mL	18 ng/mL	25 ng/mL	32 ng/mL	50 ng/mL	100 ng/mL

\*\* Semi-quantitative calibration of the Buprenorphine EIA 10 ng/mL cutoff requires additional purchase of 75 ng/mL high calibrator (Ref# 0279)