



## PRODUCT CATALOG: Drugs of Abuse Urinalysis by Catalog Number

<u>REF #</u>	<u>Product Name</u>	<u>REF #</u>	<u>Product Name</u>
0001	Universal Negative Calibrator, 5 mL	0100	<sup>510(k) Cleared</sup> Amphetamines 500 EIA Kit, small test kit
0010	<sup>510(k) Cleared</sup> Phencyclidine EIA Kit, small test kit	0101	Amphetamines 500 EIA Kit, large test kit
0011	Phencyclidine EIA Kit, large test kit	0102	Amphetamines 500 Low Calibrator (MAMP 250 ng/mL)
0012	Phencyclidine Low Calibrator (12.5 ng/mL)	0103	Amphetamines 500 Cutoff Calibrator (MAMP 500 ng/mL)
0013	Phencyclidine Cutoff Calibrator (25 ng/mL)	0104	Amphetamines 500 Intermediate Calibrator (MAMP 1000 ng/mL)
0014	Phencyclidine Intermediate Calibrator (50 ng/mL)	0105	Amphetamines 500 High Calibrator (MAMP 2000 ng/mL)
0015	Phencyclidine High Calibrator (100 ng/mL)	0107	Amphetamines 500 Level 1 Control (MAMP 375 ng/mL)
0017	Phencyclidine Level 1 Control (18 ng/mL)	0108	Amphetamines 500 Level 2 Control (MAMP 625 ng/mL)
0018	Phencyclidine Level 2 Control (32 ng/mL)	0120	<sup>510(k) Cleared</sup> Propoxyphene EIA Kit, small test kit
0020b	<sup>510(k) Cleared</sup> Opiate EIA Kit, small test kit	0121	Propoxyphene EIA Kit, large test kit
0021b	Opiate EIA Kit, large test kit	0122	Propoxyphene Low Calibrator (150 ng/mL)
0022	Opiate Low Calibrator (Morphine 150 ng/mL)	0123	Propoxyphene Cutoff Calibrator (300 ng/mL)
0023	Opiate Cutoff Calibrator (Morphine 300 ng/mL)	0124	Propoxyphene Intermediate Calibrator (600 ng/mL)
0024	Opiate Intermediate Calibrator (Morphine 600 ng/mL)	0125	Propoxyphene High Calibrator (1000 ng/mL)
0025	Opiate High Calibrator (Morphine 1000 ng/mL)	0127	Propoxyphene Level 1 Control (225 ng/mL)
0027	Opiate Level 1 Control (Morphine 225 ng/mL)	0128	Propoxyphene Level 2 Control (375 ng/mL)
0028	Opiate Level 2 Control (Morphine 375 ng/mL)	0130	<sup>510(k) Cleared</sup> Benzodiazepine EIA Kit, small test kit
0030	<sup>510(k) Cleared</sup> Cocaine Metabolite EIA Kit, small test kit	0131	Benzodiazepine EIA Kit, large test kit
0031	Cocaine Metabolite EIA Kit, large test kit		<b>BZO 200</b>
0032	Cocaine Metabolite Low Calibrator (BE 150 ng/mL)	0132	Benzodiazepine Calibrator #2/Control I (Oxazepam 100 ng/mL)
0033	Cocaine Metabolite Cutoff Calibrator (BE 300 ng/mL)	0133	Benzodiazepine Calibrator #3/Cutoff A (Oxazepam 200 ng/mL)
0034	Cocaine Metabolite Intermediate Calibrator (BE 1000 ng/mL)	0134	Benzodiazepine Calibrator #4/Control III (Oxazepam 300 ng/mL)
0035	Cocaine Metabolite High Calibrator (BE 3000 ng/mL)	0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)
0037	Cocaine Metabolite Level 1 Control (BE 225 ng/mL)		<b>BZO 300</b>
0038	Cocaine Metabolite Level 2 Control (BE 375 ng/mL)	0132	Benzodiazepine Calibrator #2 (Oxazepam 100 ng/mL)
0040	<sup>510(k) Cleared</sup> Amphetamines EIA Kit, small test kit	0133	Benzodiazepine Calibrator #3/Control II (Oxazepam 200 ng/mL)
0041	Amphetamines EIA Kit, large test kit	0134	Benzodiazepine Calibrator #4/Cutoff B (Oxazepam 300 ng/mL)
0042	Amphetamines Low Calibrator (MAMP 500 ng/mL)	0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)
0043	Amphetamines Cutoff Calibrator (MAMP 1000 ng/mL)	0136	Benzodiazepine Control IV (Oxazepam 400 ng/mL)
0044	Amphetamines Intermediate Calibrator (MAMP 1500 ng/mL)	0140	<sup>510(k) Cleared</sup> Barbiturate EIA Kit, small test kit
0045	Amphetamines High Calibrator (MAMP 2000 ng/mL)	0141	Barbiturate EIA Kit, large test kit
0047	Amphetamines Level 1 Control (MAMP 750 ng/mL)		<b>BARB 200</b>
0048	Amphetamines Level 2 Control (MAMP 1250 ng/mL)	0142	Barbiturate Calibrator #2/Control I (Secobarbital 100 ng/mL)
0070c	<sup>510(k) Cleared</sup> Cannabinoids (cTHC) EIA Kit, small test kit	0143	Barbiturate Calibrator #3/Cutoff A (Secobarbital 200 ng/mL)
0071c	Cannabinoids (cTHC) EIA Kit, large test kit	0144	Barbiturate Calibrator #4/Control III (Secobarbital 300 ng/mL)
0002c	Cannabinoids (cTHC) Negative Calibrator (0 ng/mL)	0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)
	<b>cTHC 25</b>		<b>BARB 300</b>
0072c	Cannabinoids (cTHC) Low Calibrator (12.5 ng/mL)	0142	Barbiturate Calibrator #2 (Secobarbital 100 ng/mL)
0073c	Cannabinoids (cTHC) Cutoff Calibrator (25 ng/mL)	0143	Barbiturate Calibrator #3/Control II (Secobarbital 200 ng/mL)
0007c	Cannabinoids (cTHC) Intermediate Calibrator (37.5 ng/mL)	0144	Barbiturate Calibrator #4/Cutoff B (Secobarbital 300 ng/mL)
0075c	Cannabinoids (cTHC) High Calibrator (50 ng/mL)	0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)
0006c	Cannabinoids (cTHC) Level 1 Control (18.75 ng/mL)	0146	Barbiturate Control IV (Secobarbital 400 ng/mL)
0074c	Cannabinoids (cTHC) Level 2 Control (31.25 ng/mL)	0160	<sup>510(k) Cleared</sup> Ecstasy (MDMA) EIA Kit, small test kit
	<b>cTHC 50</b>	0161	Ecstasy (MDMA) EIA Kit, large test kit
0073c	Cannabinoids (cTHC) Low Calibrator (25 ng/mL)	0162	Ecstasy Low Calibrator (100 ng/mL)
0075c	Cannabinoids (cTHC) Cutoff Calibrator (50 ng/mL)	0163	Ecstasy Cutoff Calibrator (500 ng/mL)
0076c	Cannabinoids (cTHC) Intermediate Calibrator (75 ng/mL)	0164	Ecstasy Intermediate Calibrator (750 ng/mL)
0077c	Cannabinoids (cTHC) High Calibrator (100 ng/mL)	0165	Ecstasy High Calibrator (1000 ng/mL)
0007c	Cannabinoids (cTHC) Level 1 Control (37.5 ng/mL)	0167	Ecstasy Level 1 Control (375 ng/mL)
0008c	Cannabinoids (cTHC) Level 2 Control (62.5 ng/mL)	0168	Ecstasy Level 2 Control (625 ng/mL)
	<b>cTHC 100</b>		
0075c	Cannabinoids (cTHC) Low Calibrator (50 ng/mL)		
0077c	Cannabinoids (cTHC) Cutoff Calibrator (100 ng/mL)		
0078c	Cannabinoids (cTHC) Intermediate Calibrator (150 ng/mL)		
0079c	Cannabinoids (cTHC) High Calibrator (200 ng/mL)		
0076c	Cannabinoids (cTHC) Level 1 Control (75 ng/mL)		
0009c	Cannabinoids (cTHC) Level 2 Control (125 ng/mL)		

<sup>510(k) Cleared</sup> = FDA 510(k) Cleared Assay, Calibrators, & Controls.

† Forensic Use Only

§ Product is made to order. Please allow 10 to 12 weeks of production time. Minimum purchase volume required.

\* EIA = Enzyme Immunoassay

\* All single analyte calibrators and controls are 5 mL per bottle.

\* A small test kit is 100/37.5 mL R<sub>1</sub>/R<sub>2</sub>

\* A large test kit is 1000/375 mL R<sub>1</sub>/R<sub>2</sub>



## PRODUCT CATALOG: Drugs of Abuse Urinalysis by Catalog Number

<u>REF #</u>	<u>Product Name</u>	<u>REF #</u>	<u>Product Name</u>
<b>0190b</b>	<small>510(k) Cleared</small> Methadone Metabolite (EDDP) EIA Kit, small test kit	<b>0340</b>	<small>510(k) Cleared</small> Cocaine Metabolite 150 EIA Kit, small test kit
<b>0191b</b>	Methadone Metabolite (EDDP) EIA Kit, large test kit	<b>0341</b>	Cocaine Metabolite 150 EIA Kit, large test kit
	<b>EDDPb 100</b>	<b>0342</b>	Cocaine Metabolite 150 Low Calibrator (BE 75 ng/mL)
<b>0202b</b>	Methadone Metabolite 100 Low Calibrator (50 ng/mL)	<b>0343</b>	Cocaine Metabolite 150 Cutoff Calibrator (BE 150 ng/mL)
<b>0203b</b>	Methadone Metabolite 100 Cutoff Calibrator (100 ng/mL)	<b>0344</b>	Cocaine Metabolite 150 Intermediate Calibrator (BE 300 ng/mL)
<b>0204b</b>	Methadone Metabolite 100 Intermediate Calibrator (250 ng/mL)	<b>0345</b>	Cocaine Metabolite 150 High Calibrator (BE 1000 ng/mL)
<b>0205b</b>	Methadone Metabolite 100 High Calibrator (500 ng/mL)	<b>0347</b>	Cocaine Metabolite 150 Level 1 Control (BE 112.5 ng/mL)
<b>0207b</b>	Methadone Metabolite 100 Level 1 Control (75 ng/mL)	<b>0348</b>	Cocaine Metabolite 150 Level 2 Control (BE 187.5 ng/mL)
<b>0208b</b>	Methadone Metabolite 100 Level 2 Control (125 ng/mL)		
	<b>EDDPb 300</b>	<b>0350</b>	<small>510(k) Cleared</small> Methamphetamine EIA Kit, small test kit
<b>0192b</b>	Methadone Metabolite 300 Low Calibrator (150 ng/mL)	<b>0351</b>	Methamphetamine EIA Kit, large test kit
<b>0193b</b>	Methadone Metabolite 300 Cutoff Calibrator (300 ng/mL)	<b>0352</b>	Methamphetamine Low Calibrator (250 ng/mL)
<b>0194b</b>	Methadone Metabolite 300 Intermediate Calibrator (600 ng/mL)	<b>0353</b>	Methamphetamine Cutoff Calibrator (500 ng/mL)
<b>0195b</b>	Methadone Metabolite 300 High Calibrator (1000 ng/mL)	<b>0354</b>	Methamphetamine Intermediate Calibrator (1000 ng/mL)
<b>0197b</b>	Methadone Metabolite 300 Level 1 Control (225 ng/mL)	<b>0355</b>	Methamphetamine High Calibrator (2000 ng/mL)
<b>0198b</b>	Methadone Metabolite 300 Level 2 Control (375 ng/mL)	<b>0357</b>	Methamphetamine Level 1 Control (375 ng/mL)
		<b>0358</b>	Methamphetamine Level 2 Control (625 ng/mL)
<b>0220</b>	<small>510(k) Cleared</small> Ethyl Alcohol Enzymatic Assay Kit, small test kit	<b>0370</b>	<small>510(k) Cleared</small> Carisoprodol Metabolite (SOMA) EIA Kit, small test kit <sup>§</sup>
<b>0221</b>	Ethyl Alcohol Enzymatic Assay Kit, large test kit	<b>0371</b>	Carisoprodol Metabolite (SOMA) EIA Kit, large test kit <sup>§</sup>
<b>0223</b>	Ethyl Alcohol Calibrator (100 mg/dL)	<b>0372</b>	Meprobamate Low Calibrator (50 ng/mL) <sup>§</sup>
<b>0224</b>	Ethyl Alcohol Level 1 Control (50 mg/dL)	<b>0373</b>	Meprobamate Cutoff Calibrator (100 ng/mL) <sup>§</sup>
<b>0225</b>	Ethyl Alcohol Level 2 Control (300 mg/dL)	<b>0374</b>	Meprobamate Intermediate Calibrator (200 ng/mL) <sup>§</sup>
<b>0270</b>	<small>510(k) Cleared</small> Buprenorphine EIA Kit, small test kit	<b>0375</b>	Meprobamate High Calibrator (400 ng/mL) <sup>§</sup>
<b>0271</b>	Buprenorphine EIA Kit, large test kit	<b>0377</b>	Meprobamate Level 1 Control (75 ng/mL) <sup>§</sup>
	<b>BUP 5</b>	<b>0378</b>	Meprobamate Level 2 Control (125 ng/mL) <sup>§</sup>
<b>0272</b>	Norbuprenorphine Level 1 Control (3 ng/mL)		
<b>0273</b>	Norbuprenorphine Cutoff Calibrator #1 (5 ng/mL)	<b>0380</b>	<small>510(k) Cleared</small> Hydrocodone 100 EIA Kit, small test kit
<b>0274</b>	Norbuprenorphine Level 2 Control (7 ng/mL)	<b>0381</b>	Hydrocodone 100 EIA Kit, large test kit
<b>0275</b>	Norbuprenorphine Cutoff Calibrator #2 (10 ng/mL)	<b>0382</b>	Hydrocodone 100 Low Calibrator (50 ng/mL)
<b>0277</b>	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)	<b>0383</b>	Hydrocodone 100 Cutoff Calibrator (100 ng/mL)
<b>0278</b>	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)	<b>0384</b>	Hydrocodone 100 Intermediate Calibrator (150 ng/mL)
<b>0279</b>	Norbuprenorphine High Calibrator (75 ng/mL)	<b>0385</b>	Hydrocodone 100 High Calibrator (300 ng/mL)
	<b>BUP 10</b>	<b>0387</b>	Hydrocodone 100 Level 1 Control (75 ng/mL)
<b>0273</b>	Norbuprenorphine Low Calibrator (5 ng/mL)	<b>0388</b>	Hydrocodone 100 Level 2 Control (125 ng/mL)
<b>0274</b>	Norbuprenorphine Level 1 Control (7 ng/mL)		
<b>0275</b>	Norbuprenorphine Cutoff Calibrator (10 ng/mL)	<b>0390</b>	<small>510(k) Cleared</small> Hydrocodone 300 EIA Kit, small test kit
<b>0276</b>	Norbuprenorphine Level 2 Control (13 ng/mL)	<b>0391</b>	Hydrocodone 300 EIA Kit, large test kit
<b>0277</b>	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)	<b>0392</b>	Hydrocodone 300 Low Calibrator (150 ng/mL)
<b>0278</b>	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)	<b>0393</b>	Hydrocodone 300 Cutoff Calibrator (300 ng/mL)
<b>0279</b>	Norbuprenorphine High Calibrator (75 ng/mL)	<b>0394</b>	Hydrocodone 300 Intermediate Calibrator (500 ng/mL)
		<b>0395</b>	Hydrocodone 300 High Calibrator (800 ng/mL)
<b>0290</b>	<small>510(k) Cleared</small> 6-Acetylmorphine EIA Kit, small test kit	<b>0397</b>	Hydrocodone 300 Level 1 Control (225 ng/mL)
<b>0291</b>	6-Acetylmorphine EIA Kit, large test kit	<b>0398</b>	Hydrocodone 300 Level 2 Control (375 ng/mL)
<b>0299</b>	6-Acetylmorphine Negative Calibrator (0 ng/mL)	<b>0400</b>	SPICE I (JWH-018) EIA Kit, small test kit <sup>†</sup>
<b>0292</b>	6-Acetylmorphine Low Calibrator (5 ng/mL)	<b>0401</b>	SPICE I (JWH-018) EIA Kit, large test kit <sup>†</sup>
<b>0293</b>	6-Acetylmorphine Cutoff Calibrator (10 ng/mL)	<b>0002c</b>	THC Negative Calibrator (0 ng/mL)
<b>0294</b>	6-Acetylmorphine Intermediate Calibrator (20 ng/mL)	<b>0402</b>	SPICE I (JWH-018) Low Calibrator (10 ng/mL) <sup>†</sup>
<b>0295</b>	6-Acetylmorphine High Calibrator (40 ng/mL)	<b>0403</b>	SPICE I (JWH-018) Cutoff Calibrator (20 ng/mL) <sup>†</sup>
<b>0297</b>	6-Acetylmorphine Level 1 Control (7.5 ng/mL)	<b>0404</b>	SPICE I (JWH-018) Intermediate Calibrator (35 ng/mL) <sup>†</sup>
<b>0298</b>	6-Acetylmorphine Level 2 Control (12.5 ng/mL)	<b>0405</b>	SPICE I (JWH-018) High Calibrator (50 ng/mL) <sup>†</sup>
		<b>0407</b>	SPICE I (JWH-018) Level 1 Control (15 ng/mL) <sup>†</sup>
<b>0310</b>	<small>510(k) Cleared</small> Fentanyl EIA Kit (Q), small test kit	<b>0408</b>	SPICE I (JWH-018) Level 2 Control (25 ng/mL) <sup>†</sup>
<b>0311</b>	Fentanyl EIA Kit (Q), large test kit		
<b>0313</b>	Norfentanyl Cutoff Calibrator (5 ng/mL)	<b>0410</b>	<small>510(k) Cleared</small> Tramadol EIA Kit, small test kit
<b>0317</b>	Norfentanyl Level 1 Control (3.75 ng/mL)	<b>0411</b>	Tramadol EIA Kit, large test kit
<b>0318</b>	Norfentanyl Level 2 Control (6.25 ng/mL)	<b>0412</b>	Tramadol Low Calibrator (50 ng/mL)
<b>0330</b>	<small>510(k) Cleared</small> Opiate 2000 EIA Kit, small test kit	<b>0413</b>	Tramadol Cutoff Calibrator (100 ng/mL)
<b>0331</b>	Opiate 2000 EIA Kit, large test kit	<b>0414</b>	Tramadol Intermediate Calibrator (225 ng/mL)
<b>0332</b>	Opiate 2000 Low Calibrator (Morphine 1000 ng/mL)	<b>0415</b>	Tramadol High Calibrator (400 ng/mL)
<b>0333</b>	Opiate 2000 Cutoff Calibrator (Morphine 2000 ng/mL)	<b>0417</b>	Tramadol Level 1 Control (75 ng/mL)
<b>0334</b>	Opiate 2000 Intermediate Calibrator (Morphine 4000 ng/mL)	<b>0418</b>	Tramadol Level 2 Control (125 ng/mL)
<b>0335</b>	Opiate 2000 High Calibrator (Morphine 6000 ng/mL)		
<b>0337</b>	Opiate 2000 Level 1 Control (Morphine 1500 ng/mL)		
<b>0338</b>	Opiate 2000 Level 2 Control (Morphine 2500 ng/mL)		

510(k) Cleared = FDA 510(k) Cleared Assay, Calibrators, & Controls.

<sup>†</sup> Forensic Use Only

<sup>§</sup> Product is made to order. Please allow 10 to 12 weeks of production time. Minimum purchase volume required.

\* EIA = Enzyme Immunoassay

\* All single analyte calibrators and controls are 5 mL per bottle.

\* A small test kit is 100/37.5 mL R<sub>1</sub>/R<sub>2</sub>

\* A large test kit is 1000/375 mL R<sub>1</sub>/R<sub>2</sub>



**PRODUCT CATALOG: Drugs of Abuse Urinalysis by Catalog Number**

<b>REF #</b>	<b>Product Name</b>	<b>REF #</b>	<b>Product Name</b>
0420	SPICE II (UR-144/XLR-11) EIA Kit, small test kit <sup>†</sup>	0630	<span style="border: 1px solid black; padding: 1px;">510(k) Cleared</span> Methadone II EIA Kit, small test kit
0421	SPICE II (UR-144/XLR-11) EIA Kit, large test kit <sup>†</sup>	0631	Methadone II EIA Kit, large test kit
0002c	THC Negative Calibrator (0 ng/mL)	0112	Methadone Low Calibrator (150 ng/mL)
0422	SPICE II (UR-144/XLR-11) Low Calibrator (10 ng/mL) <sup>†</sup>	0113	Methadone Cutoff Calibrator (300 ng/mL)
0423	SPICE II (UR-144/XLR-11) Cutoff Calibrator (20 ng/mL) <sup>†</sup>	0114	Methadone Intermediate Calibrator (600 ng/mL)
0424	SPICE II (UR-144/XLR-11) Intermediate Cal (35 ng/mL) <sup>†</sup>	0115	Methadone High Calibrator (1000 ng/mL)
0425	SPICE II (UR-144/XLR-11) High Calibrator (50 ng/mL) <sup>†</sup>	0117	Methadone Level 1 Control (225 ng/mL)
0427	SPICE II (UR-144/XLR-11) Level 1 Control (15 ng/mL) <sup>†</sup>	0118	Methadone Level 2 Control (375 ng/mL)
0428	SPICE II (UR-144/XLR-11) Level 2 Control (25 ng/mL) <sup>†</sup>	0640	Ketamine EIA Kit, small test kit <sup>†</sup>
0450	Ethyl Glucuronide III EIA Kit, small test kit <sup>†</sup>	0641	Ketamine EIA Kit, large test kit <sup>†</sup>
0451	Ethyl Glucuronide III EIA Kit, large test kit <sup>†</sup>	0642	Norketamine Low Calibrator (25 ng/mL) <sup>†</sup>
	<b>EtG 200 for EtG III EIA Kit</b>	0643	Norketamine Cutoff Calibrator #1 (50 ng/mL) <sup>†</sup>
0452	Ethyl Glucuronide 200 Low Calibrator (100 ng/mL) <sup>†</sup>	0644	Norketamine Cutoff Calibrator #2 (100 ng/mL) <sup>†</sup>
0453	Ethyl Glucuronide 200 Cutoff Calibrator (200 ng/mL) <sup>†</sup>	0645	Norketamine Intermediate Calibrator (250 ng/mL) <sup>†</sup>
0454	Ethyl Glucuronide 200 Intermediate Calibrator (500 ng/mL) <sup>†</sup>	0646	Norketamine High Calibrator (500 ng/mL) <sup>†</sup>
0455	Ethyl Glucuronide 200 High Calibrator (1000 ng/mL) <sup>†</sup>		<b>KET 50</b>
0457	Ethyl Glucuronide 200 Level 1 Control (150 ng/mL) <sup>†</sup>	0647	Norketamine Level 1 Control (37.5 ng/mL) <sup>†</sup>
0458	Ethyl Glucuronide 200 Level 2 Control (250 ng/mL) <sup>†</sup>	0648	Norketamine Level 2 Control (62.5 ng/mL) <sup>†</sup>
	<b>EtG 500 for EtG III EIA Kit</b>		<b>KET 100</b>
0462	Ethyl Glucuronide 500 Low Calibrator (200 ng/mL) <sup>†</sup>	0657	Norketamine Level 1 Control (75 ng/mL) <sup>†</sup>
0463	Ethyl Glucuronide 500 Cutoff Calibrator (500 ng/mL) <sup>†</sup>	0658	Norketamine Level 2 Control (125 ng/mL) <sup>†</sup>
0464	Ethyl Glucuronide 500 Intermediate Calibrator (1000 ng/mL) <sup>†</sup>	0852	<span style="border: 1px solid black; padding: 1px;">510(k) Cleared</span> Multi-Analyte Set C Low Calibrator, 15 mL <sup>§</sup>
0465	Ethyl Glucuronide 500 High Calibrator (2000 ng/mL) <sup>†</sup>	0853	Multi-Analyte Set C Cutoff Calibrator, 15 mL <sup>§</sup>
0467	Ethyl Glucuronide 500 Level 1 Control (375 ng/mL) <sup>†</sup>	0854	Multi-Analyte Set C Intermediate Calibrator, 15 mL <sup>§</sup>
0468	Ethyl Glucuronide 500 Level 2 Control (625 ng/mL) <sup>†</sup>	0855	Multi-Analyte Set C High Calibrator, 15 mL <sup>§</sup>
0520	<span style="border: 1px solid black; padding: 1px;">510(k) Cleared</span> Cotinine II EIA Kit, small test kit	0857	Multi-Analyte Set C Level 1 Control, 15 mL <sup>§</sup>
0521	Cotinine II EIA Kit, large test kit	0858	Multi-Analyte Set C Level 2 Control, 15 mL <sup>§</sup>
0522	Cotinine II Low Calibrator (100 ng/mL)		<i>Tests COC 150, MAMP, AMP 500, MTD, OPI2K, BZO 200, and BARB 200</i>
0523	Cotinine II Cutoff Calibrator (200 ng/mL)	0871	Multi-Analyte Set D Negative Calibrator, 15 mL <sup>§</sup>
0524	Cotinine II Intermediate Calibrator (400 ng/mL)	0872	Multi-Analyte Set D Low Calibrator, 15 mL <sup>§</sup>
0525	Cotinine II High Calibrator (1000 ng/mL)	0873	Multi-Analyte Set D Cutoff Calibrator, 15 mL <sup>§</sup>
0527	Cotinine II Level 1 Control (150 ng/mL)	0874	Multi-Analyte Set D Intermediate Calibrator, 15 mL <sup>§</sup>
0528	Cotinine II Level 2 Control (250 ng/mL)	0875	Multi-Analyte Set D High Calibrator, 15 mL <sup>§</sup>
0550	Fentanyl EIA Kit (SQ), small test kit <sup>†</sup>	0877	Multi-Analyte Set D Level 1 Control, 15 mL <sup>§</sup>
0551	Fentanyl EIA Kit (SQ), large test kit <sup>†</sup>	0878	Multi-Analyte Set D Level 2 Control, 15 mL <sup>§</sup>
0552	Norfentanyl Low Calibrator (2.5 ng/mL) <sup>†</sup>		<i>Tests AMP 1000, BZO 200, COC 300, OPI 300, and OXY 100</i>
0553	Norfentanyl Cutoff Calibrator (5 ng/mL) <sup>†</sup>	0891	Multi-Analyte Set E Negative Calibrator, 15 mL <sup>§</sup>
0554	Norfentanyl Intermediate Calibrator (10 ng/mL) <sup>†</sup>	0892	Multi-Analyte Set E Low Calibrator, 15 mL <sup>§</sup>
0555	Norfentanyl High Calibrator (20 ng/mL) <sup>†</sup>	0893	Multi-Analyte Set E Cutoff Calibrator, 15 mL <sup>§</sup>
0557	Norfentanyl Level 1 Control (3.75 ng/mL) <sup>†</sup>	0894	Multi-Analyte Set E Intermediate Calibrator, 15 mL <sup>§</sup>
0558	Norfentanyl Level 2 Control (6.25 ng/mL) <sup>†</sup>	0895	Multi-Analyte Set E High Calibrator, 15 mL <sup>§</sup>
0610	<span style="border: 1px solid black; padding: 1px;">510(k) Cleared</span> Oxycodone III EIA Kit, small test kit	0897	Multi-Analyte Set E Level 1 Control, 15 mL <sup>§</sup>
0611	Oxycodone III EIA Kit, large test kit	0898	Multi-Analyte Set E Level 2 Control, 15 mL <sup>§</sup>
	<b>OXY III 100</b>		<i>Tests BARB 200, BUP 10, MTD 300, and PCP 25</i>
0612	Oxycodone 100 Low Calibrator (50 ng/mL)		
0613	Oxycodone 100 Cutoff Calibrator (100 ng/mL)		
0614	Oxycodone 100 Intermediate Calibrator (150 ng/mL)		
0615	Oxycodone 100 High Calibrator (300 ng/mL)		
0617	Oxycodone 100 Level 1 Control (75 ng/mL)		
0618	Oxycodone 100 Level 2 Control (125 ng/mL)		
	<b>OXY III 300</b>		
0622	Oxycodone 300 Low Calibrator (150 ng/mL)		
0623	Oxycodone 300 Cutoff Calibrator (300 ng/mL)		
0624	Oxycodone 300 Intermediate Calibrator (500 ng/mL)		
0625	Oxycodone 300 High Calibrator (800 ng/mL)		
0627	Oxycodone 300 Level 1 Control (225 ng/mL)		
0628	Oxycodone 300 Level 2 Control (375 ng/mL)		

510(k)  
Cleared = FDA 510(k) Cleared Assay, Calibrators, & Controls.

<sup>†</sup> Forensic Use Only

<sup>§</sup> Product is made to order. Please allow 10 to 12 weeks of production time. Minimum purchase volume required.

\* EIA = Enzyme Immunoassay

\* All single analyte calibrators and controls are 5 mL per bottle.

\* A small test kit is 100/37.5 mL R<sub>1</sub>/R<sub>2</sub>

\* A large test kit is 1000/375 mL R<sub>1</sub>/R<sub>2</sub>