



PRODUCT CATALOG: Drugs of Abuse Oral Fluid by Catalog Number

<u>REF #</u>	<u>Product Name</u>	<u>REF #</u>	<u>Product Name</u>
S0005	LZI Oral Fluid Collector, 100 tubes/case	S0130c	Benzodiazepine EIA Kit, small test kit ^{†§}
S0010c	Phencyclidine EIA Kit, small test kit ^{†§}	S0131c	Benzodiazepine EIA Kit, large test kit ^{†§}
S0011c	Phencyclidine EIA Kit, large test kit ^{†§}	S0008	Oral Fluid LIS Negative Calibrator ^{†§}
S0008	Oral Fluid LIS Negative Calibrator ^{†§}	S0133c	Benzodiazepine Cutoff Calibrator (Oxazepam 20 ng/mL) ^{†§}
S0012c	Phencyclidine Low Calibrator (5 ng/mL) ^{†§}	S0135c	Benzodiazepine High Calibrator (Oxazepam 50 ng/mL) ^{†§}
S0013c	Phencyclidine Cutoff Calibrator (10 ng/mL) ^{†§}	S0137c	Benzodiazepine Low Calibrator/Level 1 Control (Oxazepam 10 ng/mL) ^{†§}
S0014c	Phencyclidine Intermediate Calibrator (15 ng/mL) ^{†§}	S0138c	Benzodiazepine Intermediate Cal/Level 2 Control (Oxazepam 30 ng/mL) ^{†§}
S0015c	Phencyclidine High Calibrator (25 ng/mL) ^{†§}	S0140b	Barbiturate EIA Kit, small test kit ^{†§}
S0017c	Phencyclidine Level 1 Control (7.5 ng/mL) ^{†§}	S0141b	Barbiturate EIA Kit, large test kit ^{†§}
S0018c	Phencyclidine Level 2 Control (12.5 ng/mL) ^{†§}	S0001	Oral Fluid Negative Calibrator ^{†§}
S0020c	Opiate EIA Kit, small test kit ^{†§}	S0143b	Barbiturate Cutoff Calibrator (Secobarbital 20 ng/mL) ^{†§}
S0021c	Opiate EIA Kit, large test kit ^{†§}	S0145b	Barbiturate High Calibrator (Secobarbital 50 ng/mL) ^{†§}
S0008	Oral Fluid LIS Negative Calibrator ^{†§}	S0147b	Barbiturate Low Calibrator/Level 1 Control (Secobarbital 10 ng/mL) ^{†§}
S0022c	Opiate Low Calibrator (Morphine 20 ng/mL) ^{†§}	S0148b	Barbiturate Intermediate Cal/Level 2 Control (Secobarbital 30 ng/mL) ^{†§}
S0023c	Opiate Cutoff Calibrator (Morphine 40 ng/mL) ^{†§}	S0220	Ethyl Alcohol Assay Kit, small test kit ^{†§}
S0024c	Opiate Intermediate Calibrator (Morphine 60 ng/mL) ^{†§}	S0221	Ethyl Alcohol Assay Kit, large test kit ^{†§}
S0025c	Opiate High Calibrator (Morphine 80 ng/mL) ^{†§}	S0001	Oral Fluid Negative Calibrator ^{†§}
S0027c	Opiate Level 1 Control (Morphine 30 ng/mL) ^{†§}	S0223	Ethyl Alcohol Assay Calibrator (50 mg/dL) ^{†§}
S0028c	Opiate Level 2 Control (Morphine 50 ng/mL) ^{†§}	S0226	Ethyl Alcohol Assay Calibrator/Control (25 mg/dL) ^{†§}
S0030d	Cocaine Metabolite EIA Kit, small test kit ^{†§}	S0227	Ethyl Alcohol Assay Calibrator/Control (75 mg/dL) ^{†§}
S0031d	Cocaine Metabolite EIA Kit, large test kit ^{†§}	S0228	Ethyl Alcohol Assay Calibrator (100 mg/dL) ^{†§}
S0007d	Oral Fluid LC Negative Calibrator ^{†§}	S0290	^{§10k} _{clear} 6-Acetylmorphine EIA Kit, small test kit [§]
S0033d	Cocaine Metabolite Low Calibrator (BE 10 ng/mL) ^{†§}	S0291	6-Acetylmorphine EIA Kit, large test kit [§]
S0034d	Cocaine Metabolite Cutoff Calibrator (BE 20 ng/mL) ^{†§}	S0298	6-Acetylmorphine Negative Calibrator (0 ng/mL) [§]
S0035d	Cocaine Metabolite Intermediate Calibrator (BE 35 ng/mL) ^{†§}	S0292	6-Acetylmorphine Low Calibrator/Level 1 Control (2 ng/mL) [§]
S0036d	Cocaine Metabolite High Calibrator (BE 50 ng/mL) ^{†§}	S0293	6-Acetylmorphine Cutoff Calibrator (4 ng/mL) [§]
S0037d	Cocaine Metabolite Level 1 Control (BE 15 ng/mL) ^{†§}	S0294	6-Acetylmorphine Intermediate Calibrator (10 ng/mL) [§]
S0038d	Cocaine Metabolite Level 2 Control (BE 25 ng/mL) ^{†§}	S0295	6-Acetylmorphine High Calibrator (20 ng/mL) [§]
S0040b	^{§10k} _{clear} Amphetamine EIA Kit, small test kit [§]	S0297	6-Acetylmorphine Level 2 Control (6 ng/mL) [§]
S0041b	Amphetamine EIA Kit, large test kit [§]	S0530	Cotinine II EIA Kit, small test kit ^{†§}
S0001	Oral Fluid Negative Calibrator [§]	S0531	Cotinine II EIA Kit, large test kit ^{†§}
S0042b	Amphetamine Low Calibrator (20 ng/mL) [§]	S0001	Oral Fluid Negative Calibrator ^{†§}
S0043b	Amphetamine Cutoff Calibrator (50 ng/mL) [§]	S0532	Cotinine II Low Calibrator (12.5 ng/mL) ^{†§}
S0044b	Amphetamine Intermediate Calibrator (100 ng/mL) [§]	S0533	Cotinine II Cutoff Calibrator (25 ng/mL) ^{†§}
S0045b	Amphetamine High Calibrator (140 ng/mL) [§]	S0534	Cotinine II Intermediate Calibrator (50 ng/mL) ^{†§}
S0046b	Amphetamine Level 1 Control (37.5 ng/mL) [§]	S0535	Cotinine II High Calibrator (75 ng/mL) ^{†§}
S0047b	Amphetamine Level 2 Control (62.5 ng/mL) [§]	S0537	Cotinine II Level 1 Control (18.75 ng/mL) ^{†§}
S0050b	^{§10k} _{clear} Methamphetamine EIA Kit, small test kit [§]	S0538	Cotinine II Level 2 Control (31.25 ng/mL) ^{†§}
S0051b	Methamphetamine EIA Kit, large test kit [§]		
S0001	Oral Fluid Negative Calibrator [§]		
S0052b	Methamphetamine Low Calibrator (20 ng/mL) [§]		
S0053b	Methamphetamine Cutoff Calibrator (50 ng/mL) [§]		
S0054b	Methamphetamine Intermediate Calibrator (100 ng/mL) [§]		
S0055b	Methamphetamine High Calibrator (140 ng/mL) [§]		
S0056b	Methamphetamine Level 1 Control (37.5 ng/mL) [§]		
S0057b	Methamphetamine Level 2 Control (62.5 ng/mL) [§]		

^{§10k}_{clear} = FDA Approved Assay, Calibrators, & Controls.

[†] Forensic Use Only (Non-FDA Approved)

[§] Product is made to order. Please allow 10 to 12 weeks of production time. Minimum purchase volume required.

* The Negative Calibrator to be used with each specific assay has been listed with each assay set.

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 75/37.5 mL R₁/R₂

* A large test kit is 750/375 mL R₁/R₂