

**PRODUCT CATALOG: Drugs of Abuse Urinalysis – EU ONLY by Catalog Number**

REF #	Product Name	REF #	Product Name
0001	Universal Negative Calibrator, 5 mL	0100 STOKK Cleared	Amphetamines 500 EIA Kit, small test kit
0010 STOKK Cleared	Phencyclidine EIA Kit, small test kit	0101	Amphetamines 500 EIA Kit, large test kit
0011	Phencyclidine EIA Kit, large test kit	0102	Amphetamines 500 Low Calibrator (MAMP 250 ng/mL)
0012	Phencyclidine Low Calibrator (12.5 ng/mL)	0103	Amphetamines 500 Cutoff Calibrator (MAMP 500 ng/mL)
0013	Phencyclidine Cutoff Calibrator (25 ng/mL)	0104	Amphetamines 500 Intermediate Calibrator (MAMP 1000 ng/mL)
0014	Phencyclidine Intermediate Calibrator (50 ng/mL)	0105	Amphetamines 500 High Calibrator (MAMP 2000 ng/mL)
0015	Phencyclidine High Calibrator (100 ng/mL)	0107	Amphetamines 500 Level 1 Control (MAMP 375 ng/mL)
0017	Phencyclidine Level 1 Control (18 ng/mL)	0108	Amphetamines 500 Level 2 Control (MAMP 625 ng/mL)
0018	Phencyclidine Level 2 Control (32 ng/mL)	0120 STOKK Cleared	Propoxyphene EIA Kit, small test kit
0020b STOKK Cleared	Opiate EIA Kit, small test kit	0121	Propoxyphene EIA Kit, large test kit
0021b	Opiate EIA Kit, large test kit	0122	Propoxyphene Low Calibrator (150 ng/mL)
0022	Opiate Low Calibrator (Morphine 150 ng/mL)	0123	Propoxyphene Cutoff Calibrator (300 ng/mL)
0023	Opiate Cutoff Calibrator (Morphine 300 ng/mL)	0124	Propoxyphene Intermediate Calibrator (600 ng/mL)
0024	Opiate Intermediate Calibrator (Morphine 600 ng/mL)	0125	Propoxyphene High Calibrator (1000 ng/mL)
0025	Opiate High Calibrator (Morphine 1000 ng/mL)	0127	Propoxyphene Level 1 Control (225 ng/mL)
0027	Opiate Level 1 Control (Morphine 225 ng/mL)	0128	Propoxyphene Level 2 Control (375 ng/mL)
0028	Opiate Level 2 Control (Morphine 375 ng/mL)	0130 STOKK Cleared	Benzodiazepine EIA Kit, small test kit
0030 STOKK Cleared	Cocaine Metabolite EIA Kit, small test kit	0131	Benzodiazepine EIA Kit, large test kit
0031	Cocaine Metabolite EIA Kit, large test kit	BZO 200	
0032	Cocaine Metabolite Low Calibrator (BE 150 ng/mL)	0132	Benzodiazepine Calibrator #2/Control I (Oxazepam 100 ng/mL)
0033	Cocaine Metabolite Cutoff Calibrator (BE 300 ng/mL)	0133	Benzodiazepine Calibrator #3/Cutoff A (Oxazepam 200 ng/mL)
0034	Cocaine Metabolite Intermediate Cal (BE 1000 ng/mL)	0134	Benzodiazepine Calibrator #4/Control III (Oxazepam 300 ng/mL)
0035	Cocaine Metabolite High Calibrator (BE 3000 ng/mL)	0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)
0037	Cocaine Metabolite Level 1 Control (BE 225 ng/mL)	BZO 300	
0038	Cocaine Metabolite Level 2 Control (BE 375 ng/mL)	0132	Benzodiazepine Calibrator #2 (Oxazepam 100 ng/mL)
0040 STOKK Cleared	Amphetamines EIA Kit, small test kit	0133	Benzodiazepine Calibrator #3/Control II (Oxazepam 200 ng/mL)
0041	Amphetamines EIA Kit, large test kit	0134	Benzodiazepine Calibrator #4/Cutoff B (Oxazepam 300 ng/mL)
0042	Amphetamines Low Calibrator (MAMP 500 ng/mL)	0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)
0043	Amphetamines Cutoff Calibrator (MAMP 1000 ng/mL)	0136	Benzodiazepine Control IV (Oxazepam 400 ng/mL)
0044	Amphetamines Intermediate Calibrator (MAMP 1500 ng/mL)	0140 STOKK Cleared	Barbiturate EIA Kit, small test kit
0045	Amphetamines High Calibrator (MAMP 2000 ng/mL)	0141	Barbiturate EIA Kit, large test kit
0047	Amphetamines Level 1 Control (MAMP 750 ng/mL)	BARB 200	
0048	Amphetamines Level 2 Control (MAMP 1250 ng/mL)	0142	Barbiturate Calibrator #2/Control I (Secobarbital 100 ng/mL)
0070c STOKK Cleared	Cannabinoids (cTHC) EIA Kit, small test kit	0143	Barbiturate Calibrator #3/Cutoff A (Secobarbital 200 ng/mL)
0071c	Cannabinoids (cTHC) EIA Kit, large test kit	0144	Barbiturate Calibrator #4/Control III (Secobarbital 300 ng/mL)
0002c	Cannabinoids (cTHC) Negative Calibrator (0 ng/mL)	0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)
cTHC 25		BARB 300	
0072c	Cannabinoids (cTHC) Low Calibrator (12.5 ng/mL)	0142	Barbiturate Calibrator #2 (Secobarbital 100 ng/mL)
0073c	Cannabinoids (cTHC) Cutoff Calibrator (25 ng/mL)	0143	Barbiturate Calibrator #3/Control II (Secobarbital 200 ng/mL)
0007c	Cannabinoids (cTHC) Intermediate Calibrator (37.5 ng/mL)	0144	Barbiturate Calibrator #4/Cutoff B (Secobarbital 300 ng/mL)
0075c	Cannabinoids (cTHC) High Calibrator (50 ng/mL)	0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)
0006c	Cannabinoids (cTHC) Level 1 Control (18.75 ng/mL)	0146	Barbiturate Control IV (Secobarbital 400 ng/mL)
0074c	Cannabinoids (cTHC) Level 2 Control (31.25 ng/mL)	0160 STOKK Cleared	Ecstasy (MDMA) EIA Kit, small test kit
cTHC 50		0161	Ecstasy (MDMA) EIA Kit, large test kit
0073c	Cannabinoids (cTHC) Low Calibrator (25 ng/mL)	0162	Ecstasy Low Calibrator (100 ng/mL)
0075c	Cannabinoids (cTHC) Cutoff Calibrator (50 ng/mL)	0163	Ecstasy Cutoff Calibrator (500 ng/mL)
0076c	Cannabinoids (cTHC) Intermediate Calibrator (75 ng/mL)	0164	Ecstasy Intermediate Calibrator (750 ng/mL)
0077c	Cannabinoids (cTHC) High Calibrator (100 ng/mL)	0165	Ecstasy High Calibrator (1000 ng/mL)
0007c	Cannabinoids (cTHC) Level 1 Control (37.5 ng/mL)	0167	Ecstasy Level 1 Control (375 ng/mL)
0008c	Cannabinoids (cTHC) Level 2 Control (62.5 ng/mL)	0168	Ecstasy Level 2 Control (625 ng/mL)
cTHC 100			
0075c	Cannabinoids (cTHC) Low Calibrator (50 ng/mL)		
0077c	Cannabinoids (cTHC) Cutoff Calibrator (100 ng/mL)		
0078c	Cannabinoids (cTHC) Intermediate Calibrator (150 ng/mL)		
0079c	Cannabinoids (cTHC) High Calibrator (200 ng/mL)		
0076c	Cannabinoids (cTHC) Level 1 Control (75 ng/mL)		
0009c	Cannabinoids (cTHC) Level 2 Control (125 ng/mL)		

STOKK Cleared = U.S. FDA 510(k) Cleared Assay, Calibrators, & Controls.^s Product is made to order. Please allow 10 to 12 weeks of production time.
Minimum purchase volume required.

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂* A large test kit is 1000/375 mL R₁/R₂

**PRODUCT CATALOG: Drugs of Abuse Urinalysis – EU ONLY by Catalog Number**

REF #	Product Name	REF #	Product Name
0190b 510(k) Cleared	Methadone Metabolite (EDDP) EIA Kit, small test kit	0340 510(k) Cleared	Cocaine Metabolite 150 EIA Kit, small test kit
0191b	Methadone Metabolite (EDDP) EIA Kit, large test kit	0341	Cocaine Metabolite 150 EIA Kit, large test kit
	EDDPb 100	0342	Cocaine Metabolite 150 Low Calibrator (BE 75 ng/mL)
0202b	Methadone Metabolite 100 Low Calibrator (50 ng/mL)	0343	Cocaine Metabolite 150 Cutoff Calibrator (BE 150 ng/mL)
0203b	Methadone Metabolite 100 Cutoff Calibrator (100 ng/mL)	0344	Cocaine Metabolite 150 Intermediate Calibrator (BE 300 ng/mL)
0204b	Methadone Metabolite 100 Intermediate Calibrator (250 ng/mL)	0345	Cocaine Metabolite 150 High Calibrator (BE 1000 ng/mL)
0205b	Methadone Metabolite 100 High Calibrator (500 ng/mL)	0347	Cocaine Metabolite 150 Level 1 Control (BE 112.5 ng/mL)
0207b	Methadone Metabolite 100 Level 1 Control (75 ng/mL)	0348	Cocaine Metabolite 150 Level 2 Control (BE 187.5 ng/mL)
0208b	Methadone Metabolite 100 Level 2 Control (125 ng/mL)		
	EDDPb 300	0350 510(k) Cleared	Methamphetamine EIA Kit, small test kit
0192b	Methadone Metabolite 300 Low Calibrator (150 ng/mL)	0351	Methamphetamine EIA Kit, large test kit
0193b	Methadone Metabolite 300 Cutoff Calibrator (300 ng/mL)	0352	Methamphetamine Low Calibrator (250 ng/mL)
0194b	Methadone Metabolite 300 Intermediate Calibrator (600 ng/mL)	0353	Methamphetamine Cutoff Calibrator (500 ng/mL)
0195b	Methadone Metabolite 300 High Calibrator (1000 ng/mL)	0354	Methamphetamine Intermediate Calibrator (1000 ng/mL)
0197b	Methadone Metabolite 300 Level 1 Control (225 ng/mL)	0355	Methamphetamine High Calibrator (2000 ng/mL)
0198b	Methadone Metabolite 300 Level 2 Control (375 ng/mL)	0357	Methamphetamine Level 1 Control (375 ng/mL)
		0358	Methamphetamine Level 2 Control (625 ng/mL)
0220 510(k) Cleared	Ethyl Alcohol Enzymatic Assay Kit, small test kit	0370 510(k) Cleared	Carisoprodol Metabolite (SOMA) EIA Kit, small test kit ^s
0221	Ethyl Alcohol Enzymatic Assay Kit, large test kit	0371	Carisoprodol Metabolite (SOMA) EIA Kit, large test kit ^s
0223	Ethyl Alcohol Calibrator (100 mg/dL)	0372	Meprobamate Low Calibrator (50 ng/mL) ^s
0224	Ethyl Alcohol Level 1 Control (50 mg/dL)	0373	Meprobamate Cutoff Calibrator (100 ng/mL) ^s
0225	Ethyl Alcohol Level 2 Control (300 mg/dL)	0374	Meprobamate Intermediate Calibrator (200 ng/mL) ^s
		0375	Meprobamate High Calibrator (400 ng/mL) ^s
0270 510(k) Cleared	Buprenorphine EIA Kit, small test kit	0377	Meprobamate Level 1 Control (75 ng/mL) ^s
0271	Buprenorphine EIA Kit, large test kit	0378	Meprobamate Level 2 Control (125 ng/mL) ^s
	BUP 5		
0272	Norbuprenorphine Level 1 Control (3 ng/mL)	0380 510(k) Cleared	Hydrocodone 100 EIA Kit, small test kit
0273	Norbuprenorphine Cutoff Calibrator #1 (5 ng/mL)	0381	Hydrocodone 100 EIA Kit, large test kit
0274	Norbuprenorphine Level 2 Control (7 ng/mL)	0382	Hydrocodone 100 Low Calibrator (50 ng/mL)
0275	Norbuprenorphine Cutoff Calibrator #2 (10 ng/mL)	0383	Hydrocodone 100 Cutoff Calibrator (100 ng/mL)
0277	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)	0384	Hydrocodone 100 Intermediate Calibrator (150 ng/mL)
0278	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)	0385	Hydrocodone 100 High Calibrator (300 ng/mL)
0279	Norbuprenorphine High Calibrator (75 ng/mL)	0387	Hydrocodone 100 Level 1 Control (75 ng/mL)
	BUP 10	0388	Hydrocodone 100 Level 2 Control (125 ng/mL)
0273	Norbuprenorphine Low Calibrator (5 ng/mL)		
0274	Norbuprenorphine Level 1 Control (7 ng/mL)	0390 510(k) Cleared	Hydrocodone 300 EIA Kit, small test kit
0275	Norbuprenorphine Cutoff Calibrator (10 ng/mL)	0391	Hydrocodone 300 EIA Kit, large test kit
0276	Norbuprenorphine Level 2 Control (13 ng/mL)	0392	Hydrocodone 300 Low Calibrator (150 ng/mL)
0277	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)	0393	Hydrocodone 300 Cutoff Calibrator (300 ng/mL)
0278	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)	0394	Hydrocodone 300 Intermediate Calibrator (500 ng/mL)
0279	Norbuprenorphine High Calibrator (75 ng/mL)	0395	Hydrocodone 300 High Calibrator (800 ng/mL)
		0397	Hydrocodone 300 Level 1 Control (225 ng/mL)
0290 510(k) Cleared	6-Acetylmorphine EIA Kit, small test kit	0398	Hydrocodone 300 Level 2 Control (375 ng/mL)
0291	6-Acetylmorphine EIA Kit, large test kit		
0299	6-Acetylmorphine Negative Calibrator (0 ng/mL)	0410 510(k) Cleared	Tramadol EIA Kit, small test kit ^s
0292	6-Acetylmorphine Low Calibrator (5 ng/mL)	0411	Tramadol EIA Kit, large test kit ^s
0293	6-Acetylmorphine Cutoff Calibrator (10 ng/mL)	0412	Tramadol Low Calibrator (50 ng/mL) ^s
0294	6-Acetylmorphine Intermediate Calibrator (20 ng/mL)	0413	Tramadol Cutoff Calibrator (100 ng/mL) ^s
0295	6-Acetylmorphine High Calibrator (40 ng/mL)	0414	Tramadol Intermediate Calibrator (225 ng/mL) ^s
0297	6-Acetylmorphine Level 1 Control (7.5 ng/mL)	0415	Tramadol High Calibrator (400 ng/mL) ^s
0298	6-Acetylmorphine Level 2 Control (12.5 ng/mL)	0417	Tramadol Level 1 Control (75 ng/mL) ^s
		0418	Tramadol Level 2 Control (125 ng/mL) ^s
0330 510(k) Cleared	Opiate 2000 EIA Kit, small test kit		
0331	Opiate 2000 EIA Kit, large test kit	0500	SPICE I (JWH-018) EIA Kit, small test kit
0332	Opiate 2000 Low Calibrator (Morphine 1000 ng/mL)	0501	SPICE I (JWH-018) EIA Kit, large test kit
0333	Opiate 2000 Cutoff Calibrator (Morphine 2000 ng/mL)	0002c	THC Negative Calibrator (0 ng/mL)
0334	Opiate 2000 Intermediate Calibrator (Morphine 4000 ng/mL)	0502	SPICE I (JWH-018) Low Calibrator (10 ng/mL)
0335	Opiate 2000 High Calibrator (Morphine 6000 ng/mL)	0503	SPICE I (JWH-018) Cutoff Calibrator (20 ng/mL)
0337	Opiate 2000 Level 1 Control (Morphine 1500 ng/mL)	0504	SPICE I (JWH-018) Intermediate Calibrator (35 ng/mL)
0338	Opiate 2000 Level 2 Control (Morphine 2500 ng/mL)	0505	SPICE I (JWH-018) High Calibrator (50 ng/mL)
		0507	SPICE I (JWH-018) Level 1 Control (15 ng/mL)
		0508	SPICE I (JWH-018) Level 2 Control (25 ng/mL)

510(k) Cleared = U.S. FDA 510(k) Cleared Assay, Calibrators, & Controls.^s Product is made to order. Please allow 10 to 12 weeks of production time.
Minimum purchase volume required.

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂* A large test kit is 1000/375 mL R₁/R₂

**PRODUCT CATALOG: Drugs of Abuse Urinalysis – EU ONLY by Catalog Number**

REF #	Product Name	REF #	Product Name
0510	SPICE II (UR-144/XLR-11) EIA Kit, small test kit	0630 510(k) Cleared	Methadone II EIA Kit, small test kit
0511	SPICE II (UR-144/XLR-11) EIA Kit, large test kit	0631	Methadone II EIA Kit, large test kit
0002c	THC Negative Calibrator (0 ng/mL)	0112	Methadone Low Calibrator (150 ng/mL)
0512	SPICE II (UR-144/XLR-11) Low Calibrator (10 ng/mL)	0113	Methadone Cutoff Calibrator (300 ng/mL)
0513	SPICE II (UR-144/XLR-11) Cutoff Calibrator (20 ng/mL)	0114	Methadone Intermediate Calibrator (600 ng/mL)
0514	SPICE II (UR-144/XLR-11) Intermediate Cal (35 ng/mL)	0115	Methadone High Calibrator (1000 ng/mL)
0515	SPICE II (UR-144/XLR-11) High Calibrator (50 ng/mL)	0117	Methadone Level 1 Control (225 ng/mL)
0517	SPICE II (UR-144/XLR-11) Level 1 Control (15 ng/mL)	0118	Methadone Level 2 Control (375 ng/mL)
0518	SPICE II (UR-144/XLR-11) Level 2 Control (25 ng/mL)	0660	Ketamine EIA Kit, small test kit ^{‡§}
0520 510(k) Cleared	Cotinine II EIA Kit, small test kit	0661	Ketamine EIA Kit, large test kit ^{‡§}
0521	Cotinine II EIA Kit, large test kit	0662	Norketamine Low Calibrator (25 ng/mL) ^{‡§}
0522	Cotinine II Low Calibrator (100 ng/mL)	0663	Norketamine Cutoff Calibrator #1 (50 ng/mL) ^{‡§}
0523	Cotinine II Cutoff Calibrator (200 ng/mL)	0664	Norketamine Cutoff Calibrator #2 (100 ng/mL) ^{‡§}
0524	Cotinine II Intermediate Calibrator (400 ng/mL)	0665	Norketamine Intermediate Calibrator (250 ng/mL) ^{‡§}
0525	Cotinine II High Calibrator (1000 ng/mL)	0666	Norketamine High Calibrator (500 ng/mL) ^{‡§}
0527	Cotinine II Level 1 Control (150 ng/mL)	KET 50	
0528	Cotinine II Level 2 Control (250 ng/mL)	0667	Norketamine Level 1 Control (37.5 ng/mL) ^{‡§}
0530	Ethyl Glucuronide III EIA Kit, small test kit [§]	0668	Norketamine Level 2 Control (62.5 ng/mL) ^{‡§}
0531	Ethyl Glucuronide III EIA Kit, large test kit [§]	KET 100	
EtG III 200		0677	Norketamine Level 1 Control (75 ng/mL) ^{‡§}
0532	Ethyl Glucuronide 200 Low Calibrator (100 ng/mL) [§]	0678	Norketamine Level 2 Control (125 ng/mL) ^{‡§}
0533	Ethyl Glucuronide 200 Cutoff Calibrator (200 ng/mL) [§]	0852 510(k) Cleared	Multi-Analyte Set C Low Calibrator, 15 mL [§]
0534	Ethyl Glucuronide 200 Intermediate Calibrator (500 ng/mL) [§]	0853	Multi-Analyte Set C Cutoff Calibrator, 15 mL [§]
0535	Ethyl Glucuronide 200 High Calibrator (1000 ng/mL) [§]	0854	Multi-Analyte Set C Intermediate Calibrator, 15 mL [§]
0537	Ethyl Glucuronide 200 Level 1 Control (150 ng/mL) [§]	0855	Multi-Analyte Set C High Calibrator, 15 mL [§]
0538	Ethyl Glucuronide 200 Level 2 Control (250 ng/mL) [§]	0857	Multi-Analyte Set C Level 1 Control, 15 mL [§]
EtG III 500		0858	Multi-Analyte Set C Level 2 Control, 15 mL [§]
0542	Ethyl Glucuronide 500 Low Calibrator (200 ng/mL) [§]	<i>Tests COC 150, MAMP, AMP 500, MTD, OPI2K, BZO 200, and BARB 200</i>	
0543	Ethyl Glucuronide 500 Cutoff Calibrator (500 ng/mL) [§]	0871	Multi-Analyte Set D Negative Calibrator, 15 mL [§]
0544	Ethyl Glucuronide 500 Intermediate Calibrator (1000 ng/mL) [§]	0872	Multi-Analyte Set D Low Calibrator, 15 mL [§]
0545	Ethyl Glucuronide 500 High Calibrator (2000 ng/mL) [§]	0873	Multi-Analyte Set D Cutoff Calibrator, 15 mL [§]
0547	Ethyl Glucuronide 500 Level 1 Control (375 ng/mL) [§]	0874	Multi-Analyte Set D Intermediate Calibrator, 15 mL [§]
0548	Ethyl Glucuronide 500 Level 2 Control (625 ng/mL) [§]	0875	Multi-Analyte Set D High Calibrator, 15 mL [§]
0560	Fentanyl EIA Kit, small test kit	0877	Multi-Analyte Set D Level 1 Control, 15 mL [§]
0561	Fentanyl EIA Kit, large test kit	0878	Multi-Analyte Set D Level 2 Control, 15 mL [§]
0562	Norfentanyl Low Calibrator (2.5 ng/mL)	<i>Tests AMP 1000, BZO 200, COC 300, OPI 300, and OXY 100</i>	
0563	Norfentanyl Cutoff Calibrator (5 ng/mL)	0891	Multi-Analyte Set E Negative Calibrator, 15 mL [§]
0564	Norfentanyl Intermediate Calibrator (10 ng/mL)	0892	Multi-Analyte Set E Low Calibrator, 15 mL [§]
0565	Norfentanyl High Calibrator (20 ng/mL)	0893	Multi-Analyte Set E Cutoff Calibrator, 15 mL [§]
0567	Norfentanyl Level 1 Control (3.75 ng/mL)	0894	Multi-Analyte Set E Intermediate Calibrator, 15 mL [§]
0568	Norfentanyl Level 2 Control (6.25 ng/mL)	0895	Multi-Analyte Set E High Calibrator, 15 mL [§]
0610 510(k) Cleared	Oxycodone III EIA Kit, small test kit	0897	Multi-Analyte Set E Level 1 Control, 15 mL [§]
0611	Oxycodone III EIA Kit, large test kit	0898	Multi-Analyte Set E Urine DAU Level 2 Control, 15 mL [§]
OXY III 100		<i>Tests BARB 200, BUP 10, MTD 300, and PCP 25</i>	
0612	Oxycodone 100 Low Calibrator (50 ng/mL)		
0613	Oxycodone 100 Cutoff Calibrator (100 ng/mL)		
0614	Oxycodone 100 Intermediate Calibrator (150 ng/mL)		
0615	Oxycodone 100 High Calibrator (300 ng/mL)		
0617	Oxycodone 100 Level 1 Control (75 ng/mL)		
0618	Oxycodone 100 Level 2 Control (125 ng/mL)		
OXY III 300			
0622	Oxycodone 300 Low Calibrator (150 ng/mL)		
0623	Oxycodone 300 Cutoff Calibrator (300 ng/mL)		
0624	Oxycodone 300 Intermediate Calibrator (500 ng/mL)		
0625	Oxycodone 300 High Calibrator (800 ng/mL)		
0627	Oxycodone 300 Level 1 Control (225 ng/mL)		
0628	Oxycodone 300 Level 2 Control (375 ng/mL)		

510(k) Cleared = U.S. FDA 510(k) Cleared Assay, Calibrators, & Controls.[§] Product is made to order. Please allow 10 to 12 weeks of production time.
Minimum purchase volume required.

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂* A large test kit is 1000/375 mL R₁/R₂

**PRODUCT CATALOG: Drugs of Abuse Urinalysis – EU ONLY by Catalog Number**

REF #	Product Name
0990 	Fentanyl III EIA Kit (Q), small test kit ^s
0991	Fentanyl III EIA Kit (Q), large test kit ^s
0992	Fentanyl Low Calibrator/Level 1 Control (0.5 ng/mL) ^s
0993	Fentanyl Cutoff Calibrator (1 ng/mL) ^s
0994	Fentanyl Intermediate Calibrator (2 ng/mL) ^s
0995	Fentanyl High Calibrator (4 ng/mL) ^s
0998	Fentanyl Level 2 Control (1.5 ng/mL) ^s
1040	AB-PINACA (K3) EIA Kit, small test kit ^s
1041	AB-PINACA (K3) EIA Kit, large test kit ^s
0002c	THC Negative Calibrator (0 ng/mL)
1042	AB-PINACA (K3) Low Calibrator (2.5 ng/mL) ^s
1043	AB-PINACA (K3) Cutoff Calibrator (5 ng/mL) ^s
1044	AB-PINACA (K3) Intermediate Calibrator (10 ng/mL) ^s
1045	AB-PINACA (K3) High Calibrator (20 ng/mL) ^s
1046	AB-PINACA (K3) Level 1 Control (3.75 ng/mL) ^s
1047	AB-PINACA (K3) Level 2 Control (6.25 ng/mL) ^s

 = U.S. FDA 510(k) Cleared Assay, Calibrators, & Controls.^s Product is made to order. Please allow 10 to 12 weeks of production time.
Minimum purchase volume required.

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂* A large test kit is 1000/375 mL R₁/R₂