

**PRODUCT CATALOG: Drugs of Abuse Urinalysis by Alphabetical Order**

REF #	Product Name	REF #	Product Name
0001	Universal Negative Calibrator, 5 mL	0270	Buprenorphine EIA Kit, small test kit
0290	6-Acetylmorphine EIA Kit, small test kit	0271	Buprenorphine EIA Kit, large test kit
0291	6-Acetylmorphine EIA Kit, large test kit		BUP 5
0299	6-Acetylmorphine Negative Calibrator (0 ng/mL)	0272	Norbuprenorphine Level 1 Control (3 ng/mL)
0292	6-Acetylmorphine Low Calibrator (5 ng/mL)	0273	Norbuprenorphine Cutoff Calibrator #1 (5 ng/mL)
0293	6-Acetylmorphine Cutoff Calibrator (10 ng/mL)	0274	Norbuprenorphine Level 2 Control (7 ng/mL)
0294	6-Acetylmorphine Intermediate Calibrator (20 ng/mL)	0275	Norbuprenorphine Cutoff Calibrator #2 (10 ng/mL)
0295	6-Acetylmorphine High Calibrator (40 ng/mL)	0277	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)
0297	6-Acetylmorphine Level 1 Control (7.5 ng/mL)	0278	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)
0298	6-Acetylmorphine Level 2 Control (12.5 ng/mL)	0279	Norbuprenorphine High Calibrator (75 ng/mL)
0100	Amphetamines 500 EIA Kit, small test kit		BUP 10
0101	Amphetamines 500 EIA Kit, large test kit	0273	Norbuprenorphine Low Calibrator (5 ng/mL)
0102	Amphetamines 500 Low Calibrator (MAMP 250 ng/mL)	0274	Norbuprenorphine Level 1 Control (7 ng/mL)
0103	Amphetamines 500 Cutoff Calibrator (MAMP 500 ng/mL)	0275	Norbuprenorphine Cutoff Calibrator (10 ng/mL)
0104	Amphetamines 500 Intermediate Calibrator (MAMP 1000 ng/mL)	0276	Norbuprenorphine Level 2 Control (13 ng/mL)
0105	Amphetamines 500 High Calibrator (MAMP 2000 ng/mL)	0277	Norbuprenorphine Intermediate Calibrator #1 (20 ng/mL)
0107	Amphetamines 500 Level 1 Control (MAMP 375 ng/mL)	0278	Norbuprenorphine Intermediate Calibrator #2 (40 ng/mL)
0108	Amphetamines 500 Level 2 Control (MAMP 625 ng/mL)	0279	Norbuprenorphine High Calibrator (75 ng/mL)
0040	Amphetamines EIA Kit, small test kit	0070c	Cannabinoids (cTHC) EIA Kit, small test kit
0041	Amphetamines EIA Kit, large test kit	0071c	Cannabinoids (cTHC) EIA Kit, large test kit
0042	Amphetamines Low Calibrator (MAMP 500 ng/mL)	0002c	Cannabinoids (cTHC) Negative Calibrator (0 ng/mL)
0043	Amphetamines Cutoff Calibrator (MAMP 1000 ng/mL)		cTHC 25
0044	Amphetamines Intermediate Calibrator (MAMP 1500 ng/mL)	0072c	Cannabinoids (cTHC) Low Calibrator (12.5 ng/mL)
0045	Amphetamines High Calibrator (MAMP 2000 ng/mL)	0073c	Cannabinoids (cTHC) Cutoff Calibrator (25 ng/mL)
0047	Amphetamines Level 1 Control (MAMP 750 ng/mL)	0007c	Cannabinoids (cTHC) Intermediate Calibrator (37.5 ng/mL)
0048	Amphetamines Level 2 Control (MAMP 1250 ng/mL)	0075c	Cannabinoids (cTHC) High Calibrator (50 ng/mL)
0140	Barbiturate EIA Kit, small test kit	0006c	Cannabinoids (cTHC) Level 1 Control (18.75 ng/mL)
0141	Barbiturate EIA Kit, large test kit	0074c	Cannabinoids (cTHC) Level 2 Control (31.25 ng/mL)
	BARB 200		cTHC 50
0142	Barbiturate Calibrator #2/Control I (Secobarbital 100 ng/mL)	0073c	Cannabinoids (cTHC) Low Calibrator (25 ng/mL)
0143	Barbiturate Calibrator #3/Cutoff A (Secobarbital 200 ng/mL)	0075c	Cannabinoids (cTHC) Cutoff Calibrator (50 ng/mL)
0144	Barbiturate Calibrator #4/Control III (Secobarbital 300 ng/mL)	0076c	Cannabinoids (cTHC) Intermediate Calibrator (75 ng/mL)
0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)	0077c	Cannabinoids (cTHC) High Calibrator (100 ng/mL)
	BARB 300	0007c	Cannabinoids (cTHC) Level 1 Control (37.5 ng/mL)
0142	Barbiturate Calibrator #2 (Secobarbital 100 ng/mL)	0008c	Cannabinoids (cTHC) Level 2 Control (62.5 ng/mL)
0143	Barbiturate Calibrator #3/Control II (Secobarbital 200 ng/mL)		THC 100
0144	Barbiturate Calibrator #4/Cutoff B (Secobarbital 300 ng/mL)	0075c	Cannabinoids (cTHC) Low Calibrator (50 ng/mL)
0145	Barbiturate Calibrator #5 (Secobarbital 1000 ng/mL)	0077c	Cannabinoids (cTHC) Cutoff Calibrator (100 ng/mL)
0146	Barbiturate Control IV (Secobarbital 400 ng/mL)	0078c	Cannabinoids (cTHC) Intermediate Calibrator (150 ng/mL)
		0079c	Cannabinoids (cTHC) High Calibrator (200 ng/mL)
0130	Benzodiazepine EIA Kit, small test kit	0076c	Cannabinoids (cTHC) Level 1 Control (75 ng/mL)
0131	Benzodiazepine EIA Kit, large test kit	0009c	Cannabinoids (cTHC) Level 2 Control (125 ng/mL)
	BZO 200	0370	Carisoprodol Metabolite (SOMA) EIA Kit, small test kit [§]
0132	Benzodiazepine Calibrator #2/Control I (Oxazepam 100 ng/mL)	0371	Carisoprodol Metabolite (SOMA) EIA Kit, large test kit [§]
0133	Benzodiazepine Calibrator #3/Cutoff A (Oxazepam 200 ng/mL)	0372	Meprobamate Low Calibrator (50 ng/mL) [§]
0134	Benzodiazepine Calibrator #4/Control III (Oxazepam 300 ng/mL)	0373	Meprobamate Cutoff Calibrator (100 ng/mL) [§]
0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)	0374	Meprobamate Intermediate Calibrator (200 ng/mL) [§]
	BZO 300	0375	Meprobamate High Calibrator (400 ng/mL) [§]
0132	Benzodiazepine Calibrator #2 (Oxazepam 100 ng/mL)	0377	Meprobamate Level 1 Control (75 ng/mL) [§]
0133	Benzodiazepine Calibrator #3/Control II (Oxazepam 200 ng/mL)	0378	Meprobamate Level 2 Control (125 ng/mL) [§]
0134	Benzodiazepine Calibrator #4/Cutoff B (Oxazepam 300 ng/mL)		
0135	Benzodiazepine Calibrator #5 (Oxazepam 1000 ng/mL)	0340	Cocaine Metabolite 150 EIA Kit, small test kit
0136	Benzodiazepine Control IV (Oxazepam 400 ng/mL)	0341	Cocaine Metabolite 150 EIA Kit, large test kit
		0342	Cocaine Metabolite 150 Low Calibrator (BE 75 ng/mL)
		0343	Cocaine Metabolite 150 Cutoff Calibrator (BE 150 ng/mL)
		0344	Cocaine Metabolite 150 Intermediate Calibrator (BE 300 ng/mL)
		0345	Cocaine Metabolite 150 High Calibrator (BE 1000 ng/mL)
		0347	Cocaine Metabolite 150 Level 1 Control (BE 112.5 ng/mL)
		0348	Cocaine Metabolite 150 Level 2 Control (BE 187.5 ng/mL)

= FDA 510(k) Cleared Assay, Calibrators, & Controls

† Forensic Use Only

§ Product is made to order. Please allow 10 to 12 weeks of production time.
Minimum purchase volume required.

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂* A large test kit is 1000/375 mL R₁/R₂

PRODUCT CATALOG: Drugs of Abuse Urinalysis by Alphabetical Order

REF #	Product Name	REF #	Product Name
0030 	Cocaine Metabolite EIA Kit, small test kit	0970 	Fentanyl III EIA Kit (Q), small test kit ^s
0031	Cocaine Metabolite EIA Kit, large test kit	0971	Fentanyl III EIA Kit (Q), large test kit ^s
0032	Cocaine Metabolite Low Calibrator (BE 150 ng/mL)	0973	Fentanyl Cutoff Calibrator (1 ng/mL) ^s
0033	Cocaine Metabolite Cutoff Calibrator (BE 300 ng/mL)	0977	Fentanyl Level 1 Control (0.5 ng/mL) ^s
0034	Cocaine Metabolite Intermediate Calibrator (BE 1000 ng/mL)	0978	Fentanyl Level 2 Control (1.5 ng/mL) ^s
0035	Cocaine Metabolite High Calibrator (BE 3000 ng/mL)		
0037	Cocaine Metabolite Level 1 Control (BE 225 ng/mL)	0380 	Hydrocodone 100 EIA Kit, small test kit
0038	Cocaine Metabolite Level 2 Control (BE 375 ng/mL)	0381	Hydrocodone 100 EIA Kit, large test kit
0520 	Cotinine II EIA Kit, small test kit	0382	Hydrocodone 100 Low Calibrator (50 ng/mL)
0521	Cotinine II EIA Kit, large test kit	0383	Hydrocodone 100 Cutoff Calibrator (100 ng/mL)
0522	Cotinine II Low Calibrator (100 ng/mL)	0384	Hydrocodone 100 Intermediate Calibrator (150 ng/mL)
0523	Cotinine II Cutoff Calibrator (200 ng/mL)	0385	Hydrocodone 100 High Calibrator (300 ng/mL)
0524	Cotinine II Intermediate Calibrator (400 ng/mL)	0387	Hydrocodone 100 Level 1 Control (75 ng/mL)
0525	Cotinine II High Calibrator (1000 ng/mL)	0388	Hydrocodone 100 Level 2 Control (125 ng/mL)
0527	Cotinine II Level 1 Control (150 ng/mL)		
0528	Cotinine II Level 2 Control (250 ng/mL)	0390 	Hydrocodone 300 EIA Kit, small test kit
0160 	Ecstasy (MDMA) EIA Kit, small test kit	0391	Hydrocodone 300 EIA Kit, large test kit
0161	Ecstasy (MDMA) EIA Kit, large test kit	0392	Hydrocodone 300 Low Calibrator (150 ng/mL)
0162	Ecstasy Low Calibrator (100 ng/mL)	0393	Hydrocodone 300 Cutoff Calibrator (300 ng/mL)
0163	Ecstasy Cutoff Calibrator (500 ng/mL)	0394	Hydrocodone 300 Intermediate Calibrator (500 ng/mL)
0164	Ecstasy Intermediate Calibrator (750 ng/mL)	0395	Hydrocodone 300 High Calibrator (800 ng/mL)
0165	Ecstasy High Calibrator (1000 ng/mL)	0397	Hydrocodone 300 Level 1 Control (225 ng/mL)
0167	Ecstasy Level 1 Control (375 ng/mL)	0398	Hydrocodone 300 Level 2 Control (375 ng/mL)
0168	Ecstasy Level 2 Control (625 ng/mL)		
0220 	Ethyl Alcohol Enzymatic Assay Kit, small test kit	0640	Ketamine EIA Kit, small test kit ^{ts}
0221	Ethyl Alcohol Enzymatic Assay Kit, large test kit	0641	Ketamine EIA Kit, large test kit ^{ts}
0223	Ethyl Alcohol Calibrator (100 mg/dL)	0642	Norketamine Low Calibrator (25 ng/mL) ^{ts}
0224	Ethyl Alcohol Level 1 Control (50 mg/dL)	0643	Norketamine Cutoff Calibrator #1 (50 ng/mL) ^{ts}
0225	Ethyl Alcohol Level 2 Control (300 mg/dL)	0644	Norketamine Cutoff Calibrator #2 (100 ng/mL) ^{ts}
		0645	Norketamine Intermediate Calibrator (250 ng/mL) ^{ts}
0450	Ethyl Glucuronide III EIA Kit, small test kit ^{ts}	0646	Norketamine High Calibrator (500 ng/mL) ^{ts}
0451	Ethyl Glucuronide III EIA Kit, large test kit ^{ts}		KET 50
EtG 200		0647	Norketamine Level 1 Control (37.5 ng/mL) ^{ts}
0452	Ethyl Glucuronide 200 Low Calibrator (100 ng/mL) ^{ts}	0648	Norketamine Level 2 Control (62.5 ng/mL) ^{ts}
0453	Ethyl Glucuronide 200 Cutoff Calibrator (200 ng/mL) ^{ts}		KET 100
0454	Ethyl Glucuronide 200 Intermediate Calibrator (500 ng/mL) ^{ts}	0657	Norketamine Level 1 Control (75 ng/mL) ^{ts}
0455	Ethyl Glucuronide 200 High Calibrator (1000 ng/mL) ^{ts}	0658	Norketamine Level 2 Control (125 ng/mL) ^{ts}
0457	Ethyl Glucuronide 200 Level 1 Control (150 ng/mL) ^{ts}		
0458	Ethyl Glucuronide 200 Level 2 Control (250 ng/mL) ^{ts}	0630 	Methadone II EIA Kit, small test kit
EtG 500		0631	Methadone II EIA Kit, large test kit
0462	Ethyl Glucuronide 500 Low Calibrator (200 ng/mL) ^{ts}	0112	Methadone Low Calibrator (150 ng/mL)
0463	Ethyl Glucuronide 500 Cutoff Calibrator (500 ng/mL) ^{ts}	0113	Methadone Cutoff Calibrator (300 ng/mL)
0464	Ethyl Glucuronide 500 Intermediate Calibrator (1000 ng/mL) ^{ts}	0114	Methadone Intermediate Calibrator (600 ng/mL)
0465	Ethyl Glucuronide 500 High Calibrator (2000 ng/mL) ^{ts}	0115	Methadone High Calibrator (1000 ng/mL)
0467	Ethyl Glucuronide 500 Level 1 Control (375 ng/mL) ^{ts}	0117	Methadone Level 1 Control (225 ng/mL)
0468	Ethyl Glucuronide 500 Level 2 Control (625 ng/mL) ^{ts}	0118	Methadone Level 2 Control (375 ng/mL)
0310 	Fentanyl EIA Kit (Q), small test kit	0190b 	Methadone Metabolite (EDDP) EIA Kit, small test kit
0311	Fentanyl EIA Kit (Q), large test kit	0191b	Methadone Metabolite (EDDP) EIA Kit, large test kit
0313	Norfentanyl Cutoff Calibrator (5 ng/mL)		EDDPb 100
0317	Norfentanyl Level 1 Control (3.75 ng/mL)	0202b	Methadone Metabolite 100 Low Calibrator (50 ng/mL)
0318	Norfentanyl Level 2 Control (6.25 ng/mL)	0203b	Methadone Metabolite 100 Cutoff Calibrator (100 ng/mL)
0550	Fentanyl EIA Kit (SQ), small test kit ^{ts}	0204b	Methadone Metabolite 100 Intermediate Calibrator (250 ng/mL)
0551	Fentanyl EIA Kit (SQ), large test kit ^{ts}	0205b	Methadone Metabolite 100 High Calibrator (500 ng/mL)
0552	Norfentanyl Low Calibrator (2.5 ng/mL) ^{ts}	0207b	Methadone Metabolite 100 Level 1 Control (75 ng/mL)
0553	Norfentanyl Cutoff Calibrator (5 ng/mL) ^{ts}	0208b	Methadone Metabolite 100 Level 2 Control (125 ng/mL)
0554	Norfentanyl Intermediate Calibrator (10 ng/mL) ^{ts}		EDDPb 300
0555	Norfentanyl High Calibrator (20 ng/mL) ^{ts}	0192b	Methadone Metabolite 300 Low Calibrator (150 ng/mL)
0557	Norfentanyl Level 1 Control (3.75 ng/mL) ^{ts}	0193b	Methadone Metabolite 300 Cutoff Calibrator (300 ng/mL)
0558	Norfentanyl Level 2 Control (6.25 ng/mL) ^{ts}	0194b	Methadone Metabolite 300 Intermediate Calibrator (600 ng/mL)
		0195b	Methadone Metabolite 300 High Calibrator (1000 ng/mL)
		0197b	Methadone Metabolite 300 Level 1 Control (225 ng/mL)
		0198b	Methadone Metabolite 300 Level 2 Control (375 ng/mL)

 = FDA 510(k) Cleared Assay, Calibrators, & Controls[†] Forensic Use Only^s Product is made to order. Please allow 10 to 12 weeks of production time.
Minimum purchase volume required.

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂* A large test kit is 1000/375 mL R₁/R₂

PRODUCT CATALOG: Drugs of Abuse Urinalysis by Alphabetical Order

REF #	Product Name	REF #	Product Name
0350 	Methamphetamine EIA Kit, small test kit	0625	Oxycodone 300 High Calibrator (800 ng/mL)
0351	Methamphetamine EIA Kit, large test kit	0627	Oxycodone 300 Level 1 Control (225 ng/mL)
0352	Methamphetamine Low Calibrator (250 ng/mL)	0628	Oxycodone 300 Level 2 Control (375 ng/mL)
0353	Methamphetamine Cutoff Calibrator (500 ng/mL)		
0354	Methamphetamine Intermediate Calibrator (1000 ng/mL)	0010 	Phencyclidine EIA Kit, small test kit
0355	Methamphetamine High Calibrator (2000 ng/mL)	0011	Phencyclidine EIA Kit, large test kit
0357	Methamphetamine Level 1 Control (375 ng/mL)	0012	Phencyclidine Low Calibrator (12.5 ng/mL)
0358	Methamphetamine Level 2 Control (625 ng/mL)	0013	Phencyclidine Cutoff Calibrator (25 ng/mL)
		0014	Phencyclidine Intermediate Calibrator (50 ng/mL)
0852 	Multi-Analyte Set C Low Calibrator, 15 mL [§]	0015	Phencyclidine High Calibrator (100 ng/mL)
0853	Multi-Analyte Set C Cutoff Calibrator, 15 mL [§]	0017	Phencyclidine Level 1 Control (18 ng/mL)
0854	Multi-Analyte Set C Intermediate Calibrator, 15 mL [§]	0018	Phencyclidine Level 2 Control (32 ng/mL)
0855	Multi-Analyte Set C High Calibrator, 15 mL [§]		
0857	Multi-Analyte Set C Level 1 Control, 15 mL [§]	0120 	Propoxyphene EIA Kit, small test kit
0858	Multi-Analyte Set C Level 2 Control, 15 mL [§]	0121	Propoxyphene EIA Kit, large test kit
<i>Tests COC 150, MAMP, AMP 500, MTD, OPI2K, BZO 200, and BARB 200</i>		0122	Propoxyphene Low Calibrator (150 ng/mL)
0871	Multi-Analyte Set D Negative Calibrator, 15 mL [§]	0123	Propoxyphene Cutoff Calibrator (300 ng/mL)
0872	Multi-Analyte Set D Low Calibrator, 15 mL [§]	0124	Propoxyphene Intermediate Calibrator (600 ng/mL)
0873	Multi-Analyte Set D Cutoff Calibrator, 15 mL [§]	0125	Propoxyphene High Calibrator (1000 ng/mL)
0874	Multi-Analyte Set D Intermediate Calibrator, 15 mL [§]	0127	Propoxyphene Level 1 Control (225 ng/mL)
0875	Multi-Analyte Set D High Calibrator, 15 mL [§]	0128	Propoxyphene Level 2 Control (375 ng/mL)
0877	Multi-Analyte Set D Level 1 Control, 15 mL [§]		
0878	Multi-Analyte Set D Level 2 Control, 15 mL [§]	0400	SPICE I (JWH-018) EIA Kit, small test kit [†]
<i>Tests AMP 1000, BZO 200, COC 300, OPI 300, and OXY 100</i>		0401	SPICE I (JWH-018) EIA Kit, large test kit [†]
0891	Multi-Analyte Set E Negative Calibrator, 15 mL [§]	0002c	THC Negative Calibrator (0 ng/mL)
0892	Multi-Analyte Set E Low Calibrator, 15 mL [§]	0402	SPICE I (JWH-018) Low Calibrator (10 ng/mL) [†]
0893	Multi-Analyte Set E Cutoff Calibrator, 15 mL [§]	0403	SPICE I (JWH-018) Cutoff Calibrator (20 ng/mL) [†]
0894	Multi-Analyte Set E Intermediate Calibrator, 15 mL [§]	0404	SPICE I (JWH-018) Intermediate Calibrator (35 ng/mL) [†]
0895	Multi-Analyte Set E High Calibrator, 15 mL [§]	0405	SPICE I (JWH-018) High Calibrator (50 ng/mL) [†]
0897	Multi-Analyte Set E Level 1 Control, 15 mL [§]	0407	SPICE I (JWH-018) Level 1 Control (15 ng/mL) [†]
0898	Multi-Analyte Set E Level 2 Control, 15 mL [§]	0408	SPICE I (JWH-018) Level 2 Control (25 ng/mL) [†]
<i>Tests BARB 200, BUP 10, MTD 300, and PCP 25</i>			
0020b 	Opiate EIA Kit, small test kit	0420	SPICE II (UR-144/XLR-11) EIA Kit, small test kit [†]
0021b	Opiate EIA Kit, large test kit	0421	SPICE II (UR-144/XLR-11) EIA Kit, large test kit [†]
0022	Opiate Low Calibrator (Morphine 150 ng/mL)	0002c	THC Negative Calibrator (0 ng/mL)
0023	Opiate Cutoff Calibrator (Morphine 300 ng/mL)	0422	SPICE II (UR-144/XLR-11) Low Calibrator (10 ng/mL) [†]
0024	Opiate Intermediate Calibrator (Morphine 600 ng/mL)	0423	SPICE II (UR-144/XLR-11) Cutoff Calibrator (20 ng/mL) [†]
0025	Opiate High Calibrator (Morphine 1000 ng/mL)	0424	SPICE II (UR-144/XLR-11) Intermediate Cal (35 ng/mL) [†]
0027	Opiate Level 1 Control (Morphine 225 ng/mL)	0425	SPICE II (UR-144/XLR-11) High Calibrator (50 ng/mL) [†]
0028	Opiate Level 2 Control (Morphine 375 ng/mL)	0427	SPICE II (UR-144/XLR-11) Level 1 Control (15 ng/mL) [†]
		0428	SPICE II (UR-144/XLR-11) Level 2 Control (25 ng/mL) [†]
0330 	Opiate 2000 EIA Kit, small test kit		
0331	Opiate 2000 EIA Kit, large test kit	0410 	Tramadol EIA Kit, small test kit [§]
0332	Opiate 2000 Low Calibrator (Morphine 1000 ng/mL)	0411	Tramadol EIA Kit, large test kit [§]
0333	Opiate 2000 Cutoff Calibrator (Morphine 2000 ng/mL)	0412	Tramadol Low Calibrator (50 ng/mL) [§]
0334	Opiate 2000 Intermediate Calibrator (Morphine 4000 ng/mL)	0413	Tramadol Cutoff Calibrator (100 ng/mL) [§]
0335	Opiate 2000 High Calibrator (Morphine 6000 ng/mL)	0414	Tramadol Intermediate Calibrator (225 ng/mL) [§]
0337	Opiate 2000 Level 1 Control (Morphine 1500 ng/mL)	0415	Tramadol High Calibrator (400 ng/mL) [§]
0338	Opiate 2000 Level 2 Control (Morphine 2500 ng/mL)	0417	Tramadol Level 1 Control (75 ng/mL) [§]
		0418	Tramadol Level 2 Control (125 ng/mL) [§]
0610 	Oxycodone III EIA Kit, small test kit		
0611	Oxycodone III EIA Kit, large test kit		
	OXY III 100		
0612	Oxycodone 100 Low Calibrator (50 ng/mL)		
0613	Oxycodone 100 Cutoff Calibrator (100 ng/mL)		
0614	Oxycodone 100 Intermediate Calibrator (150 ng/mL)		
0615	Oxycodone 100 High Calibrator (300 ng/mL)		
0617	Oxycodone 100 Level 1 Control (75 ng/mL)		
0618	Oxycodone 100 Level 2 Control (125 ng/mL)		
	OXY III 300		
0622	Oxycodone 300 Low Calibrator (150 ng/mL)		
0623	Oxycodone 300 Cutoff Calibrator (300 ng/mL)		
0624	Oxycodone 300 Intermediate Calibrator (500 ng/mL)		

 = FDA 510(k) Cleared Assay, Calibrators, & Controls[†] Forensic Use Only[§] Product is made to order. Please allow 10 to 12 weeks of production time.
Minimum purchase volume required.

* EIA = Enzyme Immunoassay

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 100/37.5 mL R₁/R₂* A large test kit is 1000/375 mL R₁/R₂